

HISTOLOGY PRODUCT GUIDE

1.800.331.2716

#### **ORDERING INFO:**

- Call Monday-Friday, 8am-6pm EST to speak directly with one of our customer service representatives.
- Closed Saturdays & Sundays (Online: Open 24/7)
- · Voice Messaging Available

#### **RETURN POLICY:**

We pride ourselves in providing quality medical products at an outstanding value, and our goal is for you to be completely satisfied with your order. To ensure quality and maintain integrity in the event of damaged or discrepant orders, please review the following guidelines carefully.

#### SHIPPING DISCREPANCIES:

Please examine your packages upon receipt to ensure your order is correct. This can be done by identifying the box with the enclosed packing slip. Please keep in mind that with multiple box shipments the packing list will be in only one box. It is your responsibility to confirm that you have received all of the items indicated on your packing slip. Failure to report shipping discrepancies within five business days prevents us from accurately problem-solving the issue; therefore variances reported after 5 business days will not be credited. If the product you order is a custom order and requires a longer lead time than 7 days, a customer service representative will notify you.


#### **RETURNS:**

All return requests must be called into our **Customer Care Center** at 1.800.331.2716 prior to return and must be assigned an **RMA** (**Return Material Authorization**) number. If applicable, a FedEx pickup will be scheduled at time of authorization. Boxes containing product that have been written on or compromised in any way will not be accepted for return. All products must be returned in the **original packaging and packed in an additional outer box** to reduce the possibility of damage during return shipping.

- Returns may be subject to a minimum
   15% restocking charge.
- Non-refrigerated product must be returned within 30 days of receipt. Unfortunately, we can <u>NOT</u> accept returns for refrigerated product.
- Any product refused or returned without an authorized RMA will be subject to a restocking fee or denial of credit or replacement.
- Mercedes Scientific reserves the right to deny a return request outside of the above guidelines. After 30 days, the RMA and return option are no longer valid and will NOT be authorized or accepted.

EMAIL: ORDERS@MERCEDESSCIENTIFIC.COM

• Capital equipment subject to additional terms.



## **TABLE OF CONTENTS**

Equipment	2
Specimen Collection	33
Grossing	37
Marking	47
Embedding	52
Staining	62
Cassettes	68
Microscopy	79
Cytology	97
Chemicals	104
Labware	125
Safety	134
Spill Clean Up	148
Weights & Measures	157
Labels	162
Apparel	169
Index	177



### **About Mercedes Scientific**

Mercedes Scientific has established itself as the leader in laboratory supplies and equipment distribution with an unparalleled commitment to customer service. Specializing in sourcing exceptional products domestically and internationally, Mercedes Scientific maintains close partnerships with manufacturers and partners to ensure customers receive quality products at some of the most competitive pricing in the industry. Founded in 1991 and located in Lakewood Ranch, Florida, Mercedes Scientific is privately held, woman-owned, and a two-time winner of Inc. 500s Fastest-Growing Companies in America. The company continues to innovate to bring new levels of service excellence to the life sciences, clinical, histology, toxicology, and drug screening markets.



Alex Miller, CEO





**800.331.2716**Weekdays, 8:00am-6:00pm ET



**800.359.8807** All Day, Every Day



**MercedesScientific.com** User-Friendly & Secure, 24/7



orders@mercedesscientific.com All Day, Every Day

While every effort has been made to ensure that the material contained herein is accurate at time of publication, Mercedes Scientific is not responsible for typographical or photographic errors, inaccuracies or misprints contained within printed materials, web site content, web-based advertisements, email advertisements and/or online catalogs. Specifications, colors, legal restrictions, pricing and availability of all products is subject to change or discontinuation at any time without prior notice. Mercedes Scientific reserves the right to change prices, or not sell an item, due to error in pricing. It may be necessary to substitute an item with a product of equal or better quality and value that varies slightly from the print and/or online description. Please contact one of our customer care representatives at 800.331.2716 to confirm pricing. All products and service marks contained within Mercedes Scientific's printed materials, web site content, web-based advertisements, email advertisements and/or online catalogs are the property of their respective owners. References that Mercedes Scientific may make to any names, marks, products or services of third parties do not necessarily constitute or imply Mercedes Scientific's affiliation with, endorsement of, or endorsement by the third party, product or service.

# **TN50**

# **Semi-Automatic Cryostat**

### **Features:**

- · Exceptional German engineering
- Ideal for Mohs surgery and traditional pathology frozen section labs
- Modern rotary microtome with motorized µm advance/return
- Spacious cooling chamber for efficient and safe workflow
- · Large anti-roll plate
- 24 easily accessible freezing positions
- Heated viewing window (glass, electrically, conductive coating, transparent)
- Intuitive, simple electronic control panel
- · Automatic defrost
- · Removable condensate bottle with level sensor
- Optional UV-C disinfection solution

# **Cryochamber Specifications**

- Material: Stainless Steel
- Dimensions: 191/2"W x 213/4"D x 153/4"H
- · Illumination: LED, White
- Temperature Range: 0°C to -35°C
- Sample Preparation Positions: 24 Positions,

Temperature Range: -40°C to -45°C









# **TN55**



**Semi-Automatic Cryostat** with Quick Freeze -60°C Cold Spot

• Dimensions (WxDxH): 26¾ " x 235% " x 43¾"

• Weight (w/o accessories): 280 lbs

• Electrical Parameters: 115V; 50/60 Hz

#### **Features:**

- 24 easily accessible freezing positions + 2 quick-freeze positions (-50°C to -55°C)
- · Exceptional German engineering
- · Ideal for Mohs surgery and traditional pathology applications
- Modern rotary microtome with motorized µm advance/return
- · Spacious cooling chamber for efficient and safe workflow
- · Large anti-roll plate
- · Heated viewing window (glass, electrically, conductive coating, transparent)
- · Intuitive electronic control panel
- · Automatic defrost
- · Removable condensate bottle with level sensor
- · Optional UV-C disinfection light

# **Cryochamber Specifications**

Material: Stainless Steel

· Illumination: LED. White

• Temperature Range: 0°C to -35°C

• Sample Preparation Positions: 24 Positions, Temperature Range: -40°C to -45°C









# TN60 Cryostat with Object Head Cooling

The **NEW** German-made **Tanner Scientific® TN60 Cryostat** features a smooth microtome for easy, intuitive operation. For extra convenience the system includes a spacious stainless steel chamber, 24 easily accessible freezing positions and a well illuminated cooling chamber that can be effectively cooled down to -35°C. For additional user safety the system is equipped with an efficient **UV-C disinfection system**.

#### **Features**

- · Exceptional German engineering
- 2-Year Warranty\*
- · UV-C disinfection
- Dual compressors for separately controlled specimen head-cooling (up to -50 °C)
- 24 easily accessible freezing positions
   + 2 quick-freeze positions
   (-50 °C to -55 °C)
- Removable Flat Mounting Surface Accessory
- · Smooth rotary microtome
- · Spacious cooling chamber
- Glass anti-roll plate
- · Heated viewing window
- Easy-to-reach electronic control panel
- · Automatic defrost
- Removable condensate bottle with level sensor

# **Cryochamber Specifications**

- Material: Stainless Steel
- Dimensions: 191/2"W x 213/4"D x 153/4"H
- Illumination: LED, White
- Temperature Range: 0°C to -35°C
- Sample Preparation Positions: 24 Positions, Temperature Range: -40°C to -45°C

### **Microtome Specifications**

- Section Thickness Range: 0.5μm to 100μm
- Section Thickness Selection:
   0.5 to 2μm in 0.5μm Steps
   2 to 20μm in 1.0μm Steps
   20 to 50μm in 2.0μm Steps
   50 to 100μm in 5.0μm Steps
- Trimming Thickness Setting Range: 0.5 to 750μm

MADE IN GERMANY

- Horizontal Specimen Feed: 28mm
- Vertical Stroke: 58mm
- Coarse Feed: Motorized (600μm/s + 3.000μm/s)
- Specimen Retraction:
- 0 200μm, Free Programmable
- Horizontal Orientation: 8°
- Vertical Orientation: 8°
- Z-Axis Orientation: 360°





## **Standard Features**

- Precision German engineering/manufacturing
- 3-component disposable blade holders
- Smooth handwheel
- Modern handwheel lock
- Specimen retraction upon backstroke
- Vibration free microtome base plate
- Spacious, integrated section waste tray
- Wide range of accessories

# **Optional Microtome Accessories**

- Choice of High or Low Profile Blade Holder
- Quick Load Cassette Clamp (Orientable)
- Universal Cassette Clamp (Orientable & Fixed)
- Standard Object Clamp (Orientable & Fixed)
- Super-Mega Cassette Clamp (Fixed)
- Foil Clamp
- Illuminated Magnifying Glass

## **Technical Specifications**

- Section Thickness Setting Range: 0.5 to 60μm
- Section Thickness Selection:
   0.5 to 2μm in 0.5μm Steps
   2 to 10μm in 1.0μm Steps
   10 to 60μm in 2.0μm Steps
- Trimming Thickness Setting Range: 10μm, 20μm, 30μm, 40μm
- Horizontal Specimen Feed: 28 mm
- Vertical Stroke: 60 mm
- Coarse Feed: TN6000 Manual, TN7000/TN8000 Automated
- Specimen Retraction: Automatic
- Horizontal Specimen Orientation: 8°
- Vertical Specimen Orientation: 8°

- Z-Axis Specimen Orientation: 360°
- **Dimensions:** 187/8"W x 24"D x 133/4"H
- Includes: Blade Holder, Waste Tray
- Weight: 73 lbs



# **TN7000**

Manual Sectioning with Automated Coarse Feed Feature

## **Individual Features**

- Multiple Cut Modes
- Universal Cassette Clamp (Orientable & Fixed)



WITH 2-YEAR WARRANTY!\*

## **Standard Features**

- Precision German engineering/manufacturing
- 3-component disposable blade holders
- · Smooth handwheel
- · Modern handwheel lock
- · Specimen retraction upon backstroke
- Vibration free microtome base plate
- · Spacious, integrated section waste tray
- Wide range of accessories
- · Acoustic signal alerts at end position
- Trimming function in 4 steps

## **Technical Specifications**

- Section Thickness Setting Range: 0.5 to 60μm
- Section Thickness Selection:
   0.5 to 2μm in 0.5μm Steps
   2 to 10μm in 1.0μm Steps
   10 to 60μm in 2.0μm Steps
- Trimming Thickness Setting Range: 10μm, 20μm, 30μm, 40μm
- Horizontal Specimen Feed: 28 mm
- Vertical Stroke: 60 mm
- Coarse Feed: Automatic
- Specimen Retraction: Automatic
- Horizontal Specimen Orientation: 8°
- Vertical Specimen Orientation: 8°
- Z-Axis Specimen Orientation: 360°

\* 2-Year Warranty:

2 years parts, 1 year labor

# **Optional Microtome Accessories**

TANNER\*

- · Choice of High or Low Profile Blade Holder
- · Quick Load Cassette Clamp (Orientable)
- · Standard Object Clamp (Orientable & Fixed)
- · Super-Mega Cassette Clamp (Fixed)
- Foil Clamp
- Illuminated Magnifying Glass
  - **Dimensions:** 181/8" W x 24" D x 133/4" H
  - Includes: Blade Holder, Waste Tray
  - Weight: 73 lbs

\* 2-Year Warranty: 2 years parts, 1 year labor

# TANNER\*

# **TN8000**

**Dual-Use Microtome** 

## **Individual Features**

- Multiple Cut Modes
- Variable Speed (Touch Slider)
- · Slow Down in Cutting Window (Feature)
- Universal Cassette Clamp (Orientable & Fixed)



## **Standard Features**

- Precision German engineering/manufacturing
- 3-component disposable blade holders
- Smooth handwheel
- Modern handwheel lock
- Specimen retraction upon backstroke
- Vibration free microtome base plate
- · Spacious, integrated section waste tray
- · Wide range of accessories
- · Acoustic signal alerts at end position
- Trimming function in 4 steps

# **Optional Microtome Accessories**

- Choice of High or Low Profile Blade Holder
- Quick Load Cassette Clamp (Orientable)
- Standard Object Clamp (Orientable & Fixed)
- Super-Mega Cassette Clamp (Fixed)
- Foil Clamp
- Illuminated Magnifying Glass

## **Technical Specifications**

- Section Thickness Setting Range: 0.5 to 60μm
- Section Thickness Selection:
   0.5 to 2μm in 0.5μm Steps
   2 to 10μm in 1.0μm Steps
   10 to 60μm in 2.0μm Steps
- Trimming Thickness Setting Range: 10μm, 20μm, 30μm, 40μm
- Horizontal Specimen Feed: 28 mm
- Vertical Stroke: 60 mm
- Coarse Feed: Automatic
- Specimen Retraction: Automatic
- Horizontal Specimen Orientation: 8°
- Vertical Specimen Orientation: 8°
- Z-Axis Specimen Orientation: 360°

- **Dimensions:** 187/8"W x 24"D x 133/4"H
- Includes: Blade Holder, Waste Tray
- Weight: 73 lbs

\* 2-Year Warranty: 2 years parts, 1 year labor



# **TN1600**

# TANNER\*

# Tissue Embedding Center (SMALLER FOOTPRINT)

### **Features:**

- · German engineered, German manufactured
- · Manual dispensing lever and foot pedal
- · Adjustable magnifier glass
- · Ergonomic hand rest
- · Spacious, heated paraffin collection tray
- 3.8 liter paraffin container
- · Cooling plate
- Two heating trays for 130 cassettes or 250 embedding molds
- · Ergonomic work area
- · LED light illumination
- · Intuitive control panel

## **Cold Plate Specifications**

- Cooled Area Temperature Range: Down to -15°C
- Cooled Area Capacity: 80 Cassettes, Anodized Aluminum, 14¾" x 10¾"
- Control: Microprocessor, Working Day, Temperature Programming, Working Time, Time and Date
- Dimensions: 15¾"W x 22½"D x 11¾"H
- Weight: 55 lbs

# **Technical Specifications**

· Temperature Range


Paraffin Container: Up to 80°C Cassette/Mold Storage Containers: Up to 80°C Working Area: Up to 70°C Cooling Spot: Down to -5°C

- Paraffin Container: 3.8 liters
- Cooling Spot: 2" x 2", Peltier Forceps Holders 2 x 3, Magnetic, Adjustable
- Illumination:

3 x LED Spots, Adjustable

- Magnifier: Plastic, Adjustable
- Paraffin Flow: Activation Manually or by Footswitch, Adjustable Paraffin Flow Rate
- Control: Microprocessor, Temperature Programming, Working Day, Working Time, Time and Date
- Work Area: 17¾" x 9½", Anodized Aluminum
- Dimensions: 173/4"W x 221/2"D x 113/4"H
- Weight: 51 lbs
- Electrical Parameters: 115V; 50-60 Hz









#### **DID YOU KNOW?**

The cold plate **TN1600B** allows for a fast and efficient cooling of paraffin blocks after embedding and prior to sectioning. The compressor-cooled cold plate, made of anodized aluminum, guarantees consistent temperature distribution and is easy to clean.





# **TN1700**



# Modular Paraffin Embedding Center (LARGER FOOTPRINT)

### **Features:**

- · German engineered, German manufactured
- · Spacious cooling plate
- · Adjustable magnifier glass
- · Ergonomic hand/arm rest
- · Programmable working hours
- · Adjustable forceps holders
- 5.5 liter paraffin container
- · Spacious and heated paraffin collection tray
- · Working area with adjustable LED illumination and integrated cold spot

**DID YOU KNOW?** 

The TN1700 features a spacious working area, large container capacity and best in class reliability. All modules are the same working height and can be arranged anyway you choose. The extendable modular 3-component system with separate temperature control elements consists of a dispenser, a preheating unit and a separate cooling plate.

#### **SPACIOUS**

Large working area with integrated cold place. The work space has a drainage system that guides the paraffin into the collection tray.

#### **UNDER CONTROL**

The larger capacity, 5.5-liter volume tank allows uninterrupted workflow in higher volume labs.



# Dispensing Module Specifications

• Temperature Range:

Paraffin Container: Up to 80°C (1°C Increments) Working Area: Up to 70°C (1°C Increments)

Cooling Spot: Down to -10°C

• Paraffin Container: 5.5 liters

• Cooling Spot: 2" x 2", Peltier

• Forceps Holders: 2 x 4, Magnetic, Adjustable

 Illumination: 3 x LED Spots, Fixed 3 x LED Spots, Adjustable

• Magnifier: Included, Easy to Orientate

 Paraffin Flow: Activation Manually or by Footswitch, Continuously Adjustable Flow Rate

Control: Microprocessor, Program Temperature,
 Working Day, Working Time, Time and Date

Work Area: 17<sup>3</sup>/<sub>4</sub>" x 9½"

Dimensions: 17¾"W x 22½"D x 11¾"H

• Weight: 47 lbs

• Electrical Parameters: 115V; 50-60 Hz



MADE IN GERMANY

# **Cooling Unit Specifications**

• Temperature Range: Down to -15°C

• Maximum Capacity: 80 Cassettes

• Control: Microprocessor, Program Temperature, Working Day, Working Time, Time, Date

• Work Area: 143/4" x 103/4"

• **Dimensions:** 15¾"W x 22½"D x 11¾"H

• Weight: 55 lbs

• Electrical Parameters: 115V; 50-60 Hz



# Prewarming Unit **Specifications**

 Temperature Range: Up to 80°C, 1°C Increments for Working Area

• Cassette Bath: 150 Cassette Capacity, Removable, 9¾" x 10¼" x 1¼"


 Control: Microprocessor, Program Temperature, Working Day, Working Time, Time and Date

 Mold Warmer: 500 Mold Capacity, 10¼" x 8¾" x 2¾"

• Dimensions: 113/4" x 221/2" x 113/4"

· Weight: 33 lbs

• **Power:** 230V; 50-60 Hz





\* 2-Year Warranty: 2 years parts, 1 year labor









ITEM # DPT DNTLOS2

# Donatello

### **Tissue Processor Series 2**

#### **Increased Productivity**

- Processes up to 405 samples
- Automated reagent management software

#### **Internet Accessible**

- · Access/troubleshoot remotely
- · Control processor from afar
- Eliminate costly service calls

#### **Safety**

- · Safe reagent parking
- Fume control
- Xylene free options

#### **Smart, Intuitive Software**

- User-friendly icon based software, with colored touch screen for quick and easy managing of the instrument
- Remote accessible
- RFID (Radio Frequency Identification) avoids mismatching human errors in reagent exchanges

#### **LED Active System**

 Green lights reveal reagent position to guide users when changing reagents

MADE IN ITALY






ITEM # DPT TIZIANO

# Tiziano

# **Automated Slide Stainer for Special Stains**

Tiziano is simple to operate, and dramatically increases efficiency and consistency with special stains.

- REDUCED REAGENT USAGE due to tank volume
- ENHANCE QUALITY: use with Special Stains Kits from Mercedes Scientific (ready-to-use, quality-controlled reagents) or any other brand
- SPACE-SAVING DESIGN:
  - Stand-alone unit can be fitted easily into the available workspace
  - Equipped with software user-friendly touch-screen tablet
- READY-TO-GO: Comes pre-programed with 40+ stain protocols




Acid Resistant Bacteria - Ziehl Nielsen acc. Fite




Colon - Alcian Blue pH 2.5 acc. Dorling



Hypophysis - P.A.S. acc. Pearse



Hotchiss-McManus



Muscle -Picromallory











ITEM # DPT GIOTTOS2

# Giotto

### **Automatic Multistainer**

- · For histology and cytology staining
- · Run different stain protocols simultaneously
- · Noiseless Robotic Arm is fast, efficient and quiet
- · Minimizes cross contamination while maximizing productivity and throughput
- · S3 (Simultaneous Staining Scheduler) software automatically determines the best reagent reservoir configuration possible for optimal efficiency
- Stain ID Technology automatically determines the protocol based on the slide rack the user loads - allowing continuous loading and eliminating errors
- Smart LED System warns operator with acoustic and visual signals
- Water Pressure System: sensor alerts for possible missing waterflow to prevent failure of hematoxylin differentiation in running water
- One of the highest throughput H&E stainers on the market

MADE IN ITALY



DIAPATH





ITEM # CES CS50F

#### **NON-GYN SLIDE PROCESSOR**

This liquid-based non-gyn cytology system features more cells per slide and does not need duplicate preps.

#### **DUAL FILTRATION TECHNOLOGY**

- Disposable filter designed to remove non-diagnostic material to optimize cell presentation
- Cells of interest are transferred to the slide in a homogeneous thin-layer cell deposit
- · Filters are available in 20mm and 14mm diameters, thereby improving cell yield and presentation

#### **CELLSOLUTIONS F50 SPECIFICATIONS**

- Throughput: 50 samples/hr. Semi-automated
- **Dimensions:** 11"W x 15"D x 8.5"H; **Weight:** 31 lbs

The CellSolutions double membrane technology yields cellular preparations with increased cellularity, presentation, and stain quality on 93.75% of the specimens tested when compared to other liquid based methods.1



See Page 100 for preservatives.

# **FINALLY** Something New After 20 Years!

## **FASTER**

Slide is ready to stain in 20 seconds.

### CHEAPER

Less expensive to acquire and operate. Costs a fraction of the others.

## BETTER

Slides will be better looking and easier to read.

1 Smith D., Espinoza, A., Vega, I. "Comparison Study of Traditional Cytology Non-Gyn Processing Techniques to Specimens Processed with CellSolutions Method". 2017 American Society of CytoPathology Meeting: Poster Publication



# **Microscope Slide Label Printer**


A complete, ready-to-print microscope slide label printer package that comes with everything you need to print histology grade slide labels — printer, software, labels, and ribbon.

The slide labels and thermal transfer ribbon are optimized for harsh chemicals and will not absorb eosin or hematoxylin coloring.

SlideID printers function with or without LIS. They can be set up for batch printing or with on-demand printing for "laboratory best practices."

#### **FEATURES:**

- Ready-to-print
- Scalable
- · Resistant to stains, solvents and xylene
- · Will not absorb eosin or hematoxylin coloring
- · Function with or without LIS
- Can be setup for batch printing or with on-demand printing



QRL SLIDEIDPRINTERPK	BASIC PACKAGE: SlideID CL-E331 Slide Label Printer, SlideID Slide Printing Software, SlideID Microscope Slide Labels, SlideID Thermal Transfer Ribbon	Each
QRL SLIDEIDPROPK	PRO PACKAGE: Includes All Contents of the Basic Package Above Plus Denso AT30Q Barcode Scanner	Each
QRL SLIDEIDLABEL	SlideID Microscope Slide Labels, 0.875" x 0.875", 1" core (1,750 Per Roll)	Roll
QRL SLIDEIDRIBBON	SlideID Special Resin Ribbon, 1.57" W x 984' (prints 12,000 +\- Slide labels)	Each
QRL SLIDEIDRIBBONZ	SlideID Thermal Transfer Ribbon For ZEBRA® Style Printers, Xylene and solvent resistant, histology grade	Each



#### PRIMERA® SIGNATURE® **SLIDE PRINTER**

#### **FEATURES & BENEFITS:**

- Reduce Inventory & Purchasing Activity: inventory of multiple colors of slides is eliminated since it permits printing colors on white frosted slides
- · Prints Directly Onto Slides: eliminates handwriting or difficultto-apply labels
- Small Footprint: fits perfectly at the workstation
- Slides Stored in Easy-to-Load Cartridges: protects from dust and other potential contaminants
- PTLab Software: template-driven software for exact information desired on slides gets printed; Pro version configures input/ output data streams to LIS systems
- Chemical-Resistant Ink: ensures reliable identification of slides now and for years in the future
- Includes: Printer, PTLab™ Software, TuffCoat Extreme Black Ribbon, Printer Driver, UL 60601-1 Power Supply, "Green Dot" Plug

PRI 78104	Primera® Slide Label Printer	Each
ACCESSORIES		
PRI 58105	Black Ribbon (5,000 Prints)	6/Box
PRI 78251	Spare Slide Holder (100-Count)	Fach



#### PRIMERA® SIGNATURE® **CASSETTE PRINTER**


#### **FEATURES & BENEFITS:**

- · Designed for efficient, hands-free operation
- Printing cassettes on-demand or in batch mode reduces errors and increases efficiency
- Taped cassettes are loaded into individual hoppers
- Level of cassettes remaining is easily viewed through the front of the hopper's transparent blue LED backlit case
- Black or color printing
- (4) Cassette Hoppers, (40) Cassettes Per Hopper
- · Prints on standard cassettes with 35° to 45° angles with either lid on or lid off, no adjustments needed
- · Prints up to 8 cassettes per minute
- · Thermal transfer printing at 300 dpi
- · Chemical-resistant ink for cassette identification
- Ribbons: Black (4,300 prints), CMYK (550 prints)
- Output Capacity: 7 Cassettes, 17 with Extension Tray
- USB Interface: ability to integrate with Laboratory Information Systems

PRI 78405	Manual Cassette Label Printer	Each
PRI 78407	Cassette Hopper Auto Loader	Each

WEB MercedesScientific.com

19







- Thermal heater ensures heat transfer over entire surface
- · Anodized black surface provides contrast with samples
- · Simple digital temperature controller
- Temperature Range: Room Temperature to 75°C

CAS XH2001	XH-2001 Slide Warmer  Dimensions: 25"L x 8"W	Each
	Capacity: Holds 66 Slides	
CAS XH2002	XH-2002 Slide Warmer	Each

Dimensions:  $10\frac{1}{4}$ "L x 7"W Capacity: Holds 23 Slides



#### PREMIERE® STEP UP SLIDE WARMER

- Unique design accommodates twice the number of slides as a flat unit of the same size
- Easy slide retrieval since end of slide isn't on flat surface
- Anodized black surface provides contrast with samples
- Simple digital temperature controller
- · Removable tray for easy cleaning
- Temperature Range: Room Temperature to 75°C
- Dimensions: 14"L x 14"W x 4"HCapacity: Holds 40 Slides

CAS XH2003 Step Up Slide Warmer Each



#### PREMIERE® DIGITAL WAX TRIMMER

This Wax Trimmer has a compact and powerful build that is perfect for medical and research laboratories. It boasts a magnetically secured tray and a disposable tray liner that collects melted paraffin, making for easy and clean disposal. It is efficient and practical for removing excess paraffin from tissue cassettes. With this wax trimmer, there is no need to scrape the collection tray.

- Digital control panel monitors changes in heat level and to adjust temperature
- · Adjustable heat setting
- Trimming of multiple cassettes simultaneously to save operator time
- Red LED light alerts users when the unit is on and heating, ensuring user safety
- Temperature Range: Room Temp. to 120°C
- Power Requirements: 150W
- Overall Size (Including Tray): 11.25 x 5 x 8"

CAS XH90D XH-2001 Slide Warmer Each
Dimensions: 25"L x 8"W

Capacity: Holds 66 Slides

#### LINISTAT™ LINEAR STAINER

- · Compact answer to H&E frozen section protocols
- · Dependable, economical way to automate staining; uses a conveyer to prepare consistently stained slides
- · Can use for hand-staining when the situation requires
- Efficient footprint fits on top of most cryostats for streamlining workflow during frozen sectioning
- Includes: (14) Staining Stations (20 Seconds Each), Running Water and 13-Slide Capacity Collection Tank
- STAINING TIME IS UNDER FIVE MINUTES
- Color-coded slide clips correspond to Thermo Scientific Cryochrome™ colored slides and colored frozen section embedding media
- **Dimensions:** 6.25"D x 25"W x 18.4"H
- Clearance Height: 9.75"; Weight: 16.3 lbs


FIS B1000200

Linistat™ Linear Stainer

Fach







#### **HISTOPRO® 414 LINEAR STAINER**

- American made H&E stainer designed specifically for the needs of a Mohs surgery laboratory
- · Fits on top of most cryostats
- · Ideal for H&E stains up to 4 slides per slide carrier
- · Continuous loading feature to keep up with lab demand
- Programmable station timing from 2 to 300 seconds
- 14 reagent stations include 3 running rinse water stations
- · Holding station for 16 processed slides
- · Only 50mL per reagent container
- Up to 900 slides per hour at 10-second station time

RUS HISTOPRO414	H&E Linear Stainer	Each
RUS 13072	Reagent Containers	6/Pack
RUS 14509	Plastic Hood	Each







#### "G SERIES" DUCTLESS FUME HOODS

- · Portable tabletop design requires no outside venting
- Systems are designed to enclose both materials and operations under negative pressure to prevent airborne contaminants from escaping into the work environment
- Contaminated air is drawn through a variety of filters to remove pollutants
- Filters are available for both vapors and particulates and they can be combined for multiple substance filtration
- Hoods are constructed of clear, thick, polycarbonate and acrylic panels with sturdy aluminum framing
- Clear vinyl curtain provides easy access while maintaining optimal airflow across the opening

#### NOTE: FILTERS SOLD SEPARATELY

AIR G24	32½"H x 24"W x 18"D	Each
AIR G30*	32½"H x 30"W x 18"D	Each
AIR G50	32½"H x 50"W x 18"D	Each



#### **BENCHTOP FUME ADSORBER**

- · For light duty histology coverslipping and staining areas
- Activated carbon media adsorbs organic odors and noxious vapors released by solvents (ie. xylene, toluene) while clean air is recirculated into the laboratory
- Acrylic canopy enclosures are shatter-resistant, low profile, completely portable and require no ductwork
- Glacier white, dry powder, epoxy-coated steel housing with perforated baffle for uniform airflow
- Built-in, motorized impeller's face velocity up to 35 fpm
- · Uses activated carbon, formaldehyde or ammonia filter
- Exterior Dimensions: 15.4"H x 23.4"D
  Interior Dimensions: 13.5"H x 14.5"D

**LAC 3955200** Interior: 131/2"H x 141/2"D Each

### **MERCEDES SCIENTIFIC® CARBON FILTERS**

- · Sold individually, this line of activated carbon filters provide a high-value, low-cost generic alternative for both Leica® and Shandon® machine filters
- Cartridges contain a high quality extruded activated carbon
- Absorb most organic and inorganic vapors including xylene, toluene, methanol, formalin, isopropyl alcohol and industrial denatured alcohol
- · Ensure and maintain safe air quality in the laboratory
- · Low cost allows more frequent changing and improves personnel safety





#### **MER 44005**

**DIMENSIONS:** 8" x 71/2" x 7/8" For Leica® CV5030™ Coverslipper QTY: Each



#### MER 440124

**DIMENSIONS:** 5.5" x 1.6" x 11.2" For Leica® ASP300™ Tissue Processor QTY: Each



#### **MER 440020**

**DIMENSIONS:** 7.8" x 2.3" x 16.5" For Sakura® E150, E300, VIP®5 QTY: Each




#### **MER 440126**

**DIMENSIONS:** 3.3" × 1.1" × 24.8" For Sakura® Tissue-Tek® DRS™ 2000 Stainers QTY: Each



#### MER 440120

**DIMENSIONS:** 3.9" x 1.37" x 18.8" For Leica® Autostainer XL ST5010. ST5020. ST4040 QTY: Each



#### MER 440127

**DIMENSIONS:** 4.7" x 1.45" x 7.9" For Sakura® Tissue-Tek® Tape 4764 Coverslipper QTY: Each



#### MER 440122

**DIMENSIONS:** 5.7" × 1.37" × 5.7" For Leica® CV5000™ Coverslipper QTY: Each



#### **MER 440128**

**DIMENSIONS:** 3.9" x 1.9" x 13.7" For Sakura® Tissue-Tek® Prisma®, VIP®6 QTY: Each



#### MER 440123

**DIMENSIONS:** 5.6" x 1.85" x 11.3" For Leica® TP1050™ Tissue Processor QTY: Each



#### **MER 447101**

**DIMENSIONS:** 5.7" x 1.37" x 11.2" For Thermo® Shandon™ Hyperclean® Pathology Workstation, Excelsior™ ES QTY: Each

#### **ASSORTED CARBON FILTERS**

FORMALDEHYDE specialty carbon filters are treated carbons that combine a chemical reaction with carbon's adsorption capacity. Compounds are broken down into the various elements for easy adsorption. Dual stage with monitoring port feature (FRC9) detects vapor breakthrough between filters; offering further layer of protection to users.

**ORGANIC CARBON** filters are used in a variety of Labconco<sup>\*</sup> enclosures which require periodic replacement.



#### **AIR FRC4**

Formaldehyde Filter **DIMENSIONS:** 10" x 17" x 2"

For G-Series Hoods, Scavengers, HS3000/HS3000A Blowers

QTY: Each



**AIR FR4** 

Activated Charcoal Filter **DIMENSIONS:** 10" x 17" x 2" For G-Series Hoods, Scavengers, HS3000, RF1000

**QTY:** Each

**ACTIVATED CHARCOAL** or coconut shell carbon is effective in adsorbing most commonly found organic vapors and odors. It can adsorb up to 50% of its own weight. Standard activated carbon does not have much capacity for some reactive gases. Where gases are chemically reactive, specially impregnated carbon can be recommended.

**POTASSIUM PERMANGANATE** filters remove high levels of dense smoke, odors and other chemicals that charcoal can not.



#### LAC 3924200

Organic Carbon Filter **DIMENSIONS:** 18" x 36" x 1"

For Protector Work Stations **QTY:** Each



#### LAC 3981400

Trace Odor Organic Carbon Filter **DIMENSIONS:** 24" x 18" x 1"

For XPert & RXPert Balance and

XPert Bulk Powder Filtered Systems **QTY:** Each



#### MOP BF007

Potassium Permanganate Filter **DIMENSIONS:** 14" x 7" x 1" For Mopec\* BF Series Grossing/ Trimming Workstations

QTY: Each



#### MOP BF009

Potassium Permanganate Filter **DIMENSIONS:** 15<sup>3</sup>/<sub>4</sub>" x 4" x 3" For Mopec\* MB, ME, O-Series Grossing Workstations

QTY: Each

#### **GRANULAR ACTIVATED CHARCOAL**

 Hard, highly adsoptive charcoal is used in a wide variety of air and water filters in numerous applications

BAC GAC48C	4lbs Bottle	Each
EKI 3056500G	500g Bottle	Each



#### FORMALIN RECYCLING SYSTEM

- · Process 3 gallons in approximately 1 hour
- Recycled product is the same concentration and volume of formalin recycled through the system
- Recycles formalin without removing the buffered salts or altering the formalin concentration
- · Also recycles alcoholic and zinc formalin
- Simply pour dirty 10% Neutral Buffered Formalin into the top and clean recycled 10% Neutral Buffered Formalin at the same volume, concentration and pH exits the bottom

CWS FRS3GAL	Floor Model System	Each
CWS FRS1GAL	Bench Top Model	Each
CWS FRRF	Replacement Filter	Each

#### ALCOHOL RECYCLING SYSTEM

- Process 3 gallons in approximately 1 hour
- Recycled product is the same concentration and volume of alcohol recycled through the system
- · Recycle alcohol with no capital outlay
- Simply pour dirty alcohol into the top and clean recycled alcohol at the same concentration exits the bottom

CWS ARS3GAL	Floor Model System	Each
CWS ARS1GAL	Bench Top Model	Each
CWS 3GARRS	3 Gal Replacement Filter	Each

#### **RECYCLING SYSTEM FEATURES**

- Guaranteed to save laboratories money, time, energy and effort all while helping out the community's environment
- Safe, simple and easy to use, all the while requiring very little lab space
- Run on gravity, can be used anywhere, no electricity or heating of formalin required
- · Customizable to meet your specific lab needs

## AIR FILTRATION SYSTEM

- Minimizes exposure from hazardous fumes and odors caused by alcohol, xylene, formalin, solvents and others
- Requires very little space and is designed to be placed on the floor or counter
- · Improves air quality in laboratories
- Works quietly and efficiently 7 days/week, 24 hours/day
- Dimensions: 23"H x 11"W

CWS AFS001	Air Filtration System	Each
CWS AFSRE	Replacement Air Filter	Fach





DIMENSIONS: 14" x 14" x 21"



**FLOOR MODEL** 

DIMENSIONS: 22" x 25"





#### **ROUND PREPARATION TISSUE BATH**

- Edges serve as drying surface for processed slides
- Black ceramic coated steel surface for improved visibility of tissue sections
- Equipped with on/off switch, pilot light, and thermostat
- Includes: Dust Cover, Thermometer Holder, Rubber Feet
- Temperature Range: 35°C to 65°C
- Dimensions: 81/4" x 2" (Chamber), 11" x 33/4" (Overall)

**BLK 14702** Round Preparation Tissue Bath

Each

#### PREMIERE® ROUND TISSUE BATH

- Low profile, black, finished interior chamber reduces glare to improve visibility of tissue specimens
- · Glazed glass lid for easy viewing of interior
- Top ledge allows space for preheating or for evaporating away any remaining water on the glass slides
- Durable fiberglass insulated body prevents heat loss and fluctuations in temperature for uniform surface heating
- Temperature Range: Room Temperature to 75°C
- **Dimensions:** 10½" x 8¼" (Overall), 8¼" x 2½" (Interior)

CAS XH1001 Round Tissue Bath

Each

#### PREMIERE® LIGHTED TISSUE BATH

- Low profile design with removable rectangular glass basin, illuminated with LED light from the side
- Easy-to-read LED display shows both set temperature and actual temperature
- Sleeker glass basin for better contrast of specimens
- Temperature Range: Room Temperature to 75°C +/- 1°C
- Dimensions: 14 x 14 x 41/2" (Overall), 10 x 61/2 x 2" (Interior)

CAS XH1003 Lighted Tissue Bath Each

#### STANDARD TISSUE FLOTATION BATH

- Eliminates wrinkles and distortion during preparation of paraffin embedded tissue sections
- Digital microprocessor control with advanced safety over temperature protection
- LED illuminated black background with transparent, removable Pyrex<sup>®</sup> dish allows specimen to be easily seen
- Nonvolatile memory retains the last temperature set point after being powered on/off
- Temperature Range: 5°C to 65°C
- Overall Dimensions: 141/4"L x 141/4"W x 41/2"H

BLK 145702 Standard Tissue Flotation Bath Each



#### **INNOVATION BIOLOGICAL MICROSCOPES**

- · Premium choice for cytologists, doctors, and researchers
- · Powerful performance, crisp clarity, rugged construction, and backed by a lifetime warranty on optical and mechanical components
- · Infinity Flat-Field Objectives: 4X, 10X, 20X, 40X and 100X (Oil) is premium for high contrast and resolution
- · Siedentopf-type head is inclined 30° for a comfortable user posture
- · Large, graphite-coated mechanical stage withstands wears and tears and accommodates two slides at once
- Ergonomic, durable design ensures dependability and comfort for hours of daily use
- · Bright, white, uniform LED illumination maintains constant daylight-color to show stained specimens in their true vibrant colors

LWS INMB04AIPL3	Binocular, 4 Objectives	Each
LWS INMT04AIPL3	Trinocular, 4 Objectives	Each
LWS INMB05AIPL3	Binocular, 5 Objectives	Each
LWS INMT05AIPL3	Trinocular, 5 Objectives	Each



#### **BIOVID CAMERA & MONITOR**

- · Essential tool for live presentations and client education
- Premium BioVID HD1080+ camera paired with the 1080p 11.6" HD Monitor provides a bright and beautiful image
- Simply mount the camera onto any trinocular C-mount for amazing resolution and fast-action video (60fps)
- · Using the included mouse and SD card, click to save and view images and videos...even side-by-side comparisons
- BioVID HD 1080+ camera features 1920x1080 resolution, increased sensitivity, low-noise and great color rendition to create stunning, high-resolution imaging on the attached 1080p HD monitor, large HDTV or projector
- Extra large pixel size on the Sony IMX236(C) CMOS chip ensures crisp images in bright field (even at high magnifications) while providing the sensitivity needed to view and capture precise images in low light applications such as dark field, phase contrast or epi-fluorescence
- Choose from HDMI and/or USB outputs or save images and videos to the included SD Card

#### **IDEAL FOR A VARIETY OF APPLICATIONS**

- Research, Teaching and Demonstrations
- · Pathology, Cytology, Oncology, Hematology
- · Semen Evaluation
- · Dental Applications
- · Industrial and Manufacturing

LWS BVC1080TVK3	Camera & Monitor	Each
LWS BVC1080CMT3	Camera Only	Each

WEB MercedesScientific.com



27

#### **14 INFINITY MICROSCOPES**

- Exceptional optical quality and expandability for top-notch performance in the laboratory
- · Rotating head adjusts the height for taller users
- Infinity Semi-Plan or Plan Objectives are anti-fungal, parfocal, parcentric and color-coded
- 30-year anti-fungal coating on all lenses
- · Cool, long-life LED "daylight" illumination
- 4 Objectives: 4X, 10X, 40XR, 100XR (Oil)
- · Features reverse quadruple nosepiece
- Perfect choice for durability, dependability and superior imaging in the laboratory

LWS I4MBN4AIPL3	Binocular, Infinity Plan	Each
LWS I4MTN4AIPL3	Trinocular, Infinity Plan	Each
LWS I4MBN4AISL3	Binocular, Semi-Plan	Each
LWS I4MTN4AISL3	Trinocular, Semi-Plan	Each

#### **NOTE: OTHER ATTACHMENTS & OPTIONS AVAILABLE**

Phase Contrast, 50X Oil Objective, Eyepiece Counter, Digital Cameras and Adapters

# PREMIERE® PROFESSIONAL BINOCULAR MICROSCOPES

- · Objectives face backwards for more space
- Base has a comfortable plastic plating over the metal for less heat conduction
- Field diaphragm adjustment built into the illuminator offers Koehler illumination for more lighting control

#### **FEATURES**

- 30° inclined binocular Seidentopf head, rotatable 360°
- Pair of 10X widefield eyepieces, interpupillary distance adjustable from 55mm to 75mm
- Abbe condenser NA 1.25 and Iris diaphragm with filter holder including blue filter
- DIN Achromatic Parfocal and Parcentered Objectives, Backward Mounted: 4X (N.A. O.10), 10X (N.A. 0.25), 40X (N.A. 0.65), 100X (Oil Immersion N.A. 1.25) 60XR (N.A. 0.90) MRP5000 ONLY
- Magnifications: 40X to 1000X
- Five objectives and interchangeable condenser system applies to bright field, phase contrast, dark field and polarized light techniques MRP5000 ONLY
- Lower position of mechanical stage knob MRP5000 ONLY

CAS MRP3000	Microscope, 4 Objectives	Each
CAS MRP5000	Microscope, 5 Objectives	Each











#### PREMIERE® PARAFFIN DISPENSER

- Built-in temperature controller and heated faucet keeps paraffin flowing with ease
- NEW LID WITH HOLE FOR THERMOMETER!
- · Industrial metal outer shell ensures longer life
- · Durable stainless steel interior is easy to clean
- Maximum Temperature: 72°C tolerance ± 2°C)
- Interior Dimensions: 9"L x 10"H
- Overall Dimensions: 16"L x 15"W x 18"H

CAS XH4002 2.3GL Large Capacity Each

#### **BOEKEL® PARAFFIN DISPENSER**

- · Stainless steel outer casing
- · Rapid warm-up with an advanced safety over temperature cut-off
- · Heated valve ensures a consistent flow of molten wax
- · Inner reservoir screen helps to prevent large particulate from clogging the dispenser valve
- Aluminum inner reservoir provides excellent thermal transfer and temperature control
- · Dispenser valve allows precise control of flowing wax and can be locked in the open position if desired
- Temperature Range: 40°C to 75°C
- Overall Dimensions: 81/2"W x 16"H

**BLK 145600** 1.25GL (5L) Capacity Each

#### **DIGITAL PARAFFIN DISPENSER**

- · For use in the Anatomical Pathology Laboratory
- · Suitable for paraffin dispensing during the preparation and inclusion of paraffin blocks
- Includes: Procedure Manual, Power Cable, Spare Fuse, Warranty Certified
- Temperature Control: Environment Temperature up to 90°C with a range of ± 2°C through Digital Thermostat
- Overall Dimensions: 295mm x 310mm x 180mm

**CAN MX90** Up to 4L Capacity Each







The Boekel® Large Wax Dispenser is ideal for larger histology laboratories that require large volumes of paraffin wax when embedding tissue samples. The Large Wax Dispenser has a maximum capacity of 28 liters (7.4 gallons) of wax and is the largest in the industry. This device can be used as the primary source of paraffin wax or to supplement tissue processors and embedding centers. The unit has a heated and filtered spigot to make sure the wax flows freely during operation. The digital control system is easy to use and has a timer based Hot Start system that ensures the wax is melted prior to operation. The unit has a stainless-steel lid and chamber to stay corrosion free for years of use.

- · Industry's largest wax capacity for high throughput use
- Quickly Melts from pellets to molten wax in as little as 4 hours!
- Hot Start feature will automatically start the warming based on set time to eliminate waiting
- Bright easy to use digital display
- · Stainless steel lid and wax chamber to resist corrosion
- · Heated and filtered wax spigot to keep the wax flowing quickly

Fach

- Two-year Warranty
- UL, CE, CSA approved

**BLK 1456XL** 7.4GL (28L) Capacity





OCY 40GC





QCY 30GC

#### **GRAVITY CONVECTION OVENS**

- Used for part drying, baking, curing, sterilizing, evaporating, heat treating, annealing, and testing
- · Heavy steel double-wall construction
- Doors open with high impact thermoplastic handles
- · High-density mineral wool insulates work space
- Exteriors have a hard, scratch-resistant hammer finish
- Interiors are made of corrosive-resistant aluminized steel
- Resistive-tubular incoloy sheathed heating elements are engineered for quick run-up and recovery times
- Control panel is equipped with an illuminated ON/OFF rocker switch and a heat-cycle pilot light
- Temperature uniformity greatly improved by perforated heat shield to absorb and distribute radiant heat evenly
- Models have temperature range to 232°C (450°F)

QCY 10GC	10GC Model, 0.7 ft <sup>3</sup>	Each
QCY 20GC	20GC Model, 1.27 ft <sup>3</sup>	Each
QCY 30GC	30GC Model, 2.0 ft <sup>3</sup>	Each
QCY 40GC	40GC Model, 3.0 ft <sup>3</sup>	Each



## SPECIMEN DISCS (W/HANDLES)

- Made of aluminum with comfort grips
- Helps to flatten and freeze samples
- · Make transporting and handling specimens easier
- COMPATIBLE WITH ALL LEICA® CRYOSTAT MODELS

MER 30	30mm	Each
MER 40	40mm	Each



#### **CRYOSTAT SPECIMEN DISCS**

- · Made of aluminum with comfort grips
- Embossed circular discs with outer o-ring for easy handling and edge protection
- Used in conjunction with heat extractor
- Helps to flatten and freeze tissue samples

#### FOR USE ON LEICA® CRYOSTATS

MES DISC20	20mm	Each
MES DISC25	25mm	Each
MES DISC30	30mm	Each
MES DISC40	40mm	Each
MES DISC55	55mm	Each

# FOR USE ON AVANTIK® QS11/QS12 CRYOSTATS AND MICROM HM520/HM525 CRYOSTATS

MES CHUCK30	30mm	Each
MES CHUCK40	40mm	Each



#### **RED O-RING, SPECIMEN CHUCK**

- Flat round metal specimen chuck for Leica CM1950 cryostat
- With orientation mark
- · Low-temp silicon O-ring RED

SPG NLK30R 30mm Each



# MERCEDES SCIENTIFIC® LUBRICATING OIL

- · For routine maintenance on microtomes
- Cleans, lubricates and protects equipment even in the most extreme temperature environments
- Fine point application tip permits users to reach all necessary friction points

MER 6675 2oz (60mL) Bottle Each

#### **DEODEROL**<sup>™</sup>

- Neutralizes odors instantly on contact instead of just masking or covering them up and no lingering scent
- · Clear, odorless, liquid formulation

CAN DO0717 8oz Spray Bottle Each

#### SHUR/LUBE™ LUBRICANT

- Lubricating oil prevents rust and extends machine's life
- For the cleaning, lubrication and protection of cryostats and microtomes in high volume laboratories and extreme temperature environments

TRX HLUB 3/40z Lubricating Oil Each

#### PRINTMATE HOT FOIL TAPE

- Single role of tape measures  $1\frac{1}{8}$ " x 400ft
- Prints approximately 20,000 cassettes
- Suitable for SHUR/MARK® Labellers and Thermo Scientific® Printmate printers


MER PM400 Hot Foil Tape 5/Pack



31



WEB MercedesScientific.com EQUIPMENT









#### **MICROTOME LEVELING DISC**

- Less than dime-size disc placed on flat part of the microtome makes reliable bubble movement and consistent sensitivity
- Straight cut across any XY direction is assured when small black dot is centered in the bubble

TNR TN10 Leveling Disc 3/Pack



#### **PARAFFIN COLLECTION TRAYS**

- Constructed of durable acrylic
- Fits all brands of microtomes
- Slides under instrument to catch paraffin

**MKL ML9844** 

11½"L x 6¼"W x 1¼"H

Each



#### **UNIVERSAL MICROTOME ALIGNER**

- · Accommodates all makes and models of microtomes
- Fits into the cassette clamp and aligns the microtome utilizing the bubble level with the "X", "Y" and "Z" axis
- Microtome must be level and specimen head clamp must be free of any debris to receive optimal performance
- Ensures microtome is ready to recut biopsies

MRE UMA Microtome Aligner Each

# SPECIMEN COLLECTION





0.045MM (1.75MIL) THICKNESS

# TANNER SCIENTIFIC® BIOHAZARD BAGS

- Three-wall construction (6" x 9") specimen bag
- Ideal for the transportation and delivery of biohazardous specimens or materials
- Resealable specimen bag contains a document pouch with a flap to ensure minimal movement of documents during transport
- Biohazard symbol <u>ABOVE</u> the tear line (weak-line) so when opened, the user only disposes the top portion of the bag, leaving more room in hazardous waste disposal

**TNR 69BIOBAG** 

Bio Bag, 6" x 9"

1,000/Case



#### **PROSTATEPAK™**

#### **SPECIMEN COLLECTION KITS**

- Ideal for the transportation, collection and identification of multiple tissue specimens from a biopsy
- Each container contains prefilled 10% neutral buffered formalin (NBF) and packaged in a convenient transport
- Site specific labels and a location diagram to keep specimen samples organized

MER PRO08	8 Vials	Kit
MER PRO10	10 Vials	Kit
MER PRO12	12 Vials	Kit
MER PRO14	14 Vials	Kit
MER PRO16	16 Vials	Kit

# MULTI-PURPOSE CONTAINERS

WITH SNAP-ON LIDS

- Translucent to view contents
- Containers are made of high density polyethylene (HDPE)
- Ideal for the collection and storage of histology and pathology specimens, fluids, powders and solids
- Disposable and stackable



GLO 271004*	4oz (120mL)	300/Case
GLO 271008	8oz Short (240mL)	100/Case
GLO 271009	8oz Tall (240mL)	100/Case
GLO 271016	16oz (480mL)	100/Case
GLO 271032	32oz (1,000mL)	100/Case
GLO 271064	64oz (2,000mL)	50/Case
GLO 271086	86oz (2,500mL)	25/Case
GLO 271172	172oz (5,000mL)	10/Case

NOTE: \*40Z CONTAINER IS MADE OF POLYPROPYLENE (PP)



### MERCEDES SCIENTIFIC®

### 10% (NBF) PREFILLED CONTAINERS WITH SNAP CLICK LIDS

- Specialized caps "click" into place to ensure airtight seal and prevent leaking
- Standard fixative used in most laboratories produces histologic artifacts that pathologists are trained to interpret professionally
- Produces consistent histologic artifacts and results when used as a primary fixative and followed by secondary fixation by alcohol

MER 102213	20mL Fill, 15mL Prefill	96/Case
MER 10240548	40mL Fill, 20mL Prefill	150/Case
MER 10240654	60mL Fill, 30mL Prefill	75/Case
MER 10240748	90mL Fill, 45mL Prefill	75/Case
MER 10240848	120mL Fill, 60mL Prefill	36/Case

### 10% NEUTRAL BUFFERED FORMALIN (NBF) PREFILLED CONTAINERS

- Safe, cost-effective way to transport, store and retrieve pathology specimens filled with 10% neutral buffered formalin to 50% of container's capacity
- Clear polypropylene jars come with Bio-Tite closure system lids that are dual-threaded to minimize leakage
- Enhances cellular staining qualities and improves cell wall rigidity resulting in firmer cutting
- · Highly visible OSHA compliant labeling

SCG 1007	7mL (3.5mL Prefill)	400/Case
SCG 1015	20mL (10mL Prefill)	100/Case
SCG 1030	40mL (20mL Prefill)	100/Case
SCG 1180	180mL (90mL Prefill)	32/Case
SCG 1480	480mL (240mL Prefill)	9/Case
SCG 2005	960mL (480mL Prefill)	5/Case
SCG 6064	64oz (32oz Prefill)	6/Case
SCG 6128	3.8L (1.9mL Prefill)	4/Case
SCG 2040	4L (2mL Prefill)	Each
SCG 6500	7.5L (3.75mL Prefill)	Each





















#### **FABRIC LOCK BOXES**

- Includes inner mesh pouch (envelope) for gel packs
- Adjustable strap for easy courier handling during pickup
- Standard royal blue regency color with silver text and red biohazard symbol printed on one side only

THY 36500	Small, 9" x 8" x 7", No Lock	Each
THY 36501	Small, Zipper Lock, 2 Keys	Each
THY 36502	Small, Zipper Lock, 2 Keys, Strap	Each
THY 36503	Large, 11" x 9" x 7", No Lock	Each
THY 36504	Large, Zipper Lock, 2 Keys	Each
THY 36505	Large, Zipper Lock, 2 Keys, Strap	Each





#### **PLASTIC LOCK BOXES**

- Made of sturdy high impact polystyrene (HIPS)
- · Made for hanging on the door with over-the-door straps
- Viewing window for timesaving convenience and hinged door opens down for ease of use
- · Includes lock and two matching keys for added security
- Standard off-white color with reflex blue text and red biohazard symbol printed on front side only

THY 36745	13½" x 10¼" x 6¾", No Strap	Each
THY 36506	13½" x 10¼" x 6¾". Straps	Each





#### STEEL LOCK BOXES

- Made from galvanized steel for sturdy construction
- · Hinges are powder-coated to be rust-resistant
- Insulation maintains temperature for sample integrity
- Protect specimens with standard lock (includes 2 keys)
- Hanging models come with rubber surface protectors in the back corners to protect doors
- Floor models have a dimpled bottom plate to prevent moisture from being trapped

MCS KF100WC	11½" x 11½" x 9", Door Clips	Each
MCS KF200	13" x 11½" x 9", Floor	Each
MCS KF300	14" x 13" x 12", Floor	Each
MCS KF4000	15" x 24" x 12", Floor	Each
MCS KF6000	24" x 24" x 12", Floor	Each
MCS KH200	8" x 12" x 3½", Hanging	Each
MCS KH300	10" x 11" x 4", Hanging	Each
MCS KH400	9" x 12" x 6", Hanging	Each
MCS KH500	11" x 14" x 7", Hanging	Each
MCS KH550	11" x 14" x 4", Hanging	Each
MCS KH600	14" x 14" x 7", Hanging	Each
MCS KH700	16" x 19½" x 7", Hanging	Each



















# GROSSING





#### REPLACEABLE TIP SCISSORS

- Stainless steel dissection scissors
- · Remove tips when dull and replace with new ones

CAN RTS555	5" Scissor Handle, 6 Tips	Each
CAN RT0005	5" Scissor Tips	6/Pack
CAN RTC015	5" Combo (1 Handle with Set of 6 Tips)	Pack
CAN RTS666	6" Scissor Handle, 6 Tips	Each
CAN RT0006	6" Scissors Tips	Each
CAN RTC016	6" Combo (1 Handle with Set of 6 Tips)	Pack
CAN RTS777	7" Scissor Handle, 6 Tips	Each
CAN RTS777 CAN RT0007	. ,	Each Each



#### **ERGONOMIC FORCEPS**

- · Made of high quality stainless steel
- Designed for embedding or delicate grossing
- · Lightweight handles and effortless squeeze-tension design for optimal embedding or grossing
- · Perforated paddle handles for secure, slip-free grip
- Reduce hand fatigue during extended procedures

CAN 992200	5½", Serrated, Angled	Each
CAN 992202	5½", Smooth, Angled	Each
CAN 992205	5½", Serrated (1.5mm)	Each
CAN 992210	5½", Smooth (3mm)	Each
CAN 992220	7", Serrated (3mm)	Each
CAN 992230	7", Smooth, Angled	Each
CAN 992240	7", Smooth (3mm)	Each
CAN 992250	7", 1 x 2 Tooth Tip	Each
CAN 992260	7", Serrated (1.5mm)	Each
CAN 992605	7", Needle Tip	Each



#### MILTEX® SERRATED TIPPED FORCEPS

MIL 62	4" Serrated Tips	Each
MIL 64	4½" Serrated Tips	Each
MIL 66	5" Serrated Tips	Each
MIL 68	5½" Serrated Tips	Each
MIL 610	6" Serrated Tips	Each
MIL 612	7" Serrated Tips	Each
MIL 614	8" Serrated Tips	Each
MIL 616	10" Serrated Tips	Each
MIL 618	12" Serrated Tips	Each



#### MILTEX® WIDE TIPPED FORCEPS

MIL 18780	4" Straight, Serrated Tips	Each
MIL 18782	4" Half Curved, Serrated Tips	Each
MIL 18784	4" Full Curved, Serrated Tips	Each



#### **PLASTIC FORCEPS**

- · Long, sturdy plastic forcep with guide pin
- · Serrated tip and grip provide easy, consistent use
- Perfect for basic dissection or classroom applications

CAS 973901 5", Plastic, Blue 50/Pack




### TANNER SCIENTIFIC® PATHPOUCH® BIOPSY BAGS

- Manufactured from solvent resistant nylon
- Suitable for processing small biopsy and histological specimens
- The thin mesh (0.2 x 0.2mm) provides excellent fluid exchange and avoids small sample loss.
- Each convenient peel-open bag allows for placing and removing specimens rapidly.
- Easy-to-tear, not stiff

Small Biopsy Bags measure 30 x 45 mm

TNR PATHBAG	White	1,000/Box
TNR PATHBAGBLU	Blue	1,000/Box





#### **HISTO·WRAP**

- · Freer flow of liquids, but not for use with acids
- Greater retention of even the smallest specimens
- · Hard, lint free surface finish recovers tiny fragments
- · Improved diffusion of processing reagents
- Superior wet strength even at 100°C
- · Available as squares, rectangles, and filter circles

OBX 225	2" x 2", Square	500/Box
OBX 235	2" x 3", Rectangle	250/Box
OBX 245	2" x 4", Rectangle	250/Box
OBX 40R	31/8", Circle	100/Pack
OBX 55R	5½", Circle	100/Pack
OBX SP22	2½" x 2½", Square	500/Box
OBX SP253	2½" x 3", Rectangle	500/Box
OBX SP33	3" x 3", Square	500/Box




#### **BIOPSY BAGS**

- Design allows unrestricted fluid movement around tissue
- Bags peel open for quick and easy specimen removal

MOR BBS10002	30 x 50mm, Small	1,000/Box
MOR BBM20002	45 x 75mm, Medium	200/Box
MOR BBL50002	75 x 95mm, Large	500/Box

#### PLEASE NOTE: CASSETTES NOT INCLUDED







### MERCEDES SCIENTIFIC® BIOPSY FOAM PADS

- Soft, highly flexible, polyester square foam pads
- Ideal for processing biopsy specimens in histopathology using standard cassettes

MER BLKSPNG	Black	1,000/Bag
MER BLUESPNG	Blue	1,000/Bag
MER GRNSPNG	Green	1,000/Bag
MER WHTSPNG	White	1,000/Bag

#### CDI° BIOPSY PADS

- Special thin non-stick foam (1" x 11/8") allows for optimal fluid transfer with virtually no carryover
- Pressure on biopsy is greatly reduced to decrease the occurrence of artifact

#### SIMPORT® BIOPSY PADS

- · Prevents biopsies from being lost during processing
- · Polyester urethane foam achieves optimum solvent flow

**SIM M4761** Blue 1.000/Pack

#### TRI-FOLD BIOPSY PAPER

- Perforated folding filter paper for containing biopsies during processing
- Special tri-fold design folds once and then one more time on end to secure biopsy completely
- Fit edge to edge in cassette ensuring biopsy security
- Easy security at time of grossing, efficient "unwrapping" during embedding

CAN DX26280	Blue	500/Case
CAN DX26290	White	500/Case

#### **BIOPSY WRAPS**

- Optimal infiltration while preventing loss of specimen
- Colors sharply contrast to help facilitate the location of small biopsy specimens after processing
- · Specialty filter paper cuts to an easy, ready-to-use size
- · Available in dispenser packs for convenience

CAN BB0023	Blue, 2" x 3"	500/Box
CAN BB0202W	White, 2" x 2"	500/Box
IPC RL2X3PM	Pink, 2" x 3"	1,000/Box
IPC RL2X3WM	White, 2" x 3"	1,000/Box

#### **PAPER TEA BAGS**

- · Oxygen bleached paper used in biopsy processing
- · Porous paper provides maximum tissue protection

**TEA BAG** 21/4" x 21/4". White 500/Bag

#### **OPTICAL LENS PAPER**

- · Suitable for cleaning glass lenses on cameras, microscopes and other optical instruments
- · 50 lint-free sheets glued and bound per booklet

**CAS 5220** 4" x 6" (Sheet) 25 Booklets/Pack

#### GENERAL DATA® BIOPSY PAPER

- · Non-absorbent porous mesh prevents reagent carryover
- · Efficient solvent flow-through during tissue processing
- · Reduces impressions made on tissue specimens
- · COLORS AVAILABLE: Black or White

TRX BP-W 21/2" x 1" (11/4" x 1" Folded) 500/Pack

#### **CELLSAFE® BIOPSY CAPSULES**

- Two-hinged interlocking frames with integral nylon mesh
- Extra fine mesh ensures high level of specimen security, zero tissue damage and low reagent carryover
- · COLORS AVAILABLE: White or Blue

**CEL 1110 Biopsy Capsules** 1,000/Bag














- · Cut-resistant, foam-core board with light blue color and measurement markings provides great color contrast
- · For routine grossing, frozen sections and photography
- · Can be used multiple times with different quadrants
- · Marked with inch and metric rulers and pre-measured circles for quick size referencing
- · Incinerate or dispose after use

**CAN BB081160** 8" x 10", Light Blue 60/Box



#### **REUSABLE CUTTING BOARDS**

- · Self seal cutting base
- · Permanently sealed to eliminate cracking or warping
- · 5-ply construction for increased longevity
- · Prevents blade run and damage to cutting edges
- · Minimal cut marks due to rejuvenating material qualities
- · Durable PVC surface allows for easy cleaning

MKI 5800W	12" x 18", White	Each
MKI 5800G	12" x 18", Gray	Each
MKI 5801W	18" x 24", White	Each
MKI 5801G	18" x 24", Gray	Each



#### **DISSECTABLE™**

- Made of heavy-duty, stain-resistant, thick polyethylene
- · Drain groove carved all around the edge contains fluids
- Large cutting area with multiple dimensional scales in inches, centimeters, millimeters and circles
- Other side features two smaller cutting boards also supplied with dimensional scales
- Heavy-duty base elevates DissecTable<sup>™</sup> to the right height and will also retain excess fluid if necessary

SIM M618	DissecTable™ Jr., 13" x 11" x 5"	Each
SIM M620	DissecTable™, 23" x 16" x 0.5"	Each
SIM M625	DissecTable™ Board Base	Each



#### **DISPOSABLE CUTTING BOARDS**

- Foam inner liner for pinning tissue
- · Standard and metric scale on cutting surface
- · Reference areas for assessing specimen size
- · Ideal for photographic purposes
- Alleviates clean-up worries when working with infected tissues by disposal and incineration
- · Excellent for teaching institutions

MKI 5700	Small, 8"W x 6²/₃"H	48/Pack
MKI 5701	Medium, 16"W x 6"H	24/Pack
MKI 5702	Large, 16"W x 13"H	12/Pack



#### **CORK DISSECTING BOARDS**

- · Specifically designed for pinning tissue
- Neoprene cork cemented to marine plywood base
- Provides a very durable cutting surface resistant to absorption while permitting thorough rinsing
- Removable, non-skid legs provide the option to raise the cutting board for a more ergonomic workstation

MOP BC004	23"L x 16"W x 1"H	Each
MOP BC005	16"L x 12"W x 1"H	Each
MOP BC023	23"L x 16"W x 4"H, Legs	Each
MOP BC024	16"L x 12"W x 4"H, Legs	Each



#### **TEASING NEEDLES**

- Solid white plastic handle eliminates absorption of contaminants such as mold, fungus or fixatives
- · Stainless steel needle eliminates rust and corrosion

CAS 971313	Straight	12/Pack
CAS 971314	Bent	12/Pack



#### PLASTIC SCALPEL HANDLE

- · Ideal for post mortem applications
- · Reusable, economical and autoclavable
- · Lightweight plastic scalpel handle with a comfort grip
- Fits autopsy blades #60XT, #70XT and #AUT22XT
- · Refined, lightweight and durable

HAV HAN8P #8, Plastic, Comfort Grip Each



#### STAINLESS STEEL BLADE HANDLES

- · Ideal for post mortem applications
- · Hand-crafted of premium quality stainless steel
- · Refined, highly polished and extremely durable
- Fit various blades sizes for ergonomic user preference

HAV HAN74	Fits #60XT, #70XT	Each
HAV HAN8	Fits #60XT, #70XT, #22XT	Each
HAV HANG3M	Fits #10 to #15C	Each
HAV HANG74	Fits #20 to #25	Each



• Pins (2" Length, %16" Head) are perfect for dissections

EIS PHTPINS Nickel Plated T-Pins 100/Pack







#### BARD PARKER® SCALPEL BLADES

- Intended for tissue separation or procedures requiring sharp surgical blades to puncture or cut
- Superior cutting advantage of stainless steel technology
- · Consistent sharpness, control and strength
- Rib-Back® blade design for greater strength and control
- Individually packaged sterile for safety and ease of use

ASP 371610	#10, Disposable, Sterile	10/Box
ASP 371615	#15, Disposable, Sterile	10/Box
ASP 371621	#21, Disposable, Sterile	10/Box
ASP 371310	#10, Rib-Back <sup>®</sup> , Non-Sterile	150/Case
ASP 371311	#11, Rib-Back®, Non-Sterile	150/Case
ASP 371315	#15, Rib-Back®, Non-Sterile	150/Case
ASP 371321	#21, Rib-Back®, Non-Sterile	150/Case
ASP 371322	#22, Rib-Back®, Non-Sterile	150/Case
ASP 371340	#60, Rib-Back®, Non-Sterile	150/Case
ASP 371150	#10, SafetyLock™ Rib-Back®	50/Box
ASP 371151	#11, SafetyLock™ Rib-Back®	50/Box
ASP 371153	#15, SafetyLock™ Rib-Back®	50/Box
ASP 371155	#21, SafetyLock™ Rib-Back®	150/Case
ASP 371156	#22, SafetyLock™ Rib-Back®	150/Case



#### SUPER SCALPEL PM BLADE

- REDESIGNED NEW STYLE super scalpel with integrated black plastic handle that is also ergonomical
- "Thumb-notch" for leverage and comfort
- No need to replace blade: dispose whole knife when dull
- · Designed for autopsies and large specimens

**CAN SS0120** #120 (PM40) Blade 20/Case



#### **FEATHER SURGICAL BLADES**

- Stainless steel, disposable blades are characterized by superior cutting quality, strength and corrosion resistant
- Individually wrapped and sealed in aluminum foil to maintain sterility and ensure long shelf-life

**ANX 2976** #22, Sterile 100/Box



#### **BLADE-SAFE® BLADE REMOVAL**

- · Stainless steel construction for durability and stability
- Safely removes blades from handles without hands and remove the risk of puncturing skin

MIL 4100 Blade-Safe® Blade Removal Each



#### **GLASSVAN® SCALPEL BLADES**

- High-quality raw material construction creates longer lasting cutting edges and superior sharpness
- · Individually wrapped sterile

MYC 2001T10	#10, Carbon Steel	100/Box
MYC 2001T11	#11, Carbon Steel	100/Box
MYC 2001T15	#15, Carbon Steel	100/Box
MYC 3001T15C	#15C, Stainless Steel	100/Box
MYC 2001T21	#21, Carbon Steel	100/Box
HAV FSC22	#22, Carbon Steel	100/Box
MYC 2001T25	#25, Carbon Steel	100/Box
HAV AUT60	#60, Carbon Steel	100/Box
MYC 2001T70GV	#70, Carbon Steel	100/Box



- Safe way to remove and contain used scalpel blades with single-hand technique: push and click
- · Minimizes risks of injury and infection
- · Built-in counter automatically shuts off to stop overfilling
- · Puncture proof sharps container

**QST QFRGEN** 100 Blade Capacity Each



#### **STEALTH™ SAFETY SCALPELS**

- Integrated rigid plastic handle with protective sheath
- Slide sheath down for use, slide up when not in use
- · Dispose when blade dulls

CAN SL0022	#22	10/Box
CAN SL0060	#60	10/Box
CAN SL0070	#70	10/Box



#### **FINGER SCALPELS**

- · Blade attached to small slip resistant handle
- Ultimate in dexterity and convenience

<b>CAN FS2222</b>	#10, Non-Sterile	100/Case
CAN FS0022	#22, Non-Sterile	100/Case
CAN FS0060	#60. Non-Sterile	100/Case



#### TECHNOCUT® SURGICAL SCALPELS

- · Superior sharpness and uniform cutting edges
- Each stainless steel blade covered with plastic protector
- · Plastic handles are marked with a millimeter scale
- · Individually wrapped and hermetically sealed

MYC 6008T10	#10	10/Box
MYC 6008T11	#11	10/Box
MYC 6008T12	#12	10/Box
MYC 6008T15	#15	10/Box
MYC 6008T21	#21	10/Box
MYC 6008T22	#22	10/Box



#### RETRACTABLE SCALPELS

- Top activated retractable feature exposes stainless steel blade and minimizes chance of injury during use
- Audible clicks indicate blade position
- Non-slip ribbed handles with metric scale
- Blade size indicated by color of slider
- · Individually packaged sterile

MYC 6008TR10	#10	10/Box
MYC 6008TR11	#11	10/Box
MYC 6008TR12	#12	10/Box
MYC 6008TR15	#15	10/Box



AMX 740001





AMX 740002

#### PERSONNA® BLADES

Single-Edge Prep Razor: Comfort-Coated™ for supersmooth cutting and specimen preparation

DermaBlade Shave Biopsy: Sterile, MicroCoat® coated stainless steel blade reduces skin and tissue trauma

Double-Edge Pathology: For use by pathologists in grossing tissue specimens prior to processing

Double-Edge Prep: Carefully ground and honed; added sharpness, edge retention and corrosion resistance

AMX 720001	DermaBlade Shave Biopsy	50/Box
AMX 720200	Double-Edge Pathology	100/Box
AMX 740002	Double-Edge Prep	100/Box
ANX D2877	Pathco Razor Blade Handle	Each





#### **RAZOR BLADES**

- · Single-edged, surgical, carbon-steel blades
- · Cutting edge on blade protected by paper shell
- · For industrial, general and scientific purposes

FIS S17302	0.009" Thick	100/Pack
FIS S65921	0.012" Thick	100/Pack





#### **MACRO PREP KNIVES**

- · Sheffield stainless steel with longer lasting honed edges
- Disposable blade with a molded, integrated, non-slip plastic handle (finger guard) for autopsy and grossing

MKI 5004	4" Cartilage Knife	10/Box
MKI 5006	6" Trimming Knife	10/Box
MKI 5000	8" Autopsy/Surgical Knife	10/Box
MKI 5010	10" Lung Knife	10/Box
MKI 5012	12" Brain Knife	10/Box



- · All-in-one blade and integrated handle
- Extremely thin blade for trimming and sectioning tissue samples during surgical, grossing or autopsy procedures
- · Dispose entirely when dull

**CAN ULT702** 6". Stainless Steel 5/Box



#### **DISSECTION KNIFE**

· Extra sharp, hollow ground high-carbon stainless steel blade cast permanently to an all aluminum handle

**CAN GK0007** 7", Stainless Steel 5/Box

## **MARKING**



#### TISSUE MARKING DYE

- Comes in 2oz. and 8oz. bottles
- 7 Colors Available: Red, Orange, Yellow, Green, Blue, Violet, Black
- · Special formulation to increase viscosity
- Ideal for use on formalin fixed, fresh or frozen tissues
- Pigments immediately bond to tissue surfaces without any drying or waiting time, and are preserved throughout tissue processing and staining
- Buy individually or in sets

Ask us about our set options - Accessories available!



#### **MERCEDES SCIENTIFIC®**

#### **20Z TISSUE MARKING DYE**

- · Special formulation and increased viscosity
- · Ideal for use on formalin fixed, fresh or frozen tissues
- Pigments immediately bond to tissue surfaces without any drying or waiting time, and are preserved throughout tissue processing and staining

MER DYE2OZBK Black Each MER DYE2OZBL Blue Each MER DYE2OZG Green Each MER DYE2OZO Orange Each MER DYE2OZR Red Each MER DYE2OZV Violet Each MER DYE2OZY Yellow Each	MER DYE2OZA	1 of Each Color	Set
MER DYE2OZG         Green         Each           MER DYE2OZO         Orange         Each           MER DYE2OZR         Red         Each           MER DYE2OZV         Violet         Each	MER DYE2OZBK	Black	Each
MER DYE2OZO Orange Each MER DYE2OZR Red Each MER DYE2OZV Violet Each	MER DYE2OZBL	Blue	Each
MER DYE2OZV Red Each MER DYE2OZV Violet Each	MER DYE2OZG	Green	Each
MER DYE2OZV Violet Each	MER DYE2OZO	Orange	Each
	MER DYE2OZR	Red	Each
MER DYE2OZY Yellow Each	MER DYE2OZV	Violet	Each
	MER DYE2OZY	Yellow	Each

### MERCEDES SCIENTIFIC° TISSUE MARKING DYE

• Sets Include: Plastic Black Tray, 50 Applicator Sticks

MER DYESET5	Includes 2oz and 8oz Bottles of Each: Blue, Black, Green, Red, Yellow	Set
MER DYESET6	Includes 2oz and 8oz Bottles of Each: Blue, Black, Green, Orange, Red, Yellow	Set
MER DYESET20Z6	Includes 2oz Bottle of Each: Blue, Black, Green, Orange, Red, Yellow	Set
MER DYESET80Z6	Includes 8oz Bottle of Each: Blue, Black, Green, Orange, Red, Yellow	Set
MER TRAY5	Tray Holds (5) 2oz Bottles and (5) 8oz Bottles	Each
MER TRAY6	<b>Tray Holds</b> (6) 2oz Bottles and (6) 8oz Bottles	Each
MER TRAY2OZ6	Tray Holds (6) 2oz Bottles	Each





### MERCEDES SCIENTIFIC® 80Z TISSUE MARKING DYE

• Same features as 2oz bottles (MER DYE2OZ)

MER DYE8OZBK	Black	Each
MER DYE8OZBL	Blue	Each
MER DYE8OZG	Green	Each
MER DYE8OZO	Orange	Each
MER DYE8OZR	Red	Each
MER DYE8OZV	Violet	Each
MER DYE8OZY	Yellow	Each



#### CDI<sup>®</sup> INK AID

An Acetic Acid mordant solution designed to be used with CDI's Tissue Dye formula. CDI's Tissue Marking Dyes are formulated to be used without a mordant, however Ink Aid may aid in adherence to particularly fatty tissue.

- Composition: >93% Water, <6% Acetic Acid, Propylene Glycol
- · May aid in adherence to particularly fatty tissue
- Acetic acid solution for use with CDI\* tissue dye formula

CAN IA4000 8oz Spray Bottle 5/Box



#### **IPC BLUE™**

#### **TOLUIDINE BLUE BIOPSY MARKING DYE**

- · Ideal for use at the grossing station
- Formula made in 10% neutral buffered formalin
- Dark blue color allows for superior visualization of biopsy during embedding and cutting
- · Convenient dispenser bottle prevents spillage or dry out
- · Can be added directly to the tissue or on the processor
- · No interference with routine H&E, IHC or special staining
- · Continues fixation while grossing
- Dye is completely replaced when the slide is stained
- Formalin-based dye won't evaporate as rapidly as similar alcohol-based alternatives

IPC RMRIB-01	1oz (30mL)	Each
IPC RMRIB-02	2oz (60mL)	Each
IPC RMRIB-16	16oz (180mL)	Each





- pens, lettering pens and airbrushesWill not chip, peel, feather or smear
- Can be diluted with water for transparent washes and to produce infinite shades
- · For use on paper, board and film

HIG 44201 1oz Bottle





#### **APPLICATOR STICKS**

- · Ideal for applying or stirring medicines
- · Select-grade wood provides a smooth finish
- · Made from close-grained white birch

DUK 9000	6" Wooden	1,000/Box
PRA 76200	6" Cotton Tipped	1,000/Box

#### POINTED WOOD APPLICATORS

- · Each wooden stick is pointed on one end
- · Made from non-sterile birch wood for multi-purpose use

CAN PT112 Pointed Stick, One End 100/Pack











#### **CDI® TISSUE MARKING DYES**

- · For marking margins and multiple biopsies
- · Define margins without bleeding, changing color or fading
- · Optimal viscosity for easier application
- Bottles sealed with flip top dispenser cap
- 10 Colors Available: Black, Blue, Green, Red, Yellow, Orange, Violet, Gold, Teal, Magenta

CAN 0728	8oz Bottle	Each
CAN 0727	2oz Bottle	Each

### CDI<sup>®</sup> TISSUE MARKING DYE MULTI-USE SQUEEZERS KITS

- 4mL prefilled dispenser "squeezers" with caps for recapping unused portion
- · Dye holds permanently through processing and staining
- · Air-tight packaging keeps dyes from thickening
- 7 Colors Available: Black, Blue, Green, Red, Yellow, Orange, Violet

CAN MD100	4mL Prefilled Squeezer 15/Bo	
CAN MD2000 6-Color Kit (3 of Each Color, Holding Tray)		Each
CAN MD3000	7-Color Kit (3 of Each Color, Holding Tray)	Each





CAN 04000

CAN 04300

#### CDI® EASY-SQUEEZE, DROPPER-TIPS

- Squeeze bottles easily dispense for large applications
- 10 Colors Available: Black, Blue, Green, Red, Yellow, Orange, Violet, Gold, Teal, Magenta

CAN 0725	2oz Easy-Squeeze Bottle	Each
CAN 51117	2oz. 7-Color Set. No Trav	Set

#### CDI® BRUSH TIP APPLICATORS

- For small biopsies where precise application is required
- · Available in 0.10 fl. oz (3mL) brush tip applicators
- 7-Color: Black, Blue, Green, Red, Yellow, Orange, Violet
- 3-Color: Gold, Teal, Magenta

CAN 0726 Brush Tip Applicator		5/Pack
CAN 04000	7-Color Set, Plastic Tray	Each
CAN 04300	3-Color Set, Plastic Tray	Each

#### CDI® TISSUE MARKING DYE KITS

• All Kits Include 50 Applicator Sticks

CAN 01000	5-Dye Kit, Wood Tray	Each
CAN 02000	6-Dye Kit, Wood Tray	Each
CAN 03000	7-Dye Kit, Wood Tray	Each
CAN 01000P	5-Dye Kit, Plastic Tray	Each
CAN 02000P	6-Dye Kit, Plastic Tray	Each
CAN 03000P	7-Dye Kit, Plastic Tray	Each





### MERCEDES SCIENTIFIC° BLACK LAB MARKERS

- Designed for use on lab equipment and consumable surfaces including plastic, glass, porcelain and metal
- Fine tip allows for detailed transcriptions on microscope slides and cassettes in permanent ink
- Practical clip prevents rolling away and can be secured on clothes and objects
- Cap can be stored on the end of the barrel to prevent risk of losing it
- · Special ink prevents bleeding through the paper

MER MARKER Black Lab Markers 12/Box



#### **EXTRA FINE PERMANENT MARKERS**

- Xylene-free formula marks permanently on most surfaces including glass, metal, wood and plastic
- Narrow-barrell marker features fast, nearly instantaneous dry time to help minimize smudging and smearing

VAR PIL44102	Black	12/Box
VAR PIL44103	Blue	12/Box
VAR PIL44104	Red	12/Box
VAR PIL44105	Green	12/Box



#### **TUNGSTEN-POINT LAB PEN**

- · Permanently mark glass slides, test tubes and labware
- Sharp, tungsten carbide tip etches glass, ceramic and plastic, and is retractable for safe storage

HEA HD1846 Tungsten-Point Lab Pen Each



#### KP MARKER PLUS®

- Laboratory marking pen is excellent for use on microscope slides and all types of cassettes
- · Can be used as a coverslip marking pen
- Ultra fine point is able to write clearly on such small items as the PCR tubes or microcentrifuge tubes
- · Quick drying time and virtually no smudging
- Withstands extremely low temperatures without smearing or fading
- · Chemically resistant to most solvents and stains
- · Ink is waterproof, smearproof and odorless

**KLI KPMARKER** KP Marker Plus®, Black 12/Box





#### **SUPER PAP PENS**

- Surface tension provided by the circle drawn with the pen ensures that only the amount of antibody solution needed for sufficient reaction will be applied
- Stable up to 129°C in microwave heat

CAN SPM0928	Mini Tip	Each
CAN SPR0905	Regular Tip	Each

#### **TISSUE CAPTURE PEN**

- For better adherence of tissue sections to slides
- Pen "charges" approximately 3,500-5,000 slides

CAN TCO417 Tissue Capture Pen Each



#### **BLACK MARKING PENCILS**

- · Smooth dark writing on slides and cassettes
- · Formulated for laboratory use

CAN N1000 Black Marking Pencil 12/Box





#### ENDURE+™ & ENDURE™

#### **MICROTOME BLADES**

ENDURE+	Profile Type	Blade Length	Blade Thickness	Blade Height	иом
MER LP100	Low Profile	76mm	0.25mm	8mm	50/Pack
MER HP200	High Profile	76mm	0.30mm	14mm	50/Pack
ENDURE					
TNR 315LP	Low Profile	80mm	0.25mm	8mm	50/Pack
TNR 554HP	High Profile	80mm	0.43mm	14mm	50/Pack

#### **MERCEDES SCIENTIFIC®**

#### MICROTOMEBLADE

	Profile Type	Blade Length	Blade Thickness	Blade Height	иом
MER MLP300	Low Profile	76mm	0.25mm	8mm	50/Pack
MER MHP500	High Profile	76mm	0.43mm	14mm	50/Pack

#### **MERCEDES SCIENTIFIC®**

#### **MICROTOMEBLADEPLUS**

	Profile Type	Blade Length	Blade Thickness	Blade Height	иом
MER 300PLUS	Low Profile	76mm	0.25mm	8mm	50/Pack
MER 500PLUS	High Profile	76mm	0.43mm	14mm	50/Pack

#### **CL STURKEY® DIAMOND**

CLS DT315DD	Low Profile, Bulk Dispenser	200/Pack
CLS DT315D50	Low Profile, Standard Dispenser	50/Pack

#### **CL STURKEY® EXTREMUS**

CLS DT315X50	Low Profile, Standard Dispenser	50/Pack
CLS DT315X	Low Profile, Box	170/Pack
CLS DT315XD	Low Profile, Bulk Dispenser	200/Pack
CLS D554X50	High Profile, Standard Dispenser	50/Pack
CLS D554X	High Profile, Box	100/Pack
CLS D554XD	High Profile, Bulk Dispenser	120/Pack

#### **CL STURKEY® GOLD**

CLS DT315G50	Low Profile, Standard Dispenser	50/Pack
CLS DT315G	Low Profile, Box	170/Pack
CLS DT315GD	Low Profile, Bulk Dispenser	200/Pack
CLS D554G50	High Profile, Standard Dispenser	50/Pack
CLS D554G	High Profile, Box	100/Pack
CLS D554GD	High Profile Bulk Dispenser	120/Pack

#### **CL STURKEY® INFINITY**

CLS DT315N50	Low Profile, Standard Dispenser	50/Pack	
CLS DT554N50	High Profile Standard Dispenser	50/Pack	



































#### **CL STURKEY® MONARCH**

**CLS DT315M50** Low Profile, Standard Dispenser 50/Pack

#### **CL STURKEY® PREMIUM**

CLS D554P50	High Profile, Standard Dispenser	50/Pack
CLS D554P	High Profile, Box	100/Pack
CLS D554PD	High Profile, Bulk Dispenser	120/Pack

#### **CL STURKEY® SELECT**

CLS DT554S50	High Profile, Standard Dispenser	50/Pack	
CLS DT315S50	Low Profile, Standard Dispenser	50/Pack	

#### **CL STURKEY® SILVER**

CLS DT315R50	Low Profile, Standard Dispenser	50/Pack
CLS DT315R	Low Profile, Box	170/Pack
CLS DT315RD	Low Profile, Bulk Dispenser	200/Pack
CLS D554R50	High Profile, Standard Dispenser	50/Pack
CLS D554RD	High Profile, Bulk Dispenser	120/Pack

#### THERMO SCIENTIFIC® DISPOSABLE MICROTOME BLADES

ITEM #	DESCRIPTION	PROFILE	ANGLE	QUANTITY
FIS 3050835	MB35 Premier	Low	35°	50/Pack
FIS 3051835	MX35 Premier	Low	34°	50/Pack
FIS 3052835	MX35 Premier+	Low	34°	50/Pack
FIS 3053835	MX35 Ultra	Low	34°	50/Pack
FIS 3050822	MB22 Premier	Low	22°	50/Pack
FIS 3150734	HP35 Coated	High	35°	50/Pack
FIS 3153735	HP35N Ultra	High	34°	50/Pack

#### **ENCORE®**

CRE 7380	Low Profile	50/Pack
CRE 7390	High Profile	50/Pack

#### **DURAEDGE®**

CRE 7223	Low Profile	50/Pack
CRE 7310	High Profile	50/Pack

#### **EDGE-RITE**™

FIS 1407060	Low Profile	50/Pack
FIS 1407061	High Profile	50/Pack

#### SHUR/SHARP®

TRX DMBHD	Heavy Duty	35/Pack
TRX DMBHP	High Profile	50/Pack
TRX DMBLP	Low Profile	50/Pack



### TANNER SCIENTIFIC® METAL BASE MOLDS

- · High-quality stainless steel construction
- · Rounded corners for optimal paraffin ribboning
- · Superior thermal exchange and conductivity
- · Highly polished surface for easy specimen removal
- · Fits most embedding rings and cassettes

TNR M07075	7 x 7 x 5mm	12/Pack
TNR M15155	15 x 15 x 5mm	12/Pack
TNR M24245	24 x 24 x 5mm	12/Pack
TNR M30245	30 x 24 x 5mm	12/Pack
TNR M37245	37 x 24 x 5mm	12/Pack



### MERCEDES SCIENTIFIC° MOLD RELEASE AGENT

- Stops paraffin and resin build-up from sticking and easily releases cooled paraffin blocks
- Cleans, conditions and extends a mold's life
- · Packaged in a non-aerosol pump spray bottle

MER 6015125ML

4oz Spray Bottle





#### **DISPOSABLE BASE MOLDS**

- Specially formulated plastic material (PVC) offers excellent thermal exchange for paraffin
- Smooth interior finish and rounded corners facilitate specimen removal

SIM M4751	7 x 7 x 5mm	1,000/Case
SIM M4752	15 x 15 x 5mm	1,000/Case
SIM M4753	24 x 24 x 5mm	1,000/Case
SIM M4754	30 x 24 x 5mm	1,000/Case
SIM M4755	37 x 24 x 5mm	1,000/Case
SIM M47510	37 x 24 x 10mm	500/Case



#### **PARAFFIN BLOCK MAILERS**

- Made of sturdy PVC and transparent to view contents
- · Specifically designed for transporting paraffin blocks
- · Attached cover easily closes but still very secure
- Suitable for all regular models of tissue/biopsy cassettes
- Use for handling and shipping up to 6 blocks as needed while keeping only 1 model of block mailer in inventory

SIM M4776 Holds 6 Blocks 50/Bag

PIN 121		18mm	4 Wells	Each
PIN 122		24mm	4 Wells	Each
PIN 123		30mm	3 Wells	Each
PIN 125	0000	24mm, 30mm	2 Wells Each	Each
PIN 126		50 x 30mm Well Bar	Section Completely Fills a Standard Slide	Each
PIN 127	0000	20 x 20 Deep Well Bar	(2) 6mm Deep 20mm Wells, (2) 9mm Deep 20mm Wells	Each
PIN 128	PT PT T X 23	24 x 24 Deep Well Bar	(2) 6mm Deep 24mm Wells, (2) 9mm Deep 24mm Wells	Each
PIN 129	(0000)	10mm Well Bar	(4) 10mm Square Wells, (2) 4mm Deep Wells, (2) 6mm Deep Wells	Each
PIN 111	<b>₽ 8888 ₽</b>	Small Embedding Well Bar Set, 18mm	(4) Small Chucks, (2) Over-Chuck Freezing Blocks, (8) Dispensing Slides	Set
PIN 112	P 6805 D	Medium Embedding Well Bar Set, 24mm	(4) Small Chucks, (2) Over-Chuck Freezing Blocks, (8) Dispensing Slides	Set
PIN 113	•	Large Embedding Well Bar Set, 30mm	(4) Large Chucks, (2) Over-Chuck Freezing Blocks, (8) Dispensing Slides	Set
PIN 114	-	24/30 Well Bar Set, 24mm/30mm Wells	<ul><li>(2) 24mm Wells, (2) 30mm Wells,</li><li>(2) Large Chucks, (2) Small Chucks,</li><li>(2) Over-Chuck Freezing Blocks,</li><li>(8) Dispensing Slides</li></ul>	Set
PIN 115		50 x 30mm Well Bar Set	(1) 50 x 30mm Well Bar, (2) 55 x 35mm Chucks	Set
PIN 116	*******	20 x 20 Deep Well Bar Set	(1) 20 x 20 Deep Well Bar, (2) Small Chucks, (2) Over-Chuck Freezing Blocks, (8) Dispensing Slides	Set
PIN 119		Freezing Griddle, Elevated Freezing Block	Epoxy-Coated Forceps, Putty Knife, Extra Elevating Bars	Set
PIN 101		Complete Set	(3) Well Bars (18, 24 & 30mm), (6) Small Chucks, (4) Large Chucks, (4) Over-Chuck Freezing Blocks, Chuck Bin, Cutting Board/Freezing Griddle, (16) Dispensing Slides, Elevated Freezing Block, Angled Embedding Forceps and Accessories	Set

PIN 131		Small Chuck, 28mm Square, Stainless Steel		Each
PIN 132		Large Chuck, 36mm Square, Stainless Steel		Each
PIN 134		X-Large Chuck, 55 x 26mm Rectangle, Stainless Steel		Each
PIN 141		Over-Chuck Freezing Block		Each
PIN 142		Chuck Bin, Stainless Steel	Loose Fitting Unattached Cover, Easy Access	Each
PIN 143	MUM	(8) Dispensing Slides		Each
PIN 171	•	Chuck Adapter	For Shandon™	Each
PIN 172	3	Chuck Adapter	For Tissue-Tek®	Each
PIN 181		Frozen Section Brushes, Flat White Boar Bristles, Prevents Curling Through Brush	(2) #1 (3/8" Angle Cut), (2) #2 (1/2" Angle Cut), Cut at Angle, Very Stiff	4/Pack





#### **EPREDIA™ HISTOPLAST**

- · Formulation comprised of paraffin and plastic polymers
- · Excellent processing yields thin, wrinkle-free sections

FIS 8330	Histoplast PE, 2lb	54° to 57°C	8/Case
FIS 8331	Histoplast Plus, 2lb	56°C	8/Case
FIS 8332	Histoplast LP. 2lb	50° to 54°C	8/Case



#### **PUREAFFIN®**

- Ultra-purified paraffin wax with additives to enhance infiltration, impart elasticity and provide easier cutting
- For open and closed tissue processors
- Blended to provide thorough tissue impregnation
- Facilitates excellent sectioning attributes

PureAffin\* I™: Formulated with DMSO for enhanced infiltration and elasticity, ideal for infiltrating fatty tissues

**PureAffin**\*  $E^{\infty}$ : Use as a two-step paraffin system with *PureAffin*\*  $i^{\infty}$  or one-step for infiltration and embedding

PureAffin® Blue™: Tinted blue provides an optimal color contrast to clear or white biopsies; for superior cutting

PureAffin® X™: Harder paraffin with a lower melting point

ITEM#	DESCRIPTION	MELTING POINT	иом
CAN EEPAR1	PureAffin® I	~56°C	8/Case
CAN EEPAR2	PureAffin® E	~56°C	8/Case
CAN EEPAR2B	PureAffin® Blue	~56°C	3/Case
CAN EEPAR3	PureAffin® X	50-53°C	3/Case



#### **EPREDIA™ PARAFFIN TYPE 1, 3, 6, 9**

- Flake form has a greater surface-to-weight ratio than pellets resulting in a faster melting time
- Type 1: contains no polymer, offers rapid infiltration
- Type 3: contains minimal polymer, offers fast infiltration
- Type 6: harder than Type 3, allows thinner sections
- Type 9: highest polymer content, extremely thin sections

ITEM#	DESCRIPTION	MELTING POINT	иом
FIS 22050124	Type 1, 2lb	55° and 57°C	10/Case
FIS 22050125	Type 3, 2lb	55° and 57°C	10/Case
FIS 22050126	Type 6, 2lb	55° and 57°C	10/Case
FIS 22050127	Type 9, 2lb	55° and 57°C	10/Case



#### **PARAPLAST®**

- Regular: paraffin and plastic polymer construction; minimizes wrinkles for clear, smooth tissue specimens; compatible with most solvents and staining techniques
- Plus: added dimethyl sulphoxide (DMSO) improves rate of tissue infiltration, easier for cutting difficult tissues
- Xtra: unique blend of low molecular weight polymers and paraffin, low viscosity permits complete infiltration

MCC 501006	Regular, 1kg	56°C	8/Case
MCC 502004	Plus, 1kg	56°C	8/Case
MCC 503002	Xtra, 1kg	50° to 54°C	8/Case





#### **PARAFFIN SCRAPERS**

- Floor Scraper: 491/2"L handle, durable 4"W blade
- Maxi Scraper: 4"W blade, comfort-grip handle, cover
- Safety Scraper: retractable blade, no-slip handle
- Mini Scraper: gets into tight spots around microtomes

MKL ML1138	Floor Scraper	Each
MKL ML1139	Blades (ML1138/ML7050)	10/Pack
MKL ML7050	Maxi Scraper	Each
MKL ML7049	Safety Scraper	2/Pack
MKL ML7048	Mini Scraper	4/Pack
MKL ML7051	Blades (ML7048/ML7049)	100/Pack





#### **PARAFILM®**

- · Waterproof film prevents both sample evaporation and contamination while protecting anhydrous materials
- Can be folded repeatedly and twisted around sharp surface without rupturing
- More flexible and softer as temperature increases
- Partially returns to original length when tension removed
- · Each roll is interlined with paper to prevent adhering
- · Ideal for use with culture media

FIS 1337416	2" x 250'	Roll
FIS 1337410	4" x 125'	Roll
FIS 1337412	4" × 250'	Roll



#### **PARAFFIN BRUSH**

- · Goat hair brush removes paraffin from microtomes and embedding compound from cryostats for better sections
- Features: Copper Ferrule, Dark Blue Wood Handle, Dark Bristles and EZ Blade Removal (Magnetic Tip)

CAN SLK1000 9" Brush Each

#### **ROUND CAMEL HAIR LACQUERING BRUSHES**

- · Large nominal sizes
- Each brush has an overall length of 8"
- · Hair set in solvent resistant epoxy
- · Maroon polished handles
- · Fine black camel hair with rustproof seamless aluminum ferrule



SHB 3163SZ1	#1, <sup>3</sup> ⁄₄" Brush	Each
SHB 3163SZ2	#2, <sup>13</sup> / <sub>16</sub> " Brush	Each
SHB 3163SZ3	#3, 15/16" Brush	Each
SHB 3163SZ4	#4, 11/16" Brush	Each
SHB 3163SZ5	#5, 11/8" Brush	Each
SHB 3163SZ6	#6, <sup>13</sup> /16" Brush	Each



#### **CAMEL HAIR BRUSHES**

- · Helpful for dusting sensitive instruments
- · Hair mounted on wooden handle with metal ferrule

CAS LB36	3⁄4" Brush: 61⁄2"L	Each
CAS LB37	1" Brush: 71⁄4"L	Each
CAS LB38	2" Brush: 81⁄4"L	Each
CAS LB29	Tips for CAS LB31	12/Pack
CAS LB30	Thick Brush: 9"L	Each
CAS LB31	Thin Brush: 7³/₄"L	12/Pack



#### NEG-50™ FREEZING MEDIA

- Water-soluble frozen medium yields high-quality sections down to -50°C without chatter
- · Leaves no residue after staining
- Buffered to ensure proper alkalinity (pH = 8.0)
- · Formulated to prevent tissue separation from the slide
- Available in **FIVE** different colors for easy identification
- Compatible with the Microm Automatic Cryo Approach (ACA) system

FIS 22110617	Colorless, 4oz Bottle	2/Case
FIS 22046512	Orange, 4oz Bottle	2/Case
FIS 22050451	Blue, 4oz Bottle	2/Case
FIS 22050452	Green, 4oz Bottle	2/Case
FIS 22050453	Pink, 4oz Bottle	2/Case
FIS 22050454	Yellow, 4oz Bottle	2/Case
FIS 22050450	1 of Each 4oz Bottle	6/Case

#### HISTOGEL® (SPECIMEN PROCESSING KIT)

- Encapsulate and retain even small, friable and viscous specimens during histological processing
- · Penetrates without allowing specimen to escape
- · Eliminates unwanted discoloration around specimens
- Remains inside the paraffin block during cutting and virtually unnoticeable during sectioning
- Kit Includes: (12) 10mL Vials, (2) WonderBlock Cooling Blocks, (1)
   Dry Bath Incubator, (500) Disposable Pipets

FIS 22110678 HistoGel® Kit 1/Case







#### **TFM™ TISSUE FREEZING MEDIA**

- · Advanced formulation supports tissue during cryotomy
- Enables precise cutting with support for cryosectioning at 3 microns
- Freezes tissue quickly without promoting crystallization even at the most extreme temperatures
- Easily remove during fixation and staining procedures without leaving behind residue
- Less curling to allow easier pick-up from the knife edge and flat sections to flow freely under the anti-roll device
- Reduced water content minimizes freeze fracture and does not cause autofluorescence
- Standard frozen sections for clinical diagnosis at 3u-6u are easy to obtain
- Available in <u>FIVE</u> colors enables histotechnologists to easily differentiate between multiple sections

TRX TFMC	Clear, 4oz	Each
TRX TFMG	Green, 4oz	Each
TRX TFMR	Red, 4oz	Each
TRX TFMB	Blue, 4oz	Each
TRX TFMY	Yellow, 4oz	Each

#### ULTRAFREEZE® FROZEN SECTION MEDIA

- Formulated to be an optimal cutting compound for frozen specimens
- Available in <u>FIVE</u> colors to optimize cutting efficiency in even the busiest labs
- Does not separate from tissue as it freezes like other frozen compounds
- · Snip top with cap for variable dispensing

CAN UF1000	Clear, 4oz	Each
CAN UF1002	Blue, 4oz	Each
CAN UF1003	Green, 4oz	Each
CAN UF1004	Pink, 4oz	Each
CAN UF1005	Yellow, 4oz	Each

MER 51024

### MERCEDES SCIENTIFIC° FREEZE SPRAY

- · Original and non-flammable options
- Non-flammable freezing compound may be used as a general purpose laboratory reagent
- Fast cooling on thermoplastic appliances or archwires
- · Evaporates quickly and leaves no residue
- Contains no ozone depleting compounds (CFC/HCFC)

 MER 51022
 Original, 10oz Can
 Each

 MER 51024
 Non-Flammable, 10oz Can
 Each

### MERCEDES SCIENTIFIC® SUB-XERO®

- Designed for immunohistochemical assays, frozen tissue specimens, tissue distribution of proteins and nucleic acid molecules at precise cutting temperatures
- · Leaves no residue on slides during staining

MER 5000 4oz (120mL) Bottle Each





MER 51022

#### CWS BHC001 (BLUE)



#### **HISTOCOOL™**

- Keeps ice frozen for long periods of time, stays frozen 6-8 hours on benchtop when kept in freezer overnight
- Freezer gel in sides and bottom keep ice frozen longer than traditional chilling methods
- Lid effectively limits surface from warm-air exposure, making ice last longer
- Space-saving design allows for stacking
- · Less trips to the freezer throughout the day

CWS BHC001	Blue, Holds 30 Blocks	Each
CWS RHC001	Red, Holds 20 Blocks	Each











#### **QUALITATIVE PLAIN CIRCLES**

- · Cellulose fiber papers with good wet strength
- Suitable for gravity, vacuum filtration systems and low suction filtration; not suitable for nitrogen analysis
- · Smooth surface with medium porosity

#### **P4 GRADE**

- For slightly slower uses, more retentive than P5 Grade
- Retains MgNH, PO,, silver sulfide and cadmium sulfide
- · Collects particulates or chromic acid mist in air monitors

FIS 098036A	4.25cm	100/Pack
FIS 098036B	5.5cm	100/Pack
FIS 098036C	7cm	100/Pack
FIS 098036D	9cm	100/Pack
FIS 098036E	11cm	100/Pack
FIS 098036F	12.5cm	100/Pack
FIS 098036G	15cm	100/Pack
FIS 098036H	18.5cm	100/Pack
FIS 098036J	20.5cm	100/Pack
FIS 098036K	24cm	100/Pack

#### **P5 GRADE**

- For clarifying and removing particles 1um and larger, and in applications not requiring low-ash or ashless paper
- For general lab use including clinical and spot testing and clarifying pharmaceuticals

FIS 09801AA	2.5cm	100/Pack
FIS 09801BB	3.5cm	100/Pack
FIS 09801J	4.25cm	100/Pack
FIS 09801A	7cm	100/Pack
FIS 09801B	9cm	100/Pack
FIS 09801C	11cm	100/Pack
FIS 09801D	12.5cm	100/Pack
FIS 09801E	15cm	100/Pack
FIS 09801F	19cm	100/Pack
FIS 09801G	24cm	100/Pack
FIS 09801H	33cm	100/Pack

#### **FLUTED QUALITATIVE CIRCLES**

#### P8 GRADE

- · Cellulose fiber papers with high wet strength
- Fluting provides more filter area, faster flow rates in conical funnel and reduces likelihood of tearing
- Suitable for coarse particles and gelatinous precipitates

FIS 0979014D	12.5cm	100/Pack
FIS 0979014E	15cm	100/Pack
FIS 0979014F	18.5cm	100/Pack
FIS 0979014G	24cm	100/Pack

#### WHATMAN™ QUALITATIVE CIRCLES

#### P4 GRADE

- Cellulose filters used in qualitative analytical techniques to determine and identify materials, and to clarify liquids
- Useful as a rapid filter for routine clean-up of biological fluids or organic extracts during analysis
- Retains course particles for monitoring air pollution
- Used when high flow rates in air pollution monitoring are required and the collection of fine particles is not critical
- Extremely fast filtering with excellent retention of coarse particles and gelatinous precipitates

FIS 09825A	70mm	100/Pack
FIS 09825B	90mm	100/Pack
FIS 09825C	110mm	100/Pack
FIS 09825D	125mm	100/Pack
FIS 09825E	150mm	100/Pack
FIS 09825F	185mm	100/Pack
FIS 09825G	240mm	100/Pack
FIS 09825H	320mm	100/Pack
FIS 09825J	42.5mm	100/Pack
FIS 09825K	55mm	100/Pack
FIS 09825L	270mm	100/Pack

PLEASE NOTE: SOME SIZES MAY TAKE LONGER TO SHIP. CALL FOR AVAILABILITY OR FOR MORE INFORMATION.

### **GLASS STAIN DISHES**

- Heavy, clear, soda-lime glass construction
- Hold a variety of different sized microscope slides

CAS SD10	Holds 10-20 Slides	Each
CAS SD20	Holds 20-Place Slide Rack (CAS SR20)	Each
CAS SD30	Holds 30-Place Slide Rack (CAS SR30)	Each
CAS SD60	Holds 60-Place Slide Rack (CAS SR60)	Each

#### STAINLESS STEEL STAIN RACKS

CAS SR20	20-Place for SD20	Each
CAS SR30	30-Place for SD30	Each
CAS SR60	60-Place for SD60	Each





#### **WHEATON™ GLASS STAINING DISHES**

- Dish and cover are made from soda-lime glass
- Nontarnishable rack withstands most staining solutions
- Hinged handle easily inserts and removes slides

#### **VERTICAL**

FIS 08815	Holds 10 Slides, Glass Cover, Rectangular Body, 5 Grooves	Each
FIS 08816	Holds 10 Slides, Plastic Cap, Rectangular Body, 5 Grooves	Each
FIS 08817	Holds 16 Slides, Glass Cover, 8 Grooves: Each Holds 2 Slides	Each
WTN 08813D	Holds 10 Slides (Immersion), Linerless Polypropylene Cap	Each
HORIZONTAL		
FIS 08813A	Holds 50 Slides, Glass Cover	Each







#### **ACRYLIC STAINING STATIONS**

- 3 flint-glass coplin jars with lids in acrylic stain station
- · Labeling surface for easy identification of stains

**MKL ML9860** 83/4"L x 6"W x 13/4"H Each



#### **20-SLIDE GLASS STAINING DISH**

- Dish, cover and rack made from soda-lime glass
- Stainless steel spring wire handle for manual staining
- · Bottom of rack is open for rapid immersion and draining

FIS 088121A	Dish & Cover	3/Pack
FIS 088121B	Slide Rack	3/Pack
FIS 088121C	Handle	6/Pack



#### **STAINING TRAYS & DISHES**

- · Made of solid stainless steel for durability
- Resistant to chemicals and corrosion and easy to clean
- · Covers are used as humidity chambers

MKL ML8944	Tray, 10-Slide, Handle	Each
MKL ML8945	Dish, 30-Slide, Cover	Each
MKL ML8946	Tray, 30-Slide, Handle	Each
MKL ML8947	Dish, 60-Slide, Cover	Each
MKL ML8948	Tray, 60-Slide, Handle	Each



#### **GLASS STAINING COPLIN JAR**

· Ground-on glass cover and wide, flat, round base

CAS SD05 Holds 10 Slides Each

#### PLASTIC STAINING COPLIN JAR

- Larger opening for easy placement and removal of slides
- Screw cap features a flat top, an easy grip knurled edge, deep thread and a seal ring to prevent leaks

BEL F44208 Holds 10 Slides Each



• 3-Piece Mesh Stainer Set: 31/8", 51/2", 77/8"

· Non-autoclavable, rust-resistant and dishwasher safe

TAR STRAIN Stainless Steel Strainers Set



#### SLIDE STAINING TRAY WITH LID

- · Clean, fresh working surface area every time
- · Reduce cross contamination and messy work areas
- · Compact size and recessed handles for easy transport
- · Pour spouts dispose unwanted liquid waste
- Dark lid protects slides for light sensitive applications

HEA HD15951A

Holds 8 Slides

4/Pack

#### EASY DIP™ SLIDE STAINING RACKS

- Holds up to 12 slides for use with slide staining jars
- Slides fit into individual slots for free passage and rapid drainage of staining fluids





#### **STAINTRAY™**

- Tough ABS plastic withstands a wide range of chemicals
- Suitable for routine staining requiring a humid chamber
- · Base holds excess solution dripping from slides
- · Rubber feet ensure greater base stability
- Units are stackable for space saving purposes
- · Clear cover allows for visual examination
- · Black cover used for fluorescent work

SIM M9181	Clear Lid, Holds 10 Slides	Each
SIM M9182	Black Lid, Holds 10 Slides	Each
SIM M9201	Clear Lid, Holds 20 Slides	Each
SIM M9202	Black Lid, Holds 20 Slides	Each



#### **SLIDETRAY**<sup>™</sup>

- · Holds 20 slides in an almost horizontal position
- · Stackable for reduced space on any shelf or lab counter
- Convenient microscope slide holder made of heavy-duty plastic to last even under the most adverse conditions
- Slides can be easily removed and inserted into position
- · Uses minimum space on any shelf or laboratory counter

SIM M75520W Slide Tray, White 10/Case

#### **EASYDIP™ SLIDE STAINING SYSTEM**

- Ideal for frozen sections, special stains and processes
- Staining jars made of resistant acetal plastic
- Resists most staining agents such as alcohol and xylene
- Wide stable base for greater stability while the inside is recessed allowing for a smaller reagent volume of 80mL
- Easy to clean and no metals to corrode

SIM M90012AS	5 Jars (1 of Each Color), 1 Rack (SIM M90512DGY)	Kit
SIM M90012B	Blue Staining Jar	6/Case
SIM M90012G	Green Staining Jar	6/Case
SIM M90012P	Pink Staining Jar	6/Case
SIM M90012W	White Staining Jar	6/Case
SIM M90012Y	Yellow Staining Jar	6/Case



#### **EASYDIP™ SLIDE STAINING KIT**

 Kit Includes: (1) Anodized Aluminum Rack, (6) Assorted Color Jars, (1) Slide Staining Rack (SIM M90512DGY)

Anodized Aluminum Holder

Each

SIM M90612AS	EasyDip" Slide Staining Kit	Kit
4		



#### **SLIDE STAINING RACK**

- · Stainless steel resists corrosion under normal use
- Adjustable to fit trays or sinks up to (21") inside
- · Adjustment is made by sliding the clips to fit the tank

**TAL 954002** 23¾"L x 3½"W Each



**SIM M906** 

# CASSETTES



## **SuperFLO™ CASSETTES**

The SuperFlo™ Cassette, a new generation emerging in the Simport® Histology Family.

- · Molded from a special high-density polymer, these cassettes keep specimens safely submerged in liquid and are resistant to the chemical action of most histological solvents.
- · Lateral vent design for a maximum fluid exchange and efficient tissue processing.
- · Sloped walls, create a basket for small tissue management.
- · Available in loose, sleeved, or taped packaging for easy printer loading.

#### **SUPERFLO™ ONE-PIECE MODEL (LOOSE)**

This model is a one-piece molded cassette so no separate lids to manage. It can be opened and closed up to five times guaranteed, securely relocking to avoid specimen loss.

	О	p	е	n	İ	n	g	

SIM M530	Routine Embedding Cassette	0.91 x 3.67mm	1,500/Case
SIM M531	Biopsy Cassette	0.64 x 0.64mm	1,500/Case
SIM M532	Microbiopsy Cassette	0.30 x 0.30mm	1,500/Case

#### **TAPED SUPERFLO™ ONE-PIECE MODEL**

Add the color number provided below after the base item#, then add a "T" at the end to indicate "Taped".

_				٠		
O	n	ρ	n	ı	n	c

SIM M430	Routine Embedding Cassette	0.91 x 3.67mm	1,500/Case
SIM M431	Biopsy Cassette	0.64 x 0.64mm	1,500/Case
SIM M432	Microbiopsy Cassette	0.30 x 0.30mm	1,500/Case

#### **SUPERFLO™ HINGED-HOOK MODEL (LOOSE)**

This model has a special hinge that holds the base and the cover together. The hinge design allows the cassette to be opened and closed as often as required.

#### **SLEEVED** SUPERFLO™ HINGED-HOOK MODEL

Add the color number provided after the base item#, then add "SL" at the end to indicate "Sleeved". \*Sleeved cassettes come 750/Cs.

О			

SIM M533	Routine Embedding Cassette	0.91 x 3.67mm 1,500/Case*
SIM M534	Biopsy Cassette	0.64 x 0.64mm 1,500/Case*
SIM M535	Microbiopsy Cassette	0.30 x 0.30mm 1,500/Case*

#### **TAPED SUPERFLO™ HINGED-HOOK MODEL**

Add the color number provided after the base item#, then add a "T" at the end to indicate "Taped".

#### Opening


SIM M433	Routine Embedding Cassette	0.91 x 3.67mm	1,500/Case
SIM M434	Biopsy Cassette	0.64 x 0.64mm	1,500/Case
SIM M435	Microbiopsy Cassette	0.30 x 0.30mm	1,500/Case

#### **CHOOSE FROM THE FOLLOWING COLORS:**

(Add the color number to the end of the base Item#.)



GRAY (9)






YELLOW (5)









#### **EMBEDDING CASSETTES**

- Suitable for holding and identifying tissue samples in processing, embedding, and sectioning procedures
- · Fit securely in microtome chuck adapters
- Resists the chemical action of most histological solvents
- Standard metal lids (SIM M481) keeps specimens in complete safety during processing
- Slanted writing surface accepts markings easily, permitting sample identification throughout all stages of embedding and long afterwards when in archives
- Each case contains (3) dispenser boxes of 500 cassettes

SIM M480	Embedding, No Lid	1,500/Case	
SIM M481	Metal Lids for M480	25/Pack	

#### **MACROSETTE™ CASSETTES**

- Designed to hold larger tissue specimens during the embedding process, as well as in a storage cabinet
- TWICE AS HIGH AS HISTOSETTE® CASSETTES (13MM)
- · One-piece integral lid eliminates the need for steel lids
- Can be opened and closed as often as necessary and always relock securely without danger of specimen loss
- Large labeling areas are located on three sides of the cassettes for convenience
- Each case contains (3) dispenser boxes of 250 cassettes

**SIM M512** Embedding, White 750/Case



#### **CHOOSE FROM THE FOLLOWING COLORS:**

(Add the color number to the end of the base Item#.)

WHITE (2)

BLUE (6)

PINK (3)

PEACH (7)

ORANGE (11)

GREEN (4)

TAN (8)

AQUA (12)

YELLOW (5)

GRAY (9)

### **HISTOSETTE® I CASSETTES**

- Molded from a special high density polymer
- · Keep specimens safely submerged in liquid
- · Resists the chemical action of most histological solvents
- Efficient flow-through slots maximize fluid exchange and ensure proper drainage
- · One-piece integral lid eliminates the need for steel lids
- Just snap apart and conveniently lock the lid into the base of the cassette
- Can be opened and closed as necessary and always relock securely without danger of specimen loss
- M490/M491: Not suitable for automated printers
- M498/M499: Anterior writing area is 45° instead of 30°
- Each case contains (3) dispenser boxes of 500 cassettes

SIM M490	Tissue, 30° Angle	1,500/Case
SIM M491	Biopsy, 30° Angle	1,500/Case
SIM M498	Tissue, 45° Angle	1,500/Case
SIM M499	Biopsy, 45° Angle	1,500/Case

### **HISTOSETTE® II CASSETTES**

- · Molded from a special high density polymer
- · Keep specimens safely submerged in liquid
- Resists the chemical action of most histological solvents
- Efficient flow-through slots maximize fluid exchange and ensure proper drainage
- Anterior writing area is slanted at a 45° angle

### M485/M486

- Most convenient for cassette labeling instruments since covers are already removed from cassettes and are packaged separately in the case
- Each case contains (2) dispenser boxes of 500 cassettes and (1) dispenser box of 1,000 covers

### M492/M493

- Can be opened and closed as necessary and always relock securely without danger of specimen loss
- One-piece integral lid eliminates the need for steel lids
- Each case contains (3) dispenser boxes of 500 cassettes

SIM M485	Tissue, Lids Separate	1,000/Case
SIM M486	Biopsy, Lids Separate	1,000/Case
SIM M492	Tissue, Lids Attached	1,500/Case
SIM M493	Biopsy, Lids Attached	1.500/Case

### **CHOOSE FROM THE FOLLOWING COLORS:**























### **SLIMSETTE™ CASSETTES**

- Unique recessed cover allows cassettes to be stacked in automatic labeling machines and in storage cabinets
- Molded from a special high density polymer
- · Keep specimens safely submerged in liquid
- Resists the chemical action of most histological solvents
- Can be opened and closed as often as necessary and always relock securely without danger of specimen loss
- · Large labeling areas for easy identification
- Anterior writing area is slanted at a 45° angle
- M509: 114 openings each measuring 1 x 5 mm
- M510: 392 openings ensure maximum fluid exchange
- M511: 4 square compartments each measuring 13mm
- Each case contains (3) dispenser boxes of 500 cassettes with covers assembled

SIM M509	Tissue, 1 Compartment	1,500/Case
SIM M510	Biopsy, 1 Compartment	1,500/Case
SIM M511	Tissue, 4 Compartments	1,500/Case

### **SWINGSETTE™ CASSETTES**

- · Molded from a special high-density polymer
- · Keep specimens safely submerged in liquid
- · Resists the chemical action of most histological solvents
- Efficient flow-through slots maximize fluid exchange and ensure proper drainage.
- Special hinge holds the base and cover together, allows cassettes to be opened and closed as often as necessary
- Cover can be removed and re-inserted easily without danger of specimen loss
- · Large tab for convenient and easy opening of lid
- · Anterior writing area is slanted at a 45° angle

### M515/M516

• Each case contains (3) dispenser boxes of 500 cassettes

### M517/M518

- · Covers are packaged separately in the case
- Each case contains (2) dispenser boxes of 500 cassettes and (1) dispenser box of 1,000 covers

SIM M515	Tissue, Lids Attached	1,500/Case
SIM M516	Biopsy, Lids Attached	1,500/Case
SIM M517	Tissue, Lids Separate	1,000/Case
SIM M518	Biopsy, Lids Separate	1,000/Case

### **CHOOSE FROM THE FOLLOWING COLORS:**

(Add the color number to the end of the base Item#.)



### **UNISETTE™ CASSETTES**

- GREAT ALTERNATIVE TO SAKURA UNI-CASSETTE®
- · Molded from a special high density polymer
- · Keep specimens safely submerged in liquid
- · Resists the chemical action of most histological solvents
- Efficient flow-through slots maximize fluid exchange and ensure proper drainage
- · One-piece integral lid eliminates the need for steel lids
- Snap-latch and hinge-lock design stop early separation of base and lid and allows for one-hand operation
- Can be opened and closed as often as necessary and always relock securely without danger of specimen loss
- · Anterior writing area is slanted at a 35° angle
- Each case contains (3) dispenser boxes of 500 cassettes

SIM M505	Tissue	1,500/Case
SIM M506	Biopsy	1,500/Case

### MICROMESH™ BIOPSY CASSETTES

- NO BIOPSY PADS REQUIRED
- Offers 1,676 square openings (0.38mm) allowing for a greatly improved fluid exchange and drainage without having to use biopsy pads
- Large anterior slots in both cassette and cover ensure the Micromesh biopsy cassette will sink rapidly
- Large square compartment measuring 27mm is perfect even for needle biopsies
- Cover does not protrude above the cassette to save space allowing more cassettes to be stacked in automatic labeling machines and tissue processors
- · Molded from a special high density polymer
- · Keep specimens safely submerged in liquid
- Resists the chemical action of most histological solvents
- Can be opened and closed as often as necessary and always relock securely without danger of specimen loss
- Anterior writing area is slanted at 45° to make cassette more suitable to use with automated cassette printers
- Each case contains (4) dispenser boxes of 250 cassettes with covers assembled

SIM M507	Single Compartment	1,000/Case
SIM M508	Four Compartments	1,000/Case

### **CHOOSE FROM THE FOLLOWING COLORS:**


(Add the color number to the end of the base Item#.)


WHITE (2)
 BLUE (6)
 LILAC (10)
 PINK (3)
 PEACH (7)
 ORANGE (11)
 GREEN (4)
 TAN (8)
 AQUA (12)
 YELLOW (5)
 GRAY (9)


















### UNISETTE™ II (FOR PRIMERA® PRINTERS)

### IN QUICKLOAD™ STACKS

- Disposable plastic cassettes hold tissue specimens during processing and embedding as well as in storage
- Molded from a special high-density polymer
- · Keep specimens safely submerged in liquid
- Resists the chemical action of most histological solvents
- Efficient flow-through slots maximize fluid exchange and ensure proper drainage
- Can be opened or closed as required and always relock securely without danger of specimen loss

### **MANUAL FEED PRINTER**

3 Dispenser Boxes of 500 Cassettes

SIM M405	Unisette™ II, Tissue	1,500/Case
SIM M406	Unisette™ II, Biopsy	1,500/Case

### **AUTOMATED FEED PRINTER**

25 Stacks of Taped Cassettes

SIM M405T	Unisette™ II, Tissue	1,000/Case
SIM M406T	Unisette™ II, Biopsy	1,000/Case

### MICROMESH™ (FOR PRIMERA® PRINTERS)

### IN QUICKLOAD™ STACKS

- Square openings (0.38mm) allow for a greatly improved fluid exchange without having to use biopsy pads
- Large anterior and posterior slots in both cassette and cover ensure the cassette will be submerged rapidly
- Mesh ensures efficient fluid exchange and drainage
- One-piece integral lid eliminates the need for steel lids
- Recessed cover allows more cassettes to be stacked in automatic labeling machines and storage drawers

### **MANUAL FEED PRINTER**

4 Dispenser Boxes of 250 Cassettes

SIM M407 Micromesh™, Biopsy 1,000/Case


### **ROBOTIC FEED PRINTER**

50 Stacks of Taped Cassettes

**SIM M407T** Micromesh™, Biopsy 2,000/Case

### CHOOSE FROM THE FOLLOWING COLORS:

(Add the color number to the end of the base Item#.)



SPECIALLY DESIGNED FOR AUTOMATED & MANUAL PRIMERA® CASSETTE PRINTERS

### **SLEEVED CASSETTES**

- ESPECIALLY FOR USE WITH THERMOFISHER PRINTERS
- Suitable for hoppers; QuickLoad™ sleeve will load in cassette labeling instruments in one simple operation
- · Transparent sleeve allows viewing of cassettes in order to confirm there is no jam in sleeve during printing
- Molded from a special high density acetal polymer
- · Keep specimens safely submerged in liquid
- Totally resistant to the chemical action of most solvents used in histology laboratories
- · Efficient flow-through slots maximize fluid exchange and ensure proper reagent drainage
- Each case contains 10 sleeves of 75 cassettes
- SWINGSETTE™: 10 bags of 75 lids packaged separately

SIM M480SL	Embedding	750/Case
SIM M485SL	Histosette® II, Tissue	750/Case
SIM M486SL	Histosette® II, Biopsy	750/Case
SIM M507SL	Micromesh™, Biopsy	750/Case
SIM M509SL	Slimsette™, Tissue	750/Case
SIM M510SL	Slimsette™, Biopsy	750/Case
SIM M517SL	Swingsette <sup>™</sup> , Tissue	750/Case
SIM M518SL	Swingsette <sup>™</sup> , Biopsy	750/Case



are compatible with ThermoFisher Printers!

### PICTURED (RIGHT) THERMO SCIENTIFIC° PRINTMATE AS CASSETTE PRINTER

- Eliminate waste and streamline labeling process by printing cassettes on-demand at grossing station
- Small footprint allows the printer to fit easily next to grossing station
- Eliminate manual handwriting processes
- Print only the cassettes needed, right at the point of use
- · Prints each cassette in 10 seconds or less
- · Print barcodes on cassettes to automate downstream processes



### **CHOOSE FROM THE FOLLOWING COLORS:**

(Add the color number to the end of the base Item#.)


- WHITE (2)
- BLUE (6)
- ORANGE (11)
  - LILAC (10)

- PINK (3) GREEN (4)
- TAN (8)

PEACH (7)

AQUA (12)

- YELLOW (5)
- GRAY (9)



































### **TAPED CASSETTES**

- SPECIFIC FOR LEICA® & SAKURA® INK JET PRINTERS
- Just load the QuickLoad  $\mbox{}^{\text{\tiny{M}}}$  stack in the magazine and you are ready for printing
- Packaged in stacks of 40; each case contains 50 stacks for a total of 2,000 cassettes

SIM M480T	Embedding	2,000/Case
SIM M492T	Histosette® II, Tissue	2,000/Case
SIM M493T	Histosette® II, Biopsy	2,000/Case
SIM M505T	Unisette™, Tissue	2,000/Case
SIM M506T	Unisette™, Biopsy	2,000/Case
SIM M507T	Micromesh <sup>™</sup> , Biopsy	2,000/Case
SIM M509T	Slimsette™, Tissue	2,000/Case
SIM M510T	Slimsette™, Biopsy	2,000/Case
SIM M517T	Swingsette <sup>™</sup> , Tissue	2,000/Case
SIM M518T	Swingsette™, Biopsy	2,000/Case

### **CHOOSE FROM THE FOLLOWING COLORS:**

(Add the color number to the end of the base Item#.)

WHITE (2) BLUE (6)

PINK (3)

GREEN (4)

YELLOW (5)

PEACH (7)

— TAN (8)

GRAY (9)

76

LILAC (10)

AQUA (12)

ORANGE (11)

### **MODULAR STORAGE DRAWERS**

- · Made of high impact resistant plastic
- Permanent storage and identification of up to 165 embedding rings or 250 cassettes per drawer
- Stackable to any convenient height courtesy of interlocking ridges on top and bottom
- · Identification labels included
- Dimensions: 151/8"L x 91/8"W x 21/8"H

SIM M4956

**Modular Drawers** 

6/Case





TISSUE CASSETTE BASKETS













### **TISSUE INDEX MARKERS**

- See when a cassette has been removed from storage
- 13/4" markers stand 1/8" above cassettes for easy visibility
- · Features case number, date and initials of person fields
- 4 vertical spaces for other patient information
- · COLORS: Yellow, Pink, Green, Blue, Lilac, White

LSS TIM100 Tissue Index Marker 1,000/Pack

### **TISSUE CASSETTE BASKETS**

- Designed to accommodate most major cassette brands
- · Organize cassettes with removable dividers
- Deluxe size includes hinged lid does lock into place while small and medium sizes include lid with handle that does <u>NOT</u> lock into place
- Dimensionally identical to Sakura® cassette baskets

TOB SB01	Slide Basket Adapter	Each
TOB B050	Holds 50 Cassettes, Small, 5½"L x 4½"W x 1¼"H	Each
TOB B100	Holds 100 Cassettes, Large, 9"L x 5½"W x 1¼"H	Each
TOB B150	Holds 150 Cassettes, Deluxe, 9"L x 6"W x 21/8"H	Each



# TANNER SCIENTIFIC® CASSETTE STORAGE SYSTEM

· For the storage of paraffin embedded cassettes

• Includes: Dividers, (8) Drawers

• Capacity: 500 Cassettes

• Dimensions: 12" x 12" x 4"

· Shipped flat, easy to assemble

• MAY ALSO BE USED TO STORE SLIDES

**TNR CASBX20** 

Cassette Storage System

Each

TANNERS



CASSETTE STORMS

### TISSUE CASSETTE STORAGE UNITS

- · Holds approximately 600 cassettes per unit
- Holds approximately 3,600 cassettes per case of 6 units
- Economical, lightweight, compact design
- SHIPPED FULLY ASSEMBLED
- Dimensions: 81/4"W x 141/2"D x 41/2"H



### **INCLUDES:**

- (12) High-Impact Polystyrene Plastic Trays
- (12) Sturdy Corrugated Drawers with Front Labeling
- (6) Corrugated Drawer Enclosures for Dust-Proof Protection
- (1) Outer Carton for Long Term Storage

LSS CB21

Regular Storage Unit

6/Case



LSS ACT2

### INCLUDES:

- (12) Corrugated Drawers with Dividers 11/4" Apart
- (6) Corrugated Drawer Enclosures for Dust-Proof Protection
- (1) Outer Carton For Long-Term Storage

LSS ACT2

Corrugated Storage Unit

6/Case



# TN60 WITH UV-C DISINFECTION CRYOSTAT

### item # TNR TN60

Dimensions (WxDxH): 26¾" x 23-5/8" x 43¾"

Weight (w/o accessories): 280 lbs


Electrical Parameters: 115V; 50/60 Hz





# MERCEDES SCIENTIFIC® MICROSCOPE SLIDES

- Pathological and cytological applications in diagnostics and biological-medical research of cells and tissues
- Imported from one of the finest glass manufacturing facilities in Germany
- · Exceptional clarity and cleanliness
- Packaged in thick foil or "tropical packaging" to protect from humidity even in very warm climates
- Slides are packaged 72 per box, 20 boxes per case, 10 gross per case (1,440 slides), unless noted



### **45° BEVELED SLIDES**

MER 7200/45/BL	Blue	1,440/Case
MER 7200/45/GR	Green	1,440/Case
MER 7200/45/OR	Orange	1,440/Case
MER 7200/45/PK	Pink	1,440/Case
MER 7200/45/WH	White	1,440/Case
MER 7200/45/YL	Yellow	1,440/Case

### **45° BEVELED ADHESIVE SLIDES**

MER 7255/45/BL	Blue	1,440/Case
MER 7255/45/GR	Green	1,440/Case
MER 7255/45/OR	Orange	1,440/Case
MER 7255/45/PK	Pink	1,440/Case
MER 7255/45/WH	White	1,440/Case
MER 7255/45/YL	Yellow	1,440/Case

### 90° STRAIGHT SLIDES

MER 7200/90/BL	Blue	1,440/Case
MER 7200/90/GD	Gold	1,440/Case
MER 7200/90/GR	Green	1,440/Case
MER 7200/90/LG	Light Green	1,440/Case
MER 7200/90/LV	Lavender	1,440/Case
MER 7200/90/OR	Orange	1,440/Case
MER 7200/90/PK	Pink	1,440/Case
MER 7200/90/RD	Red	1,440/Case
MER 7200/90/SL	Silver	1,440/Case
MER 7200/90/SM	Salmon	1,440/Case
MER 7200/90/WH	White	1,440/Case
MER 7200/90/YL	Yellow	1,440/Case

### 90° STRAIGHT ADHESIVE SLIDES

MER 7255/90/GD         Gold         1,440/Case           MER 7255/90/GR         Green         1,440/Case           MER 7255/90/LV         Lavender         1,440/Case           MER 7255/90/OR         Orange         1,440/Case           MER 7255/90/PK         Pink         1,440/Case           MER 7255/90/SL         Silver         1,440/Case           MER 7255/90/SM         Salmon         1,440/Case           MER 7255/90/WH         White         1,440/Case           MER 7255/90/YL         Yellow         1,440/Case	MER 7255/90/BL	Blue	1,440/Case
MER 7255/90/LV         Lavender         1,440/Case           MER 7255/90/OR         Orange         1,440/Case           MER 7255/90/PK         Pink         1,440/Case           MER 7255/90/SL         Silver         1,440/Case           MER 7255/90/SM         Salmon         1,440/Case           MER 7255/90/WH         White         1,440/Case	MER 7255/90/GD	Gold	1,440/Case
MER 7255/90/OR         Orange         1,440/Case           MER 7255/90/PK         Pink         1,440/Case           MER 7255/90/SL         Silver         1,440/Case           MER 7255/90/SM         Salmon         1,440/Case           MER 7255/90/WH         White         1,440/Case	MER 7255/90/GR	Green	1,440/Case
MER 7255/90/PK         Pink         1,440/Case           MER 7255/90/SL         Silver         1,440/Case           MER 7255/90/SM         Salmon         1,440/Case           MER 7255/90/WH         White         1,440/Case	MER 7255/90/LV	Lavender	1,440/Case
MER 7255/90/SL         Silver         1,440/Case           MER 7255/90/SM         Salmon         1,440/Case           MER 7255/90/WH         White         1,440/Case	MER 7255/90/OR	Orange	1,440/Case
MER 7255/90/SM         Salmon         1,440/Case           MER 7255/90/WH         White         1,440/Case	MER 7255/90/PK	Pink	1,440/Case
MER 7255/90/WH White 1,440/Case	MER 7255/90/SL	Silver	1,440/Case
	MER 7255/90/SM	Salmon	1,440/Case
<b>MER 7255/90/YL</b> Yellow 1,440/Case	MER 7255/90/WH	White	1,440/Case
	MER 7255/90/YL	Yellow	1,440/Case

NOTE: LIGHT GREEN AND RED ARE NOT AVAILABLE COLORS FOR THE PLATINUMLINE  $^\circ$  STARFROST  $^\sim$  ADHESIVE SLIDES.



90° STRAIGHT COLOR SLIDES



90° STRAIGHT ADHESIVE SLIDES





CLIPPED CORNER SLIDES

### **CLIPPED CORNER SLIDES**

MER 7200/90/WH/CC	White, 90°	1,440/Case
MER 7200/90/BL/CC	Blue, 90°	1,440/Case
MER 7200/90/GR/CC	Green, 90°	1,440/Case
MER 7200/90/PK/CC	Pink, 90°	1,440/Case
MER 7200/90/OR/CC	Orange, 90°	1,440/Case
MER 7200/90/YL/CC	Yellow, 90°	1,440/Case

### **ADHESIVE CLIPPED CORNER SLIDES**

MER 7255/90/WH/CC	White, 90°	1,440/Case
MER 7255/90/BL/CC	Blue, 90°	1,440/Case
MER 7255/90/GR/CC	Green, 90°	1,440/Case
MER 7255/90/PK/CC	Pink, 90°	1,440/Case
MER 7255/90/OR/CC	Orange, 90°	1,440/Case
MER 7255/90/YL/CC	Yellow, 90°	1,440/Case

### **ADVANCED ADHESIVE SLIDES**

- · Hydrophilic, coated slides were specially developed for use in Immunohistochemistry and harsher applications
- Stronger adherence means no tissue falls off anymore
- Improved light transmission, no autofluorescence
- · Made to a new low iron, soda-lime glass formulation

MER 8888 White, 90°, Adhesive 1.440/Case

### **PATIENT CONTROL SLIDES**

- · Ideal for staining procedures that require a control
- Spaces for the patient sample and the control material
- 2 samples can be separated by the red line, which is printed on the reverse side for reagent movement
- · Coated to provide maximum tissue adhesion

**MER 9999** Adhesive, 90°, Coated 1.440/Case

### **1.1MM THICK SLIDES**


MER 7800	White, 90°	1,440/Case
MER 7850	Twin Frost, 45°	1,440/Case
MER 7860	Twin Frost, 90°	1.440/Case

### **1.2MM THICK SLIDES**


- · Highly resistant to corrosion and moisture
- · Made of high-quality glass with edges that are ground smooth to help eliminate sharp cutting surfaces

MER 7200/90/WH/1.2MM White, 90° 1.440/Case











### **PLAIN SLIDES**

- Clear, precleaned, premium-quality, water-white glass
- · Corrosion resistant glass with smoothly polished edges
- Cost-effective slides for general-purpose applications
- Unmatched transparency and absence of glass color tint

MER 7500	45° Edge	1,440/Case
MER 7515	90° Edge	1,440/Case

### **TWIN FROST SLIDES**

- · Non-corrosive, non-fog, precleaned, water-white glass
- · For blood smears or routine microscopy applications
- · Brilliant white opaque coating retains pen markings
- Superior readability for easier slide identification

MER 7550	45° Edge	1,440/Case
MER 7760	90° Edge	1,440/Case

### **CIRCLE SLIDES**

- Standard diagnostic slides with 1 or 2 printed white rings serve as a microscope to find centrifuged cells easier
- · Suitable for use with cytocentrifuges

### WHITE, 90°

**MER K177** 

MER 7200/90/WH/1C	1 Circle	1,440/Case
MER 7200/90/WH/1C/BX	1 Circle	72/Box
MER 7200/90/WH/2C	2 Circle	1,440/Case
WHITE, 90°, ADHESIVE		
MER 7255/90/WH/1C	1 Circle	1,440/Case
MER 7255/90/WH/1C MER 7255/90/WH/1C/BX	1 Circle	1,440/Case 72/Box

White, 90°, Adhesive

### **ETCHED SLIDES**


- · Double layer frosted end of two contrasting colors
- · Available with wettable or coated glass surface
- Used for accurate and repeated identification of cultured or manipulated cells growing on surface


### WHITE (TOP), BLACK (BOTTOM)

MER 7900/90/BW	90°, Regular	1,440/Case
MER 7955/90/BW	90°, Adhesive	1,440/Case
WHITE (TOP), BLUE (BOTTOM)		
MER 7900/90/WB	90°, Regular	1,440/Case
MER 7955/90/WB	90°, Adhesive	1,440/Case

### **THERMAL PAINT SLIDES**

**MER 7255/90/TP/PATH** 90° Edge 1,440/Case





72/Box





# TANNER SCIENTIFIC® MICROSCOPE SLIDES

- Ideal for everyday use in histology, cytology, hematology and other laboratory environments
- Offer exceptional quality at a price that fits any budget

### 45° SLIDES

TNR BLU45	Blue	1,440/Case
TNR GLD45	Gold	1,440/Case
TNR GRN45	Green	1,440/Case
TNR LAV45	Lavender	1,440/Case
TNR LTGRN45	Light Green	1,440/Case
TNR ORG45	Orange	1,440/Case
TNR PNK45	Pink	1,440/Case
TNR SALMON45	Salmon	1,440/Case
TNR SILVER45	Silver	1,440/Case
TNR WHT45	White	1,440/Case
TNR YEL45	Yellow	1,440/Case

### 90° SLIDES


TNR BLU90	Blue	1,440/Case
TNR GLD90	Gold	1,440/Case
TNR GRN90	Green	1,440/Case
TNR LAV90	Lavender	1,440/Case
TNR LTGRN90	Light Green	1,440/Case
TNR ORG90	Orange	1,440/Case
TNR PNK90	Pink	1,440/Case
TNR SALMON90	Salmon	1,440/Case
TNR SILVER90	Silver	1,440/Case
TNR WHT90	White	1,440/Case
TNR YEL90	Yellow	1,440/Case


### 90° CLIPPED CORNER SLIDES

TNR BLU90CC	Blue	1,440/Case
TNR GLD90CC	Gold	1,440/Case
TNR GRN90CC	Green	1,440/Case
TNR LAV90CC	Lavender	1,440/Case
TNR LTGRN90CC	Light Green	1,440/Case
TNR ORG90CC	Orange	1,440/Case
TNR PNK90CC	Pink	1,440/Case
TNR SALMON90CC	Salmon	1,440/Case
TNR SILVER90CC	Silver	1,440/Case
TNR WHT90CC	White	1,440/Case
TNR YEL90CC	Yellow	1,440/Case

### 90° WHITE SLIDES (NO MARKINGS)

TNR WHT90NM	Plain	1,440/Case
TNR WHT90NMAD	Adhesive	1,440/Case
TNR WHT90CCNM	Clipped Corner	1,440/Case





### 45° ADHESIVE SLIDES

TNR BLU45AD	Blue	1,440/Case
TNR GLD45AD	Gold	1,440/Case
TNR GRN45AD	Green	1,440/Case
TNR LAV45AD	Lavender	1,440/Case
TNR LTGRN45AD	Light Green	1,440/Case
TNR ORG45AD	Orange	1,440/Case
TNR PNK45AD	Pink	1,440/Case
TNR SALMON45AD	Salmon	1,440/Case
TNR SILVER45AD	Silver	1,440/Case
TNR WHT45AD	White	1,440/Case
TNR YEL45AD	Yellow	1,440/Case

### 90° ADHESIVE SLIDES

TNR BLU90AD	Blue	1,440/Case
TNR GLD90AD	Gold	1,440/Case
TNR GRN90AD	Green	1,440/Case
TNR LAV90AD	Lavender	1,440/Case
TNR LTGRN90AD	Light Green	1,440/Case
TNR ORG90AD	Orange	1,440/Case
TNR PNK90AD	Pink	1,440/Case
TNR SALMON90AD	Salmon	1,440/Case
TNR SILVER90AD	Silver	1,440/Case
TNR WHT90AD	White	1,440/Case
TNR YEL90AD	Yellow	1,440/Case

### 90° CLIPPED CORNER ADHESIVE SLIDES

TNR BLU90ADCC	Blue	1,440/Case
TNR GLD90ADCC	Gold	1,440/Case
TNR GRN90ADCC	Green	1,440/Case
TNR LAV90ADCC	Lavender	1,440/Case
TNR LTGRN90ADCC	Light Green	1,440/Case
TNR ORG90ADCC	Orange	1,440/Case
TNR PNK90ADCC	Pink	1,440/Case
TNR SALMON90ADCC	Salmon	1,440/Case
TNR SILVER90ADCC	Silver	1,440/Case
TNR WHT90ADCC	White	1,440/Case
TNR YEL90ADCC	Yellow	1,440/Case

### PLAIN SLIDES

TNR PLN45	45° Plain	1,440/Case
TNR PLN90	90° Plain	1,440/Case

### SINGLE FROSTED SLIDES

TNR SGLFROST45	45° Single Frost	1,440/Case
TNR SGLFROST90	90° Single Frost	1,440/Case

### **DOUBLE FROSTED SLIDES**

TNR DFROST45	45° Double Frost	1,440/Case
TNR DFROST90	90° Double Frost	1,440/Case



# MERCEDES SCIENTIFIC® ADVANCED ADHESIVE SLIDES

- Hydrophilic, coated slides were specially developed for use in immunohistochemistry and harsher applications
- · Stronger adherence means no tissue falls off anymore
- Improved light transmission, no autofluorescence
- · Made to a new low iron, soda-lime glass formulation

**MER 8888** 90°, White, Adhesive 1,440/Case

# MERCEDES SCIENTIFIC® IHC+ MICROSCOPE SLIDES

- Dimensions: 12 x 75 x 1mm
- · Fortified adhesion brought about by new coat material
- High transmittance for fluorescence observation
- Designed to be compatible with immunohistostainers

**MER 1111** 45° Adhesive 100/Box





### **TOMO® HYDROPHILIC ADHESION SLIDES**

- Superior tissue adhesion, staining performance, and slide surface fluidics
- Manufactured using optically clear, exceptionally flat, low fluorescence glass
- · Slides are finished with clipped
- Proprietary slide coating ensures consistent performance, even after exposure to heat stress and humid environmental conditions
- Recommended for use on Ventana Benchmark staining platform and automated IHC platforms

TOM 11	45° White, 12 mm x 75 mm Thickness 1.0mm	1,000/Case
TOM 12	45° Yellow, 25 mm x 75 mm, Thickness 1.0 mm	1,000/Case
TOM 13	45° Pink, 25 mm x 75 mm, Thickness 1.0 mm	1,000/Case
TOM 14	45° Blue, 25 mm x 75 mm, Thickness 1.0 mm	1,000/Case
TOM 15	45° Green, 25 mm x 75 mm, Thickness 1.0 mm	1,000/Case



# TANNER SCIENTIFIC® COVERGLASS

- Suitable for a host of laboratory applications: cytology, histopathology, microbiology and life sciences
- High-quality, thin, disposable coverglass provides the perfect clarity techs want and need while their non-stick treatment makes them easy to handle
- Packaged 1oz per box, 10oz per case
- Thickness: #1 = 0.13-0.17mm, #1.5 = 0.16-0.19mm
- Note: 2oz per box, 20 boxes per case\*

TNR 2222	22 x 22, #1	10oz/Case
TNR 2440	24 × 40, #1	10oz/Case
TNR 2450	24 x 50, #1	10oz/Case
TNR 24501.5	24 x 50 #1.5	10oz/Case
TNR 2450AC*	24 x 50, #1	20oz/Case
TNR 2455	24 x 55, #1	10oz/Case
TNR 2460	24 x 60, #1	10oz/Case
TNR 24601.5	24 x 60, #1.5	10oz/Case



# TANNER SCIENTIFIC® COVERSLIPPING FILM

- Length: 75 Meters, Width: 24mm
- Usable in all automatic film coverslippers
- Free of aromatic solvents like Xylol and Benzol
- "Vario Texture" surface structure allows for a better distribution of solvent drops without any effects on the high transparency of the film after coverslipping
- Can use for a minimum of 12 months after manufacturing

TNR TAPE Coverslipping Film Each

# MERCEDES SCIENTIFIC® COVERGLASS


- Made of a water-white colorless borosilicate glass with an optimum hydrolytic resistance
- Ultra thin coverglass German manufactured
- Performs exceptionally well on all automated glass coverslippers and is Leica\* recommended
- Very fine glass particles powder the surface to prevent sticking caused by stacking flat coverglass
- Thickness: #1 = 0.13-0.17mm, #1.5 = 0.16-0.19mm











### **ROUND COVERGLASS**

MER R0010	10mm, #1	5oz/Case
MER R0011	11mm, #1	5oz/Case
MER R0012	12mm, #1	5oz/Case
MER R00121.5	12mm, #1.5	5oz/Case
MER R0013	13mm, #1	5oz/Case
MER R0015	15mm, #1	5oz/Case
MER R0018	18mm, #1	5oz/Case
MER R00181.5	18mm, #1.5	5oz/Case
MER R00221.5	22mm, #1.5	10oz/Case
MER R0025	25mm, #1	10oz/Case
MER R00251.5	25mm, #1.5	10oz/Case
MER R0028	28mm, #1	10oz/Case


### **#1 THICK COVERGLASS**


MER R1515	15 x 15mm	10oz/Case
MER R1818S	18 x 18mm	10oz/Case
MER R2222	22 x 22mm	10oz/Case
MER R2240	22 x 40mm	10oz/Case
MER R2250	22 x 50mm	10oz/Case
MER R2430	24 x 30mm	10oz/Case
MER R2440	24 x 40mm	10oz/Case
MER R2445	24 x 45mm	10oz/Case
MER R2450	24 x 50mm	10oz/Case
MER R2455	24 x 55mm	10oz/Case
MER R2460	24 x 60mm	10oz/Case
MER R2525	25 x 25mm	10oz/Case
MER R2530	25 x 30mm	10oz/Case
MER R2540	25 x 40mm	10oz/Case
MER R2550	25 x 50mm	10oz/Case


### **#1.5 THICK COVERGLASS**


MER R22221.5	22 x 22mm	10oz/Case
MER R22401.5	22 x 40mm	10oz/Case
MER R22501.5	22 x 50mm	10oz/Case
MER R24301.5	24 x 30mm	10oz/Case
MER R24401.5	24 x 40mm	10oz/Case
MER R24501.5	24 x 50mm	10oz/Case
MER R24551.5	24 x 55mm	10oz/Case
MER R24601.5	24 x 60mm	10oz/Case
MER R25251.5	25 x 25mm	10oz/Case
MER R25301.5	25 x 30mm	10oz/Case


ITEM #	DESCRIPTION	QTY	TISSUE SOURCE	TISSUE TYPE
HIS 110	Acid Fast Bacteria	A=10, B=20, C=100	Mouse	Lung
AME CSA0725P	Adeno Carcinoma, QC w/PSA Antibody	25/Pack	Human	Prostate
AME CSA0925P	Adrenal Gland	25/Pack	Human	Adrenal Gland
HIS 120	Alcian Blue	A=10, B=20, C=100	Human	Umbilical Cord
HIS 1055	Alpha Feto Protein	A=10, B=20, C=100	Human	Liver
AME CSA0225P	Alzheimer's	25/Pack		
AME CSA0325P	Amoeba	25/Pack	Human	Intestine
HIS 130	Amyloid	A=10, B=20, C=100	Human	Heart
AME CSA0425P	Amyloid Alpha	25/Pack		
AME CSP0225P	Anterior Pituitary	25/Pack		
AME CSS0625P	Antrum Stomach	25/Pack	Human	Stomach
AME CSA0825P	Aorta	25/Pack	Human	Aorta
AME CSA1025P	Appendix	25/Pack	Human	Appendix
HIS 135	Argentaffin	B=20, C=100	Human	Intestine
AME CSA0525P	Argent-Argyrophil	25/Pack		
HIS 140	Argyrophil	A=10, B=20, C=100	Human	Tumor
AME CSA0625P	Aspergillus	25/Pack	Human	Lung
HIS 160	Bile	A=10, B=20, C=100	Human	Liver
AME CSB0325P	Bone Marrow	25/Pack	Human	Marrow
AME CSB0425P	Bladder	25/Pack	Human	Bladder
AME CSB0525P	Breast, EGFR	25/Pack	Human	Breast
AME CSB0625P	Breast Carcinoma, PR+ Controls (QC)	25/Pack	Human	Breast
AME CSB0725P	Breast Carcinoma, c-erb-2+	25/Pack	Human	Breast
AME CSB0825P	Breast Carcinoma, QC w/Ki67 Antibody	25/Pack	Human	Breast
AME CSB0925P	Breast Carcinoma, QC w/ER Antibody	25/Pack	Human	Breast
AME CSC0125P	Calcium	25/Pack		
AME CSC0225P	Cat Scratch Fever	25/Pack		
HIS 1103	CD3	A=10, B=20, C=100		
HIS 1105	CD5	A=10, B=20, C=100		
HIS 1108	CD8	A=10, B=20, C=100		
HIS 1115	CD15 (LEU-M1)	A=10, B=20, C=100		
HIS 1120	CD20 (L26)	A=10, B=20, C=100		
HIS 1130	CD30 (KI-1)	A=10, B=20, C=100		
HIS 1143	CD43	A=10, B=20, C=100		
AME CSC0725P	Colon Carcinoma	25/Pack	Human	Colon
AME CSC1025P	Colon	25/Pack	Human	Colon
AME CSC0325P	Copper	25/Pack	Human	Colon
AME CSC0425P	Cryptococcus	25/Pack		
HIS 1235	Cytokeratin-Monoclonal	A=10, B=20, C=100	Human	Skin
HIS 1236	Cytokeratin-Monocional	A=10, B=20, C=100	Human	Skin
HIS 1240	Cytomegalovirus	A=10, B=20, C=100 A=10, B=20, C=100	Human	Lung
HIS 1255	Desmin	A=10, B=20, C=100	Human	Leiomyoma
HIS 210	Elastic	A=10, B=20, C=100	Human	Skin
AME CSE0125P	Epstein-Barr Virus	25/Pack	Tiuliidii	JKIII
HIS 1285	Epithelial Membrane Antigen	A=10, B=20, C=100	Human	Kidney
HIS 1285	Factor 8	A=10, B=20, C=100	Human	Lymphoid
AME CSF0225P		25/Pack	Hulliali	Lymphold
	Ferric Iron  Fungus (Histopath)		Mouso	Lung
HIS 230	Fungus (Histopath)	A=10, B=20, C=100	Mouse	Lung
HIS 310	Gram	A=10, B=20, C=100	Mouse	Lung
HIS 1310	Glucagon	A=10, B=20, C=100	Human	Pancrease
AME CSG0125P	Glycogen	25/Pack		
AME CSH0825P	Hairy Cell Leukemia	25/Pack		
AME CSH0725P	Heart	25/Pack	Human	Heart
HIS 375	Helicobacter Pylori	A=10, B=20, C=100	Mouse	Intestine
HIS 376	Human Helicobacter	A=10, B=20, C=100		
HIS 380	Hepatitis	A=10, B=20, C=100	Human	Liver
HIS 1360	Hepatitis B, Core Antigen	B=20, C=100	Human	Liver














ITEM #	DESCRIPTION	QTY	TISSUE SOURCE	TISSUE TYPE
HIS 1365	Hepatitis Surface Antigen	A=10, B=20, C=100	Human	Liver
AME CSH0925P	Hepatoma, QC w/ AFP Antibody	25/Pack	Human	Liver
HIS 1380	Herpes Simplex Virus, Type 1	A=10, B=20, C=100	Human	Lung
HIS 1385	Herpes Simplex Virus, Type 2	A=10, B=20, C=100	Human	Lung
HIS 1390	Herpes Simplex Virus, Type 1 & 2	A=10, B=20, C=100	Human	Lung
HIS 1370	HMB-45	A=10, B=20, C=100	Human	Skin
AME CSH1025P	Hodgkin's Lymphoma	25/Pack		
HIS 1420	Insulin	A=10, B=20, C=100	Human	Pancrease
HIS 410	Iron	A=10, B=20, C=100	Human	Liver
HIS 1500	Kappa	A=10, B=20, C=100	Human	Lymphoid
AME CSK0225P	Kaposi Sarcoma	25/Pack		
AME CSK0125P	Kidney	25/Pack	Human	Kidney
HIS 1540	Lambda	A=10, B=20, C=100	Human	Lymphoid
HIS 1545	LCA Tonsil	A=10, B=20, C=100	Human	Lymphoid
AME CSL0125P	Legionnaires	25/Pack		
AME CSL0325P	Lipofuscin Pigment	25/Pack		
AME CSM0625P	Mantle Cell Lymphoma	25/Pack		
AME CSM0125P	Mast Cell	25/Pack		
HIS 530	Melanin	A=10, C=100	Human	Skin
AME CSM0225P	Melanoma HMB-45	25/Pack	Human	Skin
AME CSM0525P	Mesothelioma	25/Pack		
HIS 550	Mucin	A=10, B=20, C=100	Human	Intestine
AME CSC0625P	Myelin (Central Nervous System)	25/Pack		
AME CSN0325P	Neuroblastoma	25/Pack		
HIS 1610	Neuron Specific Enolase	A=10, B=20, C=100	Human	Pancrease
AME CSN0225P	Nodular Lymphoma	25/Pack		
AME CSO0125P	Ovarian Cancer	25/Pack	Human	Ovary
HIS 1640	P53	A=10, B=20, C=100	Human	Colon
HIS 610	Periodic Acid Schiff's (PAS)	A=10, B=20, C=100	Human	Normal Liver
HIS 620	Periodic Acid Schiff's w/Digestion	A=10, B=20, C=100	Human	Normal Liver
AME CSP0325P	Placenta	25/Pack	Human	Placenta
AME CSP0125P	Pneumocystis	25/Pack	Human	Lung
AME CSP0425P	Prostate, Normal	25/Pack	Human	Prostate
HIS 665	Ptah	A=10, B=20, C=100	Human	Muscle
AME CSR0225P	Renal Cell Carcinoma	25/Pack		
HIS 720	Reticulum	A=10, B=20, C=100	Human	Liver
HIS 1740	S-100	A=10, B=20, C=100	Human	Nerve
AME CSS0725P	Skeletal Muscle	25/Pack		
AME CSS0325P	Small Intestine	25/Pack	Human	Intestine
AME CSS0425P	Smooth Muscle	25/Pack		
HIS 1760	Somatostatin	A=10, B=20, C=100	Human	Pancrease
HIS 755	Spirochete	B=20, C=100	Rabbit	Testis
SCY TCS002025	Spirochete	25/Pack		
AME CSS0525P	Spleen	25/Pack	Human	Spleen
AME CSS0225P	Staphylococcus, Inoculated	25/Pack		
HIS 1770	Synaptophysin	A=10, B=20, C=100	Human	Pancrease
HIS 1780	T-Cell CD45RO	A=10, B=20, C=100	Human	Lymphoid
AME CST0425P	Testicular Cancer	25/Pack	Human	Testical
AME CST0725P	Thymus	25/Pack		
AME CST0525P	Thyroid Adenoma	25/Pack	Human	Thyroid
HIS 790	Trichrome	A=10, B=20, C=100	Human	Normal Liver
AME CST0225P	Tonsil	25/Pack	Human	Tonsil
AME CST0325P	Toxoplasma Gondii	25/Pack		
AME CST0625P	TTF-1, Thyroid	25/Pack	Human	Thyroid
AME CSU0225P	Uric Acid	25/Pack		
AME CSU0325P	Uterus	25/Pack	Human	Uterus
HIS 1870	Vimentin	A=10, B=20, C=100	Human	Tumor





### **UNIMAILER**<sup>™</sup>

### SINGLE PLASTIC SLIDE MAILERS

- · Accommodates slides with or without cover glass
- · Secured stacking by an innovative locking mechanism
- · Attach an ID label or tamper-evident tie through holes
- · Avoids any vibrations or breakage during transport
- · Available in pastel colors for easy identification

SIM M800100B	Blue	200/Case
SIM M800100G	Green	200/Case
SIM M800100P	Pink	200/Case
SIM M800100W	White	200/Case
SIM M800100Y	Yellow	200/Case

# SIM T345

### **LOCKMAILER™**

### MAGENTA MICROSCOPE SLIDE JARS

- Tamper evident multi purpose container for mailing, staining or storing microscope slides
- · Constructed of extra-strong, clear polypropylene
- · Holds up to 4 standard (1" x 3") slides vertically
- · Inside channels slotted to keep slides safely separated
- Unique tamper evident leakproof screw cap prevents someone from manipulating your slides during transport or storage situations without your knowledge
- · Can use without the tamperproof locking mechanism
- Perfect also for slide conveyors and specimen slide transport between the doctor's office and the lab
- Designed for maximum stability on a bench top while having an internal volume of only 12mL
- CapInsert™ (SIM T345 Series) may be inserted on top of closure for color coding purposes (10 Colors Available)

SIM M9504MA	Holds 4 Slides	500/Case

### **CAPINSERT™**

The Capinsert™ is used to color-code a multitude of Simport® products according to your specific needs. It is inserted on top of the closure and has a write-on frosted area for sample identification.

- · Made of polypropylene
- · Write-on frosted area for sample identification
- COLORS AVAILABLE: Blue, Gray, Green, Lilac, Orange, Pink, Red, Violet, White, Yellow, Assorted

SIM T345 CapInsert<sup>™</sup> 500/Case



### 2-PLACE PLASTIC SLIDE MAILERS

- · Ideal for secure shipping and storing of slides
- Polypropylene protectors hold 2 standard (1" x 3") slides
- · Large, easy-to-secure snaps for safe storage
- Array of colors brings some fun to the laboratory and is useful for color coding
- ASSORTED COLORS: Red, Yellow, Blue, Purple, Orange

HEA HD15983A	Assorted	25/Pack
HEA HD15983B	Natural	25/Pack

### **CLEAR PLASTIC SLIDE MAILERS**

- Flat plastic slide mailers produced from polypropylene
- Snap-close lid keeps slides secure during transport

GLO 513031	Single	100/Box
GLO 513041	Double	100/Box
GLO 513045	Triple	100/Box

### **SCREW CAP SLIDE MAILERS**

- · Inner slots keep slides separate
- · Leak tight red screw cap secures slides

GLO 513068 Quadruple 500/Case



### **5-PLACE PLASTIC SLIDE MAILERS**

- · Ideal for secure shipping, storing, or staining slides
- Polypropylene containers hold 5 standard (1" x 3") slides
- Inside wall features slotted channels to keep slides separated and to make removal easier
- · Lid closes securely to keep slides safe
- Side-Opening Mailers measure 81 x 19 x 31mm
- End Opening Mailers measure 43 x 24 x 88mm
- ASSORTED COLORS: Red, Yellow, Blue, Purple, Orange


HEA HD15986	End Open, Natural	25/Pack
HEA HD15982	Side Open, Natural	25/Pack
HEA HD15983G	End Open, Assorted	25/Pack

### **CARDBOARD SLIDE MAILERS**

- Securely hold slides with or without a cover-glass
- Features an inner finger groove for effortless slide removal along with a cover that opens for viewing
- Flat surface area is ideal for labeling or writing any necessary information
- · Lightweight and reusable to minimize cost

GLO 513001	Single	100/Box
GLO 513002	Double	100/Box
GLO 513003	Triple	50/Box
GLO 513004	Quadruple	50/Box













### **CARDBOARD SLIDE FOLDERS**

- Convenient solution for safely working, transporting and shortterm storing of slides
- Easily organize slides with the index on the front cover that corresponds with the number grid inside
- Recessed horizontal compartments and inside covers protect specimens when closed to reduce contamination
- · Front protection flaps attached by tape for repeated use
- Thumb cut in slide compartments enables easy removal of slides with slight pressure at one end
- Fold flaps flat under the slide holder to save space
- Multiple colors for increased coding and organization

MER MSL26	20-Place, Blue	Each
MER MSL27	20-Place, Green	Each
MER MSL28	20-Place, Red	Each
MER MSL29	20-Place, Yellow	Each
CAS 9000	20-Place, Black	Each
HEA HD9930	30-Place, Black	Each



### **20-PLACE SLIDE TRAY RACKS**

- Holds up to (5) 20-place slide folders (NOT INCLUDED)
- Ideal for microbiology, histology and pathology fields
- · Eliminates the worry of breakage during retrieval
- Matching clasps provide area for identification labels and doubles as a handle (EACH RACK INCLUDES 5 CLASPS)
- Portable and refrigerator safe for cold sample storage
- Stack multiple units to create your very own library
- **Dimensions:** 226mm x 340mm x 125mm

**HEA HD99040A** 

Slide Tray Rack, Blue

Each







### **VIEW-PACK™ SLIDE HOLDER**

- Standard (9" x 12") 3-ring binder with ten (8½" x 10½") vinyl microscope slide holder pages (160 slides total)
- Each page holds (16) 1" x 3" slides in individual pockets with a white back and clear front
- · Center flap holds slides securely
- Fold the page at the center crease to retrieve slides

BEL 441700000	Slide Holder Binder	Each
BEL 441710000	Slide Holder Pages	10/Pack

### **SLIDE HOLDERS**

- · Hold slides during evaluation, cleaning and drying
- · Made from chemical-resistant polyethylene
- · Solid, strong construction with rubber feet for stability

HEA HD15981A	Holds 40 Slides	Each
HEA HD15981B	Holds 78 Slides	Each
HEA HD15981C	Holds 15 Slides	Each



### **20-PLACE PLASTIC SLIDE FOLDERS**

- Base holds slides horizontally with numbered spaces
- Transparent hinged covers allows easy reading of ID labels with or without an optical bar code reader
- Removal of slides is made easy simply by pressing on one end, which will automatically lift the other end
- Stackable units take minimum space on tables or shelves
- Resist temperatures between -80°C and +80°C

SIM M75020B	Blue	10/Case
SIM M75020G	Green	10/Case
SIM M75020GY	Gray	10/Case
SIM M75020W	White	10/Case
SIM M75020Y	Yellow	10/Case
SIM M75020AS	2 of Each Color	10/Case

### SPONGE BLOCK SPACERS

- · Made of high-quality neoprene foam rubber
- · Keep slides standing upright when the drawer isn't full
- Compresses against the walls and holds itself and the slides in place due to it being wider than the drawer

**BLK 901-0007** 

Sponge Blocks, Black

12/Pack





# TANNER SCIENTIFIC® 6-DRAWER SLIDE CABINET

THE PERFECT SOLUTION TO A LABORATORY'S MICROSCOPE SLIDES ORGANIZATIONAL NEEDS

- Holds up to 5,000 (1" x 3") slides
- Solid steel construction and chemical resistant finish ensure many years of safe and secure slide storage
- Features black pull knob drawers and metal dividers
- · Stackable up to 10 cabinets high
- · All six drawers slide smoothly and easily
- BASE OPTION AVAILABLE
- **Dimensions:** 153/4"W x 19"D x 5"H

TNR 2000	6-Drawer Cabinet	Each
TNR 2001	Base for TNR 2000	Each

## EBERBACH® 2800 SERIES MICROSCOPE SLIDE CABINETS

- Safe, short-term and long-term storage of prepared slides can easily be indexed for retrieval
- · Hinged doors ensure long operation and a close fit
- Sturdy aluminum slide trays, numbered in sets of 25, are removable and interchangeable between cabinet systems
- Individual compartments for laying slides flat, maintain cover glasses in place, preserving your work year round
- Nominal capacity can be doubled, as two micro slides can be stacked in each compartment
- Flat sides and disappearing door of the unit system provide for easy stacking and/or banking
- Cabinet system will initially handle a larger number of slides with drawer space at the bottom for a card index
- · Replacement slide trays sold separately

EBE 2800	Holds 500 Slides	Each
EBE 2820	Holds 1,000 Slides	Each
EBE 2840	Holds 3,000 Slides	Each
EBE 2870	Metal Slide Tray	Each





### **CORK LINED, PLASTIC SLIDE BOXES**

- · Made from sturdy ABS plastic and designed for stacking
- · Index paper on inside lid for documented slides
- Cork lined bottom with embossed numerals on bottom
- Self locking lid (25-Place), hinged lid (50-Place), and hinge and rust-resistant clasp (100-Place) close securely
- SOLD INDIVIDUALLY

COLOR	25-PLACE	50-PLACE	100-PLACE
White	CAS 970025	CAS 970055	CAS 970125
Blue	CAS 970026	CAS 970056	CAS 970126
Green	CAS 970027	CAS 970057	CAS 970127
Red	CAS 970028	CAS 970058	CAS 970128
Yellow	CAS 970029	CAS 970059	CAS 970129



### **TALL PREMIUM SLIDE BOX**

- Sturdy ABS material for rugged use
- · Grooved slots hold individual slides in place
- · Soft cork lining ensures slides will not slip
- · Numbered identification card, lining and grooved slots
- Additional space for larger (2" x 3") microscope slides

HEA HD15979A Blue, Holds 100 Slides Each

### TRUE NORTH® SLIDE BOX

- Durable polycarbonate and stainless steel clasp
- Store in ultra-low temperature ranges down to -80°C
- Uniform grooved slots separate slides to stop contact
- Includes temperature-resistant numbered foam lining and identification sheet

HEA HD15988H Purple, Holds 100 Slides Each

### STANDARD SLIDE BOXES

- · Constructed of durable polypropylene with no liner
- · Hinged lids held securely closed by latch lock
- · Stacking feet maximize stability
- · Includes identification card

HEA HD15990A	Blue, Holds 25 Slides	Each
HEA HD15990B	Red, Holds 25 Slides	Each
HEA HD15990C	White, Holds 25 Slides	Each
HEA HD15991A	Blue, Holds 100 Slides	Each
HEA HD15991B	Red, Holds 100 Slides	Each
HEA HD15991C	White, Holds 100 Slides	Each



# TANNER SCIENTIFIC® SLIDE STORAGE SYSTEM

• Designed for the storage of microscope slides

• Includes: Dividers, (4) Drawers

• Capacity: 2,000 Slides

• Dimensions: 12"  $\times$  12"  $\times$  4"

• Shipped flat, easy to assemble

• MAY ALSO BE USED TO STORE CASSETTES

TNR SLIDBX20 Slide Storage System Each



### MICROSCOPE SLIDE STORAGE UNITS

- · Holds approximately 2,200 slides per unit
- Holds approximately 13,200 slides per case of 6 units
- SHIPPED FULLY ASSEMBLED

### **REGULAR**


Includes: (6) High-Impact Polystyrene Plastic Trays,(6) Sturdy Corrugated Drawers with Front Labeling, (6) Corrugated Drawer Enclosures for Dust-Proof Protection, (1) Outer Carton for Long Term Storage

**LSS MS11** 81/4"W x 141/2"D x 41/2"H 6/Case

### ALL CORRUGATED

Includes: (6) Corrugated Drawers with Dividers 11/8" Apart, (6) Corrugated Drawer Enclosures for Dust-Proof Protection, (1) Outer Carton for Long Term Storage

**LSS ACS1** 81/4"W x 141/2"D x 41/2"H 6/Case



### **SLIDE INDEX MARKERS**

- Very unique design makes process easy to load slides
- · Contains spaces for case number, initials and date
- Four vertical spaces available for other information
- Easy color coding provides easy slide retrieval
- Dimensions: 1"W x 3½"H

MER SIMBLANK	Blank	1,000/Box
MER SIMYELLOW	Yellow	1,000/Box
MER SIMPINK	Pink	1,000/Box
MER SIMGREEN	Green	1,000/Box
MER SIMBLUE	Blue	1,000/Box
MER SIMLILAC	Lilac	1,000/Box
MER SIMWHITE	White	1,000/Box









### **SLIDEFILE**<sup>™</sup>

- Made of high impact polystyrene
- · Slides are stored upright for easier insertion and removal
- · Read bar codes without removing slides from tray
- Two index cards numbered 1 to 100 supplied per box
- Colors for slide classification and reduce sample mix-up
- Resists temperatures between -80°C and +80°C
- Simply remove the tray and line up 400 slides in 3 rows for long term storage and maximum storage space

SIM M700100B	Blue	10/Case
SIM M700100G	Green	10/Case
SIM M700100P	Pink	10/Case
SIM M700100W	White	10/Case
SIM M700100Y	Yellow	10/Case

### SLIDEFILE™ JR.

 Simply remove the tray and line up 200 slides in 3 rows for long term storage and maximum storage space

SIM M70050B	Blue	10/Case
SIM M70050G	Green	10/Case
SIM M70050P	Pink	10/Case
SIM M70050W	White	10/Case
SIM M70050Y	Yellow	10/Case

### **DRAINRACK**<sup>™</sup>

- Made of high impact polystyrene
- Rugged tray used as a drain rack can hold up to 200 microscope slides in 100 individual numbered slots
- Slides are stored upright for easier insertion and removal
- Read bar codes without removing slides from the box
- Easy-to-read numbers from 1 to 100 identifying each slot
- Simply tilt slide forward and backward with one finger to easily and rapidly pick up

SIM M710100B	Blue	10/Case
SIM M710100G	Green	10/Case
SIM M710100P	Pink	10/Case
SIM M710100W	White	10/Case
SIM M710100Y	Yellow	10/Case

### DRAINRACK™ JR.

 Can hold up to 100 microscope slides in 50 individual numbered slots for long term storage

SIM M71050B	Blue	10/Case
SIM M71050G	Green	10/Case
SIM M71050P	Pink	10/Case
SIM M71050W	White	10/Case
SIM M71050Y	Yellow	10/Case





MER 1042 Staining Dipper

# MERCEDES SCIENTIFIC® QUIK-DIP® STAIN

The Quik-Dip® Stain is a high-quality, quick-time staining kit that serves numerous purposes in today's cost-conscious lab setting. It can be used for fine needle aspirations, cytology, blood smears, bone marrow biopsies as well as bone marrow smears, touch preps, and can be used to detect H.Pylori. The kit contains a fixative for air-dried cell suspensions and/or touch prep slides.

### FREE SAMPLES AVAILABLE!

MER 1001	3 x Gallon Set	Set
MER 1002	3 x 32oz Set	Set
MER 1003	3 x 4oz Set	Set
MER 1001A	1GL Green Bottle	Each
MER 1001B	1GL Red Bottle	Each
MER 1001C	1GL Blue Bottle	Each
MER 1002A	32oz Green Bottle	Each
MER 1002B	32oz Red Bottle	Each
MER 1002C	32oz Blue Bottle	Each
MER 1003A	4oz Green Bottle	Each
MER 1003B	4oz Red Bottle	Each
MER 1003C	4oz Blue Bottle	Each

### **STAINING SETS**

- Suitable for histological, cytological and hematological staining procedures with total immersion for slides
- Metal staining racks made of stainless steel, lightweight for mobility

MER 1011	12-Place Metal Rack, Staining Dipper, (12) Green Stain Dishes with Lids	Set
MER 1012	12-Place Metal Rack, Staining Dipper, (12) Clear Stain Dishes with Lids	Set
MER 1021	3-Place Metal Rack, Staining Dipper, (3) Green Stain Dishes with Lids	Set

### STAINING DISHES

- Dimensions: 4"L x 2"W x 4"H; Capacity: 300mL
- · Lids reduce reagent evaporation or contamination
- · Resistant to acids and solvents

MER 1032	Clear, Plastic	Each
MER 1030	Green, Plastic	Each

### STAINING DIPPER


· Plastic handle attached for total immersion

MER 1042	Holds 25 Slides	Each



- Throughput: 50 samples/hr. Semi-automated
- **Dimensions:** 11"W x 15"D x 8.5"H; **Weight:** 31 lbs

The CellSolutions® double membrane technology yields cellular preparations with increased cellularity, presentation, and stain quality on 93.75% of the specimens tested when compared to other liquid based methods.1



# **FINALLY** Something New After 20 Years!

### **FASTER**

Slide is ready to stain in 20 seconds.

### CHEAPER

Less expensive to acquire and operate. Costs a fraction of the others.

### **BETTER**

Slides will be better looking and easier to read.

1 Smith D., Espinoza, A., Vega, I. "Comparison Study of Traditional Cytology Non-Gyn Processing Techniques to Specimens Processed with CellSolutions Method". 2017 American Society of CytoPathology Meeting: Poster Publication

The **CellSolutions**® preservatives and the **CellSolutions**® **F50** provide a high quality and efficient method for liquid-based nongyn cytology processing.

### We Have Enhanced Preservatives:

- Effectively lyses red blood cells and removes hemoglobin and proteins in solution
- · Maintains cellular morphology
- · Offers ability to perform ancillary molecular testing
- Classified as non-hazmat; no additional shipping charges, processing or specimen transport fees

### **Cytology Preservatives**

### **F50 GENERAL CYTOLOGY PRESERVATIVE**

CES CF101500	10mL Vial	500/Kit
CES CF101200	10mL Vial	200/Kit
CES CF10125	10mL Vial	25/Kit
CES CF10100	30mL Prefill; 50mL Conical	100/KIt
CES CF101L	1L Bottle	Each
CES CF101G	4 x 1L Bottles	4/Case

### RED LYTIC GENERAL CYTOLOGY PRESERVATIVE

CES CR10225	90mL Cup	25/Case
CES CR102L	1L Bottle	Each
CES CR102G	4 x 1L Bottles	4/Case

### **BLUE PRESERVATIVE (URINE)**

CES CB10225	90mL Cup	25/Case
CES CB102L	1L Bottle	Each
CES CB102G	4 x 1L Bottles	4/Case

### F50 Slide/Filter Kits


CES CS400FS	Small, 14mm	400/Kit
CES CS400E	Large 20mm	400/Kit



### **Red Lytic Preservative**

- Preservative fluid formulated for the preservation of non-cervicovaginal (non-gyn) cells in suspension
- Thin-layer cytology slides are processed from these cell suspensions using the CellSolutions® Automated and GluCyte® Manual Methods for cytology slide preparation
- Slide preparations are evaluated for the presence of cancer or its precursor lesions by cytotechnologists and pathologists trained to evaluate CellSolutions® prepared slides.































### M964: CYTOSEP™ FUNNELS

FOR SHANDON° CYTOSPIN° 4 CYTOCENTRIFUGE

SIM M96410FW	Single Funnel with White Filter & Cap	50/Pack
SIM M96410FT	Single Funnel with Tan Filter & Cap	50/Pack
SIM M9641	Single Funnel Only	500/Case
SIM M96420FW	Double Funnel with White Filter & Cap	500/Case
SIM M9641D	Double Funnel Only	500/Case

### M965: CYTOSEP™ FUNNELS

FOR SHANDON° CYTOSPIN° 4 CYTOCENTRIFUGE

SIM M96510FW	Single Funnel with White Filter & Cap	40/Pack
SIM M96510FT	Single Funnel with Tan Filter & Cap	40/Pack
SIM M9651	Single Funnel Only	50/Pack
SIM M96520FW	Double Funnel with White Filter & Cap	480/Case
SIM M9651D	Double Funnel Only	50/Pack

### M964 & M965 FUNNEL ACCESSORIES

SIM M965B	Base Only	50/Pack
SIM M965C	Cap Only	50/Pack
SIM M965FW	White Filter Paper for Single Funnel	200/Case
SIM M965FT	Tan Filter Paper for Single Funnel	200/Case
SIM M965FWD	White Filter Paper for Double Funnel	200/Case

### M963: CYTOSEP™ FUNNELS

FOR SAKURA® CYTO-TEK® CYTOCENTRIFUGE

SIM M9631	Fluid Chamber Only, 1mL	200/Case
SIM M963B	Base Holder Only	200/Case
SIM M963C	Cap Only	200/Case
SIM M963FW	White Filter Card Only	200/Case

### M967: CYTOSEP™ FUNNELS

FOR ELITECH CYTOPRO® CYTOCENTRIFUGE

SIM M96710FW	Single Chamber with White Filter Paper & Cap	48/Pack
SIM M967FW	White Filter Paper for Single Funnel	100/Pack
SIM M96720FW	Double Chamber with White Filter Paper & Cap	48/Pack
SIM M967FWD	White Filter Paper for Double Funnel	100/Pack

## CYTOSEP™ CYTOLOGY FUNNEL CLIPS

- Reusable and autoclavable
- Holds funnel disposable sample chamber against slides
- Metal will not rust or corrode during decontamination

SIM M964B Stainless Steel Clips



















6/Case

### PAPANICOLAOU STAIN SOLUTIONS

- Imparts a characteristic range of coloration to aid in the identification and classification of exfoliative cells of vaginal, cervical, prostatic and other body secretions
- · Alcohol-based counterstains used in cytology for routine cytoplasmic staining using the Papanicolaou Method

### **OG-6 CYTOLOGY STAIN**


• Used to detect the presence of keratin in the cytoplasm

EKI 8037500ML	16oz Bottle	Each
EKI 80371GL	1GL Bottle	Each

### **EA-50 CYTOLOGY STAIN**

· Contains dye strengths optimized for traditional use

EKI 8033500ML	16oz Bottle	Each
EKI 80331GL	1GL Bottle	Each





PAPANICOLAOU STAIN

ANX 9300101P

# THY 36827G

### THY 36825G

THY 36848G

### **CYTOLOGY BRUSHES**

THY 36825G	Rovers® Cervex-Brush®	25/Pack
THY 36827G	Rovers® Cervex-Brush® Combi	25/Pack
THY 36840G	Rovers® EndoCervex-Brush®	50/Pack
THY 36848G	Rovers® EndoCervex-Brush®-S	50/Pack



### **CERVICAL SCRAPERS**

- · Sterile, wooden and naturally biodegradable
- · Ayres style with splinter free edges
- · Cervical specimen collection for pap smears

SOL 59200	7", Non-Sterile	250/Box
SOL 59600	7", Sterile, Ayres	100/Pack

### SLIDE-RITE® FIXATIVE SPRAY

- · Water soluble fixative for exfoliative cytology
- · Forms a clear, protective coating over specimens and cytological smears and fixates cells immediately
- Protects samples during transportation and in storage
- Dries in as little as 5 minutes

FIS 0435588 2oz Spray Bottle Each

### SAFETEX® FIXATIVE SPRAY

- Water and alcohol soluble polyethylene glycol formulation for affixing specimens
- Tough thin film for maximizing cell preservation and ensuring optimum staining environment
- · No hazardous polymerization will occur

ANX 930010IP	1oz Spray Bottle	12/Case
ANX 930022	4oz Spray bottle	12/Case





### NO-TOUCH® PAP SMEAR KITS

- · Provides maximum slide specimen protection
- · Ventilation folds allow wet slides to dry in courier pack

ANX 230080	Slide, Brush, Scraper, 1oz Fixative Spray	500/Case
ANX 230100	Slide, Brush, Scraper, Swab, Fixative Pouch	500/Case
ANX 230180	Slide, Brush, Scraper	500/Case



### SAMPLE STORAGE TRAYS

- Sturdy corrugated trays hold 25 samples
- · Preprinted on all sides for vial identification
- · Stackable and easy to assemble
- · Packed flat for convenient storage

MKL ML3488	SurePath™ Storage Tray, 8½"L x 8½"W x 2¼"H	10/Pack
MKL ML3487	ThinPrep® Storage Tray, 101/8"L x 101/8"W x 3"H	10/Pack

### PAP-FIX PLUS™ PAP JARS

- · Outstanding performance in preserving and transporting pap smear slides
- Denatured 95% ethyl alcohol preserves pap smear samples
- Rigid polypropylene construction provides fortress-like protection for up to four slides in transit



### PREFILLED PAP JARS (ETHYL ALCOHOL)

ANX 3085500G8K	Green Caps	100/Pack
ANX 3085500B8K	Blue Caps	100/Pack
ANX 3085500K8K	Fuchsia Caps	100/Pack
ANX 3085500L8K	Lavender Caps	100/Pack
ANX 308550008K	Natural Caps	100/Pack
ANX 3085500E8K	Orange Caps	100/Pack
ANX 3085500R8K	Red Caps	100/Pack
ANX 3085500Y8K	Yellow Caps	100/Pack

### **EMPTY PAP JARS, BULK PACK**

ANX 2405400G8K	Green Caps	1,000/Pack
ANX 2405400B8K	Blue Caps	1,000/Pack
ANX 2405400K8K	Fuchsia Caps	1,000/Pack
ANX 2405400L8K	Lavender Caps	1,000/Pack
ANX 240540008K	Natural Caps	1,000/Pack
ANX 2405400E8K	Orange Caps	1,000/Pack
ANX 2405400R8K	Red Caps	1,000/Pack
ANX 2405400Y8K	Yellow Caps	1,000/Pack
ANX 2405400X8K	Assorted Caps	1,000/Pack



### **CENTRIFUGE TUBES**

- · Manufactured and packaged using an improved process, while maintaining standards
- · Graduations are printed in bold, easy-toread black ink and are accurate to +/-2%
- Large white marking area provides excellent contrast for identification
- Easy to write on flat surface polyethylene (HDPE) caps feature integrated leak-proof seal
- · Tubes include graduations molded into the conical bottom, ideal for pellet viewing and small volume quantification
- · Convenient zip-closure bag is a practical alternative to folding and taping opened bags

Easily refill reusable racks with tubes packed in bags

CTS 229413	15mL, 25/Bag, Sterile	500/Case
CTS 229435	50mL, 25/Bag, Sterile	500/Case

### FOAM RACK, STERILE

Foam rack is a sturdy alternative to disposable paperboard racks

CTS 229412	15mL, 50/Rack, Sterile	300/Case
CTS 229422	50mL, 50/Rack, Sterile	300/Case









# MERCEDES SCIENTIFIC® ALCOHOLS & SOLVENTS

Mercedes Scientific offers a complete product line of alcohols and solvents that are ideal for tissue processing, and all histology and cytology staining procedures or applications in the laboratory environment. Review all warnings, safety measures and instructions before use.

MER 1200	100% Reagent Alcohol, 1GL	Each
MER 1210	95% Reagent Alcohol, 1GL	Each
MER 1202	80% Reagent Alcohol, 1GL	Each
MER 1203	70% Reagent Alcohol, 1GL	Each
MER 1220	Isopropyl Alcohol 99%, 1GL	Each

**ALCOHOLS:** pure filtered and laboratory tested to ensure a consistent high quality for all formulating, processing and staining procedures in histology and cytology

**METHANOL:** for general tissue processing as a substitute for ethanol or as a cytology fixative

**ACETONE:** use as a mounting medium remover, mordant or cleaning solvent; for hardening and dehydrating tissues

XYLENE: highly refined to deliver consistent performance, high solvency ensures rapid and complete de-waxing and clearing; for optimal tissue infiltration and deparaffinization

MER 1230	Methanol, 1GL	Each
MER 1240	Acetone, 1GL	Each
MER 1250	Xylenes, 1GL	Each

# XYLENE BASED



### **MERCEDES SCIENTIFIC® OPTIC·MOUNT I™ (XYLENE)**

- · Antioxidant included to prevent the fading of standard hematoxylin and eosin stains
- · Miscible with all routine clearing reagents
- · Optimum refractive index for viewing and photo microscopy, and will not yellow with age

MER 7722	4oz Bottle	Each
MER 7723	500mL Bottle	Each

### SHUR/MOUNT®

- · Xylene based routine coverslipping medium for histology and cytology slides
- · Acrylic medium provides superior optical clarity and has a refractive index of 1.4
- · Recommended for automated coverslipping instruments that use both glass and tape coverslips
- · Contains an antioxidant to prevent fading and an UV absorbent to prevent yellowing with long term storage
- · Low-viscosity coverslipping medium dries quickly for viewing, storage and filing

TRX LCA 16oz Bottle 4/Case

### CYTOSEAL® XYL

- · Xylene-based rapid drying allows microscopic examination of slides soon after application
- · Inhibit the fading or yellowing of stained specimens

FIS 83124 4oz Bottle Each

### **BALSAM BOTTLE**

- · Ideal for applying media used in mounting slides
- Soda-lime glass bottle is sealed by a ground-glass cap fitted to ground-glass neck of the bottle
- Includes: Glass Applicator Rod, Glass Cap

FIS 030821 100mL Bottle 6/Case

### SHANDON CONSUL-MOUNT™

- · Permanent, organic, xylene-based mountant
- · Clear liquid with medium viscosity and refractive index
- For Thermo Scientific® Shandon Consul™ Coverslipper

FIS 9990440	Histology Media, 4oz	4/Case
FIS 9990441	Cytology Media, 4oz	4/Case

# TOLUENE BASED



## MERCEDES SCIENTIFIC° OPTIC·MOUNT I™ (TOLUENE)

- Mounting medium designed to be used with staining procedures using xylene as a clearing agent
- · Miscible with all routine clearing reagents
- Formulated with antioxidant to prevent fading of standard hematoxylin and eosin stains
- Optimum refractive index for photo microscopy

MER 7720	4oz Bottle	Each
MER 7724	500mL Bottle	Each

#### CYTOSEAL® 60

- Low viscosity medium dries quickly and allows for rapid, even spreading to virtually eliminate air bubbles
- Prevent cracking, fading and discoloration of specimens

FIS 83104	4oz Bottle	Each
FIS 831016	16oz Bottle	Each



#### **PERMOUNT®**

- · Gives a thinner adhesive layer for better optical quality
- Refractive index near that of fixed protein keeps image free of optical imperfections such as diffraction halos
- Suitable for mounting cover slips to slides with thick or thin specimens for long-term storage
- Won't shift the colors of biological dyes, or lose adhesive quality, cohesiveness or color with age
- Preserves most biological dyes with little or no fading when slides are stored in darkness

FIS SP15100	100mL Bottle	Each
FIS SP15500	500mL Bottle	Each

#### **SHANDON-MOUNT™**

- Performs quick, high-quality conventional coverslipping
- · Overnight-drying liquid cover slip with low viscosity

FIS 1900333	4oz Bottle	Each
FIS 9990435	500mL Bottle	2/Case

#### RICHARD-ALLAN® MOUNTING MEDIA

- · Rapidly dries and optically clear
- · Store slides for prolonged periods without stains fading
- Miscible with clearing agents including xylene substitutes

FIS 22110610	4oz Bottle	Each
FIS 22050102	1pt Bottle	2/Case

Mercedes Scientific

## WATER BASED













SIG F4680

AME MMC0226

AME MMC0243

#### **AQUASLIP™ AQUEOUS**

- Aqueous mounting medium prevents stained tissue chromogens from fading, while allowing permanent storage of stained slides
- · Doesn't require coverslipping
- Compatible with chromogens such as AEC, DAB, Fast Red, BCIP/NBT, and BCIP/INT

AME MMC0619 30mL Bottle Each

#### **CC/MOUNT™**

- · Aqueous mounting medium with high refractive index
- Permanent preservation of tissue sections stained with peroxidase and alkaline phosphatase based systems as well as with various fluorescent dyes
- Permits the long-term storage of antigens localized using organic insoluble chromogen substrates

SIG C9368 30mL Bottle Each

#### **AQUA-MOUNT™**

- Formulated for immunostaining and other techniques
- Non-fluorescing medium with low viscosity for maximum specimen absorption
- · Resists fading and polymerizes quickly for slide handling

**THS 41799008** 20mL Bottle Each

#### **FLUOROMOUNT™ AQUEOUS**

- Preserves fluorescence labeling in tissue sections
- Aqueous based mounting medium for cover slipping tissues stained with various fluorescent dyes

SIG F4680 25mL Bottle Each

#### SHUR/MOUNT®

- Water soluble, non-fluorescing formulation for mounting sections from aqueous solutions
- Use for frozen sections and special procedures that will not tolerate alcohols or solvents
- · Coverslip directly from water or buffer using standard methods

TRX LCW 16oz Bottle 4/Case

#### **COVERSAFE™**

- Fast-drying, non-toxic, low-hazard and optically clear
- · Works well with all clearing agents
- Non-volatility prevents problems with bubbles, solvent evaporation spaces or mounting medium flow
- · Retains excellent refractivity indefinitely
- · Used in automatic coverslipping machines
- Does **NOT** contain xylene, toluene, dioxane or benzene

AME MMC0243	1PT Bottle	Each
AME MMC0226	100mL Bottle	Each

#### **MERCEDES SCIENTIFIC® BOND-IT ADHESIVE**

- · Blended solution promotes a strong adhesive bond between tissue sections and microscopic slides
- · Use as an additive to waterbaths or as a slide coating application to reduce the loss of tissue sections
- · Suited to all immunohistochemical procedures

MER 3225500ML	16oz Bottle	Each
MER 32251GL	1GL Bottle	Each

#### MERCEDES SCIENTIFIC® 10% FORMALIN FIXATIVE (NBF)

- · Neutral buffered formalin solution is used in proper tissue fixation for accurate histopathologic evaluation
- · No dilution or addition of other agents is necessary
- · Penetrates quickly, but fixes slowly

MER 44991GL	1GL Bottle	Each
MER 449920L	20L Cube	Each

#### **FORMALIN WITH EOSIN TINT**

- · Intended for use as a general histology fixative
- · Addition of eosin to the formulation allows more easy visuals of small fragments of tissue during the transfer of the tissue from the fixative to the tissue cassette
- · AVAILABLE SIZES: 1GL, 2.5GL Cube, 5GL Cube

**CAN FX10** Formalin w/Fosin Tint Fach

#### **MERCEDES SCIENTIFIC®**

#### LIQUIMOUNT™ (ALCOHOL BASED)

- · Simultaneously clears and mounts coverslips in one step
- · Eliminate the need to use hazardous xylene as a clearing agent with a xylene-free formulation
- No more final step required of clearing and coverslip directly from absolute isopropyl alcohol
- Suited to manual and automated coverslipping processes

**MER 773440Z** 4oz Bottle Fach

#### FORM-C™ FORMALIN

- · Pre-buffered formalin concentrate
- MIX IN 30 SECONDS: simply add 1GL of Form-C to 4GL of water to make 5GL of 10% neutral buffered formalin
- · High purity, true neutral pH and non-clogging
- · Reduces tissue processor maintenance
- · Stable in freezing temperatures

**ABS FORM04** 1GL Bottle 4/Case

#### GROSSING-AID™

- · Alcoholic formalin fixative with enhancing additives renders tissues with high DNA content white or a lighter color than the surrounding tissues
- Rapidly fixes adipose tissue and connective tissues
- · Causes lymph nodes to appear white against a darker background of connective tissues and adipose tissues

**BBC 0280** 1GL Rottle Each



MER 3325500ML











ABS FORM04

BBC 0280

## MERCEDES SCIENTIFIC° XYLENE SUBSTITUTE

- · Nonhazardous, nonflammable, odorless with low toxicity
- Clearing and deparaffinization agent for processing and staining histological or cytological specimens

MER 3310 1GL Bottle Each

#### **AMERICLEAR™**

- Citrus clearing solvent and non-flammable solution
- · Ideal for processing, staining and removing paraffin
- · Replaces xylene for automated or manual systems
- · Low-toxicity, all-organic compound
- Higher flash point that removes the necessity of being stored in a flammable safety cabinet

ANX 72060044 1GL Bottle Each

#### **CITRA-CLEAR™ XYLENE SUBSTITUTE**

- 100% d-Limonene for processing and staining
- · Mild citrus odor
- · Exceptionally high purity
- Flash point 122°F

BBC 64001GL 1GL Bottle 4/Case

#### **HISTO-CLEAR™**

#### HISTO-CLEAR I

- Directly substitute for xylene and yields excellent results in automated tissue processing
- · Superior results with low citrus odor
- · Compatible with most popular mounting media
- Leaves tissue less hard and brittle than xylene for cutting thin sections and prolonging microtome blade life
- Enhances clarity and vibrance of acidophilic stains

#### **HISTO-CLEAR II**

- · Used in dry mounting procedures and aqueous mounting
- Greatly reduced citrus odor compared to other citrus-based histological clearing agents
- · Completely biodegradable to reduce disposal costs

NAT HS200	Histo-Clear™ I, 1GL Bottle	Each
NAT HS202	Histo-Clear™ II, 1GL Bottle	Each

#### **SLIDEBRITE™**

- · Xylene substitute for fast slide drying time
- · Gentle, yet very effective deparaffinization
- · Non-flammable, non-toxic, biodegradable, odorless
- · No vapor masking or disguising additives
- · Enhances nuclear detail
- · No tissue shrinkage or morphology change

ABS SB04	SlideBrite™, 1GL Bottle	4/Case
ABS SMM04	S-Mount®, 4oz Bottle	12/Case









3310 BBC 64001GL

NAT HS200

ABS SB04



**DELTACAL LITE** 

#### **DELTACAL**<sup>™</sup>

- Used for the decalcification of routine, bone marrow core and immunohistochemical specimens
- **BENEFITS:** safe and simple, very effective, reproducible, economical and eco-friendly
- FEATURES: convenient packaging, superior results, better detail in decalcification, biodegradable, non-hazardous, safe and drain disposable

#### **DELTACAL LITE**

· UNIVERSAL effective decalcifying agent

DEX 22002	1L Bottle	2/Pack
DEX 22004	4L Bottle	1/Pack
DEX 22006	4L Bottle	6/Pack
DEX 22012	1L Bottle	12/Pack

#### **DELTACAL**

Very RAPID decalcifying agent

DEX 21002	1L Bottle	2/Pack
DEX 21004	4L Bottle	1/Pack
DEX 21006	4L Bottle	6/Pack
DEX 21012	1L Bottle	12/Pack

#### **DELTAFORM**

· SLOW, MILD decalcifying agent

DEX 23002	1L Bottle	2/Pack
DEX 23004	4L Bottle	1/Pack
DEX 23006	4L Bottle	6/Pack
DEX 23012	1L Bottle	12/Pack

#### **RDO DECALCIFIERS**

#### RAPIF

- Gently decalcifies specimens in a fraction of the time
- · Offers crisp nuclear staining and definition
- · Gives distinguishable chromatin
- · Shows great surface decalcification
- Safe enough to be disposed down lab sinks and drains

#### GOLD

- · Use when biopsies contain minimal amounts of calcium
- Less aggressive decalcification on large and small tissue samples so the nucleic acids within cells are not altered
- · Superb, sharp and crisp nuclear hematoxylin staining
- · No filtration required since it's clear and nonprecipatite
- · Stellar morphological representation of bone biopsies

APX RDO	RDO Rapid, 1GL	Each
APX RDOGOLD	RDO Gold, 1GL	Each

#### RAPID·CAL<sup>™</sup>

- · Specially formulated to decalcify bone rapidly
- Gentle enough to be used as a routine decalcifying solution for bone marrow biopsies
- Usual biopsy obtained from an 11-gauge needle requires about 45-60 minutes to decalcify; biopsies from a larger 8-gauge needle require approximately 60-90 minutes
- Compatible with mercury-containing fixatives such as B-5, Zenker's, Helley's and 10% neutral buffered formalin

BBC 6029	1PT Bottle	Each
BBC 6030	1QT Bottle	Each
BBC 6040	1GL Bottle	Each







APX RDOGOLD

APX RDO

BBC 6040



## MERCEDES SCIENTIFIC° SCOTT'S TAP WATER SUBSTITUTE

#### **BLUING REAGENT**

- Buffered, slightly alkaline solution used in histology and cytology after hematoxylin staining
- Ensures thorough nuclei bluing

MER 11160 1GL Bottle Each





## MERCEDES SCIENTIFIC° HEMATOXYLIN STAINS

- Progressive nuclear stain available in <u>THREE</u> different strengths for cytology and histology applications
- **GILLS I:** Ideal for both gynecological and non-gynecological cytology preparations
- GILLS II: Used as a counterstain in IHC procedures, a frozen section nuclear stain, a routine hematoxylin and eosin (H&E) stain, and for more intense cytology staining
- GILLS III: Primarily for tissue sections when a darker nuclear stain is preferred with shorter staining time

#### **GILLS HEMATOXYLIN**

MER 14796500ML	Gills I, 16oz	Each
MER 147961L	Gills I, 1L	Each
MER 147961GL	Gills I, 1GL	Each
MER 24796500ML	Gills II, 16oz	Each
MER 247961L	Gills II, 1L	Each
MER 247961GL	Gills II, 1GL	Each
MER 34796500ML	Gills III, 16oz	Each
MER 347961L	Gills III, 1L	Each
MER 347961GL	Gills III, 1GL	Each

#### **MODIFIED HARRIS HEMATOXYLIN**

- · No filtering required, formulation contains no mercury
- · Enhanced with glacial acetic acid and ethylene glycol

MER 4799500ML	16oz Bottle	Each
MER 47991L	1L Bottle	Each
MER 47991GL	1GL Bottle	Each

#### **HEMATOXYLIN**

- · Specific nuclear staining of high precision and clarity
- · Vivid, high clarity results with minimal blue background

MER MA1010 500mL Bottle 4/Case

#### RICHARD-ALLAN SCIENTIFIC®

#### **GILLS HEMATOXYLIN**

• Progressive nuclear stains with three strengths available

FIS 22050199	Gills 1, 1PT	4/Case
FIS 22050200	Gills 1, 1GL	Each
FIS 22050201	Gills 2, 1PT	4/Case
FIS 22050202	Gills 2, 1GL	Each
FIS 22050203	Gills 3, 1PT	4/Case
FIS 22050204	Gills 3, 1GL	Each

#### MODIFIED HARRIS HEMATOXYLIN

- · Regressive nuclear stain for both histology and cytology
- More stable than conventional Harris Hematoxylin

FIS 22050206 1GL Bottle Each



## EPREDIA™ SIGNATURE SERIES™ CLARIFIER™ 1, 2

- Enhance the nuclear definition and clarity of Signature Series hematoxylins with Epredia™ Clarifier™ 1, 2
- Improves contrast between Nuclei and Cytoplasm
- Removes background staining caused by excessive tissue adhesive
- Removes unwanted hematoxylin staining from adhesive without affecting nuclear staining
- Clarifier 1 is recommended with Hematoxylin 7211 and Hematoxylin 1, Clarifier 2 is recommended with Hematoxylin 7212 and Hematoxylin 2

FIS 22050116	Clarifier 1™, 1GL Bottle	Each
FIS 22050117	Clarifier 2™, 1GL Bottle	Each

#### **MERCEDES SCIENTIFIC®**

#### **EOSIN Y (1%, W/V IN ALCOHOL)**

- Key component of an Eosin Y Working Solution in the routine hematoxylin and eosin stain and can be used in either manual or automated staining platforms
- Distinguishes between the cytoplasm of different types of cells by staining cytoplasmic components differing shades and intensities of pink to red

MER 4131500ML	16oz Bottle	Each
MER 41311L	1L Bottle	Each
MER 41311GL	1GL Bottle	Each

#### EOSIN

- Alcoholic eosin allows for exceptional differentiation and contrast between cytoplasmic components and nuclei
- · Rinse tissue sections in alcohol prior to staining

MER MA1015 500mL Bottle 4/Case

#### **EOSIN Y ALCOHOLIC** (WITH PHLOXINE B)

- Appropriate amount of acid for proper counterstaining
- Excellent differentiation of cell cytoplasms and other non-nuclear tissue components from bright visuals
- Ideal 3-color spectrum is seen erythrocytes (red), smooth muscle (dull pink), connective tissue (bright pink) - with proper dehydration and differentiation

BBC 3630 1GL Bottle Each

#### PROTOCOL® DEFINE® CLARIFIER

- Removes residual stain that may bind to slide adhesives
- · Use after staining with hematoxylin
- · Use only with progressive staining

FIS 310350 1GL Bottle Each

#### RICHARD-ALLAN SCIENTIFIC®


#### SIGNATURE SERIES EOSIN-Y

- · Bright, vivid cytoplasmic staining
- · Offers excellent delineation of cytoplasmic components
- · Unmatched color contrast between nuclei and cytoplasm
- Various shades of red and pink differentiate between muscle, red blood cells and connective tissue


**FIS 22050110** 1PT Bottle 4/Case

























#### **RAPID·FIX™**

- Fast-acting, specially buffered, formaldehyde based fixative for use at the grossing (cutting) station
- Fixes fatty tissues approximately 50% faster than 10% Neutral Buffered Formalin, and is still safe enough for use on all tissue specimens
- Compatible with all other primary fixatives and does not interfere with any special stain or immunohistochemical procedures
- Improves overall fixation, dehydration, clearing and wax infiltration when used routinely

BBC 0172 1GL Bottle Each

#### FINEFIX FORMALIN-FREE FIXATIVE

- Molecular fixative for optimum H&E, HC and IHC staining as well as fully recovery of nucleic acids
- · Removes formalin from laboratories
- · Storage and transport of specimens in cold climate

FIS 22310293 1GL Bottle Each

#### PARA/GARD® PARAFFIN REPELLENT

- · Xylene-free and non-toxic paraffin repellent
- Quickly and easily removes paraffin from adhering to laboratory surfaces and equipment
- Creates a protective coating to prevent paraffin buildup
- · Safe for use on wood, metal and plastic

TRX PG	4oz Bottle	Each
TRX PG16	16oz Bottle	Each

#### **CLEANEROL™ SPRAY**

- Stain and tissue dye remover, paraffin cleaner and repellent for microtomes, countertops and more
- Biodegradable, earth-friendly, non-toxic, pleasant smell
- Removes wax, ink, blood, grease, mold, adhesives, etc

CAN CL2525 4oz Bottle 4/Case

#### **SOFT BLOCK**

- Softens bone, uterus, nails and other difficultto-cut tissues without the need for acids
- Odorless, technician- and environment-friendly formula won't destroy antigen sites and won't interfere with IHC or In-Situ, special stains, or routine H&E stains
- CAN BE POURED DOWN DRAIN
- WON'T DAMAGE MICROTOME
- Works in 5-15 minutes depending on tissue size and type

MKL ML7052 500mL Bottle Each

#### **NAIL PREP NAIL SOFTENER**

- Exclusive formulation designed for pre-processing nail preparation
- Fast, efficient, and ready to use, with two concentrations
- · Softens brittle nail clippings
- · Makes specimens easier to cut and process
- · Helps to keep tissue from washing off the slide

BBC NP500 500mL Bottle Each

#### A

ITEM #	DESCRIPTION	VOLUME(S)
EKI 1101	Acetic Acid, Glacial	16oz, 2.3L, 1GL
EKI 1109	Acetic Acid Solution, 1%	16oz, 1L, 1GL
EKI 1113	Acetic Acid Solution, 2%	16oz, 1L, 1GL
EKI 1115	Acetic Acid Solution, 3%	4oz, 16oz, 1L, 1GL
EKI 1116	Acetic Acid Solution, 4%	16oz, 1L, 1GL
EKI 1125	Acetic Acid Solution, 5%	16oz, 1L, 1GL
EKI 1127	Acetic Acid Solution, 10%	16oz, 1L, 1GL
EKI 1183	Acetone-Xylenes, 1+1	16oz
EKI 1188	Acid-Alcohol Decolorizer, 1%	16oz, 1GL
EKI 1187	Acid-Alcohol Decolorizer, 3%	16oz, 1L, 1GL
EKI 1184	Acid Fuchsin, Powder	5g, 25g
EKI 1107	Acid Rinse Clarifier/Enhancer	1GL
EKI 4494	Acid Zinc Formalin Fixative	1L, 1GL, 20L
EKI 1193	Alcian Blue 8GX, Powder	5g, 25g
EKI 1198	Alcian Blue Stain, 1%, pH 1.0	8oz, 16oz
EKI 1195	Alcian Blue Stain, 1%, pH 2.5	8oz, 16oz
EKI 1192	Alizarin Red S Powder	25g, 100g
EKI 1191	Alizarin Red S Stain, 1%	16oz, 1L, 1GL
EKI 1222	Aluminum Ammonium Sulfate	500g, 2.5kg
EKI 1260	Aluminum Chloride, Hexahydrate	500g
EKI 1290	Aluminum Chloride, 6 Normal	4oz, 16oz
EKI 1420	Aluminum Potassium Sulfate	500g
EKI 1425	Aminoalkylsilane	8oz, 16oz, 1L
EKI 1556	Ammonium Chloride, Purified	500g
EKI 1625	Ammonium Hydroxide, 28%	16oz, 2.3L
EKI 1956	Aniline Blue, Powder	25g
EKI 1958	Aniline Blue Stain, 2.5%	16oz
EKI 1989	Ascorbic Acid	100g
EKI 1992	Auramine O, Powder	5g, 25g
EKI 1996	Azure A, Powder	25g

#### В

ITEM #	DESCRIPTION	VOLUME(S)
EKI 4493	B-5 Fixative Kit	1L, 1GL
EKI 4495	B-5 Zinc Chloride Fixative	1L, 1GL
EKI 2198	Basic Fuchsin, Powder	10g, 25g, 100g
EKI 2196	Basic Fuchsin Stain, 0.25%	8oz, 16oz
EKI 2197	Basic Fuchsin Stain, 0.5%	8oz
EKI 2199	Benedict Reagent	16oz, 1L, 1GL, 20L
EKI 2248	Biebrich Scarlet, Powder	10g, 25g
EKI 2250	Biebrich Scarlet-Acid Fuchsin Stain	16oz, 1L
EKI 2335	Borax Solution, 5% in Ethanol	16oz, 1L
EKI 2268	Boric Acid, Purified	500g, 2.5kg
EKI 2277	Bouin Fixative	1L, 1GL, 20L

#### В

ITEM #	DESCRIPTION	VOLUME(S)
EKI 2286	Brilliant Cresyl Blue, Powder	10g, 25g
EKI 2288	Brilliant Cresyl Blue, 1% Aqueous	4oz, 16oz
EKI 2289	Brilliant Cresyl Blue, 1% in Saline	4oz, 16oz
EKI 2344	Brilliant Green, Powder	10g, 25g
EKI 2495	Buffer Solution, pH 2.00	16oz
EKI 2500	Buffer Solution, pH 4.01	16oz, 1GL
EKI 2508	Buffer Solution, pH 5.00	16oz
EKI 2512	Buffer Solution, pH 6.00	1L
EKI 2520	Buffer Solution, pH 7.00	16oz, 1GL
EKI 2528	Buffer Solution, pH 8.00	16oz
EKI 2540	Buffer Solution, pH 10.00	16oz, 1GL

#### C

ITEM #	DESCRIPTION	VOLUME(S)
EKI 2009	Canada Balsam, Natural	4oz, 16oz
EKI 2850	Carbol Fuchsin Stain, Kinyoun	16oz, 1L
EKI 2851	Carbol Fuchsin Stain, Ziehl-Neelsen	8oz, 16oz, 1L
EKI 1270	Carmine Alum Lake, Powder	25g
EKI 1272	Carminic Acid, Powder	5g, 25g
EKI 2870	Carnoy Fixative	1L, 1GL
EKI 3040	Cedarwood Oil, Purified	4oz, 16oz
EKI 3056	Charcoal, Activated, Granular	500g
EKI 3100	Chloroform	16oz, 4L
EKI 3267	Chromic Acid Solution, 5%	8oz, 16oz
EKI 3266	Chromic Acid Solution, 10%	16oz
EKI 3260	Chromium Trioxide	125g, 500g
EKI 3270	Chromotrope 2R, Powder	10g, 25g
EKI 3296	Colloidal Iron Stock, Mowry	16oz, 1L
EKI 3314	Congo Red, Powder	5g, 25g
EKI 3450	Congo Red Stain Kit, Bennhold	16oz
EKI 3350	Copper Sulfate, S.G. 1053	16oz, 1L, 1GL, 5GL
EKI 3489	Crystal Violet, Powder	10g, 25g, 100g

#### D

ITEM #	DESCRIPTION	VOLUME(S)
EKI 3600	Davidson Fixative	1GL, 10L, 20L
EKI 3620	Decalcification Endpoint Kit	8oz
EKI 3621	Decalcifier, EDTA-Sucrose	1L, 1GL
EKI 3623	Decalcifier, Formic Acid, 5%	1L, 1GL
EKI 3622	Decalcifier, Formic Acid-Citrate	1L, 1GL
EKI 3626	Decalcifier, Hydrochloric Acid-EDTA	1L, 1GL

Ε

ITEM #	DESCRIPTION	VOLUME(S)
EKI 3995	Ehrlich Aldehyde Reagent	4oz, 16oz
EKI 4132	Eosin-Phloxine Stain Kit	16oz, 1L
EKI 4020	Eosin Y, Powder	25g, 100g
EKI 4133	Eosin Y Stain, 1% Aqueous	4oz, 16oz, 1L, 1GL
EKI 4158	Eriochrome Black T, Powder	5g, 25g
EKI 4044	Ethyl Acetate	16oz, 4L
EKI 4115	Ethylene Glycol Reagent	16oz, 4L

F

ITEM #	DESCRIPTION	VOLUME(S)
EKI 4175	Fast Green FCF, Powder	10g, 25g
EKI 4185	Ferric Chloride, Hexahydrate	500g
EKI 4193	Ferric Chloride Solution, 1%	8oz, 16oz
EKI 4191	Ferric Chloride Solution, 10%	16oz
EKI 4194	Ferric Chloride Solution, 29%	16oz, 1L
EKI 4422	Formaldehyde, 37%	16oz, 1GL, 20L
EKI 4465	Formic Acid, 88%	16oz
EKI 4464	Formic Acid, 96%	16oz
EKI 4467	Fouchet Reagent	8oz, 16oz

G

ITEM #	DESCRIPTION	VOLUME(S)
EKI 1098	Gallego Solution	1L
EKI 4488	Gelatin, 300 Bloom	500g, 2.5kg
EKI 4634	Gentian Violet, 1% Aqueous	4oz, 16oz
EKI 4635	Gentian Violet, 2% Aqueous	4oz, 16oz
EKI 4636	Gentian Violet, 2% Methanol	4oz, 16oz
EKI 4490	Giemsa Stain, Powder	10g, 25g
EKI 8377	Giemsa Stain Phosphate Buffer	16oz, 1L
EKI 4639	Giemsa Stock Stain Solution	16oz, 1L
EKI 4475	Glutaraldehyde, 25%	16oz, 4L
EKI 4476	Glutaraldehyde, 50%	16oz, 4L
EKI 4500	Glycerin	16oz, 4L
EKI 4560	Gold Chloride, Trihydrate	1g
EKI 4561	Gold Chloride Solution, 1%	100mL
EKI 4563	Gold Chloride Solution, 0.2%	100mL, 8oz
EKI 6593	Gomori Trichrome Decolorizer	8oz, 16oz
EKI 4647	Gomori Trichrome Stain	4oz, 16oz
EKI 1185	Gram Decolorizer	16oz, 1L, 1GL
EKI 4633	Gram Gentian Violet Stain	16oz, 1L, 1GL
EKI 6050	Gram Iodine Stain	4oz, 16oz, 1L, 1GL
EKI 11120	Gram Safranine Stain	16oz, 1L, 1GL
EKI 6053	Gram Stain Kit	8oz

Н

ITEM #	DESCRIPTION	VOLUME(S)
EKI 4877	Haupt Gelatin Adhesive	16oz
EKI 4711	Hayem Diluent	4oz, 16oz, 1L
EKI 4794	Hematoxylin, Mayer Modified	16oz, 1L, 1GL
EKI 4880	Hematoxylin, Powder	10g, 25g, 100g
EKI 4791	Hematoxylin, Weigert Stain	16oz
EKI 4642	Hollande Fixative	1L, 1GL
EKI 4653	Hyaluronidase Buffer	1L
EKI 4660	Hydrochloric Acid, 36.5-38%	16oz, 2.3L
EKI 5101	Hydrochloric Acid, 2%	16oz
EKI 5102	Hydrochloric Acid, 10%	1L, 1GL

ı

ITEM #	DESCRIPTION	VOLUME(S)
EKI 4815	Immersion Oil, Microscopy, Type A (Low Viscosity)	1oz, 4oz, 16oz
EKI 4816	Immersion Oil, Microscopy, Type A-B (Medium Viscosity)	1oz, 4oz, 16oz
EKI 4817	Immersion Oil, Microscopy, Type B (High Viscosity)	1oz, 4oz, 16oz
EKI 6055	Iodine, D'Antoni	4oz, 16oz
EKI 6057	Iodine, Masson Trichrome	8oz
EKI 6054	lodine, Lugol Dilute	16oz, 1L, 4L

J

ITEM #	DESCRIPTION	VOLUME(S)
EKI 6404	Jenner Stain, Powder	10g, 25g
EKI 6405	Jenner Blood Stain	16oz, 1L

Κ

ITEM #	DESCRIPTION	VOLUME(S)
EKI 4922	Karnovsky Fixative	1L
EKI 4915	Kovac Aldehyde Reagent	4oz, 16oz

L

ITEM #	DESCRIPTION	VOLUME(S)
EKI 4953	Lactophenol-Cotton Blue Stain	4oz, 16oz
EKI 5114	Leishman Stain, Powder	25g, 100g
EKI 5115	Leishman Stock Blood Stain	16oz, 1L
EKI 6720	Light Green SF Yellowish Powder	10g, 25g
EKI 6721	Light Green SF Yellowish Stain	16oz
EKI 5137	Lithium Carbonate, Saturated	1L, 1GL

#### M

ITEM #	DESCRIPTION	VOLUME(S)
EKI 5345	Magnesium Carbonate	500g, 2.5kg
EKI 5559	Malachite Green Oxalate, Powder	25g
EKI 5560	Malachite Green Stain, 1%	16oz
EKI 6955	Mayer Egg Albumin Adhesive	4oz, 8oz
EKI 5830	Metanil Yellow, Powder	25g
EKI 5831	Metanil Yellow Stain, 0.25%	16oz
EKI 5875	Methylene Blue, Powder	10g, 25g, 100g
EKI 5876	Methylene Blue Stock Solution	8oz, 16oz, 1L
EKI 5868	Methyl Green Stain, 2%	8oz, 16oz
EKI 7310	Methyl Violet 2B, Powder	25g, 100g
EKI 5870	Michel Transport Medium	16oz, 1L, 1GL
EKI 5871	Michel Wash Solution	16oz, 1L, 1GL
EKI 5877	Millonig Buffer Solution	16oz, 1L
EKI 7391	Mucicarmine Stock, Southgate	4oz, 16oz, 1L

#### Ν

ITEM #	DESCRIPTION	VOLUME(S)
EKI 7436	Neutral Red, Powder	25g
EKI 7439	Neutral Red Indicator, 0.1%	16oz
EKI 7437	Neutral Red Stain, 1%	16oz
EKI 7441	New Methylene Blue N, Powder	5g, 10g
EKI 7443	New Methylene Blue N Stain	4oz, 8oz, 16oz, 1L
EKI 6205	Nitric Acid, 70%	16oz, 2.3L
EKI 7709	Nuclear Fast Red, Powder	5g
EKI 7710	Nuclear Fast Red Stain	8oz, 1L

#### C

ITEM #	DESCRIPTION	VOLUME(S)
EKI 6218	Oil Red O, Powder	25g
EKI 6219	Oil Red O Stain in Isopropanol	16oz, 1L
EKI 6209	Oil Red O Stain in Propylene Glycol	16oz
EKI 6220	Orange G, Powder	10g, 25g, 100g
EKI 7695	Orange G Stain, 1%	16oz, 1L

# DON'T SEE THE CHEMICALS YOU'RE LOOKING FOR?

Call a Sales Representative today! at 800.331.2716!

#### Р

ITEM #	DESCRIPTION	VOLUME(S)
EKI 8031	Papanicolaou Stain, EA-36	16oz, 1GL, 10L
EKI 8033	Papanicolaou Stain, EA-50	16oz, 1GL, 10L
EKI 8035	Papanicolaou Stain, EA-65	16oz, 1GL, 10L
EKI 8037	Papanicolaou Stain, OG-6	16oz, 1GL, 10L
EKI 6329	Pararosaniline Hydrochloride, Powder	10g, 25g
EKI 6330	Periodic Acid	25g, 100g
EKI 6315	Periodic Acid Solution, 0.5%	8oz, 16oz
EKI 6316	Periodic Acid Solution, 1%	4oz, 16oz
EKI 6466	Phenol, 88-90%, Liquified, Purified	16oz
EKI 6464	Phenol, 88-90%, Liquified, Reagent	16oz
EKI 6499	Phloxine B, Powder	5g, 25g
EKI 6501	Phosphomolybdic Acid	25g, 100g, 500g
EKI 6503	Phosphomolybdic Acid, 10%	8oz, 16oz
EKI 6510	Phosphomolybdic Phosphotungstic Acid Solution	16oz
EKI 6592	Phosphotungstic Acid, Reagent	25g, 100g, 500g
EKI 8544	Picric Acid, 0.05% in Acetone	16oz
EKI 8540	Picric Acid, 0.1% in Acetone	16oz
EKI 8541	Picric Acid, 1% in Acetone	16oz
EKI 8539	Picric Acid, Saturated, Aqueous	16oz, 1L, 1GL
EKI 6612	Poly L-Lysine Slide Adhesive	8oz, 16oz
EKI 6828	Potassium Bromide	500g
EKI 6845	Potassium Bromide, Purified	500g
EKI 7143	Potassium Ferricyanide, 1%	8oz, 16oz
EKI 7144	Potassium Ferricyanide, 10%	16oz
EKI 7181	Potassium Ferrocyanide	500g, 2.5kg
EKI 7182	Potassium Ferrocyanide, 10%	16oz, 1GL
EKI 7225	Potassium Hydroxide, Pellets	500g
EKI 10080	Potassium Hydroxide, 10%	4oz, 16oz, 1L, 1GL
EKI 10090	Potassium Hydroxide, 20%	4oz, 16oz, 1GL
EKI 7305	Potassium Iodide	125g, 500g
EKI 7422	Potassium Permanganate	125g, 500g
EKI 10745	Potassium Permanganate, 1%	16oz
EKI 7729	Propylene Glycol	16oz, 4L

#### R

ITEM #	DESCRIPTION	VOLUME(S)
EKI 7996	Rhodanine, Powder	5g
EKI 7997	Rhodanine Stock, 0.2% in Alcohol	8oz, 16oz

#### S

ITEM #	DESCRIPTION	VOLUME(S)
EKI 3299	Saccomanno Fixative	1GL, 5GL
EKI 7755	Safranine O, Powder	10g, 25g, 100g
EKI 11121	Safranine O Stain, 0.25%	16oz
EKI 11115	Safranine O Stain, 1%	4oz, 16oz
EKI 11157	Schiff Stain Solution, McManus	4oz, 16oz, 1L
EKI 11160	Scott Tap Water Substitute	16oz, 1L, 1GL
EKI 7844	Silver Nitrate	25g, 100g, 500g
EKI 11450	Silver Nitrate Solution, 5%	16oz
EKI 11452	Silver Nitrate Solution, 10%	16oz
EKI 7940	Sodium Bicarbonate	500g, 2.5kg
EKI 7943	Sodium Bicarbonate Solution, 5%	16oz
EKI 7974	Sodium Borate Solution, 0.5%	16oz
EKI 8087	Sodium Carbonate Solution, 0.1%	16oz, 1L, 1GL
EKI 8165	Sodium Chloride	500g, 2.5kg
EKI 12440	Bloodbank Saline, 0.85%, with Preservative	10L, 20L
EKI 8242	Sodium Citrate Solution, 3%	4oz, 16oz
EKI 8404	Sodium Hydroxide, Pellets	500g, 2.5kg
EKI 12940	Sodium Hydroxide, 10%	4oz, 16oz
EKI 12980	Sodium Hydroxide, 20%	16oz, 1L
EKI 13080	Sodium Hydroxide, 0.100 Normal	16oz, 1L, 1GL, 20L
EKI 13200	Sodium Hydroxide, 1.00 Normal	16oz, 1L, 1GL
EKI 13210	Sodium Hydroxide, 10.00 Normal	1L, 1GL, 10L
EKI 8765	Sodium Phosphate, Dibasic, Heptahydrate	500g, 2.5kg, 12kg
EKI 8680	Sodium Phosphate, Monobasic, Monohydrate	125g, 500g, 2.5kg, 12kg
EKI 8848	Sodium Sulfate Solution, 1%	16oz
EKI 9040	Sodium Thiosulfate, Pentahydrate	500g, 2.5kg
EKI 14007	Sodium Thiosulfate Solution, 2.5%	16oz
EKI 14006	Sodium Thiosulfate Solution, 5%	16oz, 1L
EKI 9330	Sudan Black B, Powder	25g
EKI 9331	Sudan III, Powder	25g, 100g
EKI 14210	Sudan III, Saturated in Alcohol	16oz
EKI 9332	Sudan IV, Powder	25g, 100g
EKI 14286	Sudan IV Stain, Herxheimer	4oz, 16oz
EKI 14300	Sulfosalicylic Acid Solution, 3%	4oz, 16oz, 1L, 1GL
EKI 14302	Sulfosalicylic Acid Solution, 5%	4oz, 16oz, 1L
EKI 14320	Sulfosalicylic Acid Solution, 20%	4oz, 16oz, 1L, 1GL
EKI 9401	Sulfuric Acid	16oz, 2.3L
EKI 14746	Sulfuric Acid Solution, 1%	16oz

#### Т

ITEM #	DESCRIPTION	VOLUME(S)
EKI 14946	Toluidine Blue O, Powder	25g
EKI 14949	Toluidine Blue Stain, 0.1%	8oz, 16oz
EKI 14950	Trichloroacetic Acid, 3%	16oz, 1L
EKI 14951	Trichloroacetic Acid, 5%	16oz, 1L
EKI 14952	Trichloroacetic Acid, 10%	16oz, 1L
EKI 14953	Trichloroacetic Acid, 20%	4oz, 16oz
EKI 14955	Trichloroacetic Acid, 100%	4oz, 16oz
EKI 9670	Trichrome Fecal Stain	16oz

#### V

ITEM #	DESCRIPTION	VOLUME(S)
EKI 9744	Van Giesen Counterstain	16oz, 1L
EKI 9745	Verhoeff-Van Gieson	16oz
	Elastic Tissue Stain Kit	

#### W

ITEM #	DESCRIPTION	VOLUME(S)
EKI 15335	Water, Deionized, High Purity	1GL, 10L, 20L
EKI 9800	Wright Stain, Powder	10g, 25g, 100g
EKI 15385	Wright Stain Solution	16oz, 1L, 1GL
EKI 15388	Wright Stain, Buffered	4oz, 16oz, 1L, 1GL
EKI 15384	Wright Stain Buffer Solution, pH 6.4, Giordano	16oz, 1L, 1GL
EKI 15386	Wright Stain Buffer, pH 6.8	16oz, 1L, 1GL

### X

ITEM #	DESCRIPTION	VOLUME(S)
EKI 9835	Xylenes-Peanut Oil, 2+1	8oz

#### Z

ITEM #	DESCRIPTION	VOLUME(S)
EKI 15395	Zamboni Fixative	1L
EKI 15404	Zenker Fixative	1L, 1GL
EKI 4496	Zinc Formalin Fixative	1L, 1GL, 10L, 20L


	(4.77)		
ACID FAST BACTERIA			
Kit		ES) OR MER FAB2 (100 SLIDES)	
Component	MER AAL500	Acid Alcohol (0.5%) Solution	500mL
	MER CFZ125	Carbol Fuchsin Solution	125mL
	MER LGA125	Light Green Solution	125mL
ALCIAN BLUE (PH 1.0)	)		
Kit	MER AFT1 (800 SLID	ES) OR MER AFT2 (100 SLIDES)	
Component	MER AAG500	Acetic Acid Solution	500mL
	MER ANA250	Alcian Blue (pH 1.0) Solution	250mL
	MER NSF250	Nuclear Fast Red (Enhanced Stability)	250mL
ALCIAN BLUE (PH 2.5	5)		
Kit	MER AFR1 (800 SLID	ES) OR MER AFR2 (100 SLIDES)	
Component	MER AAG500	Acetic Acid Solution	500mL
	MER ANC250	Alcian Blue (pH 2.5) Solution	250mL
	MER NSF250	Nuclear Fast Red (Enhanced Stability)	250mL
ALCIAN BLUE (PAS)			
Kit	MER APS1 (800 SLID	ES) OR MER APS2 (100 SLIDES)	
Component	MER AAG125	Acetic Acid (3%) Solution	2 x 125mL
	MER ANC250	Alcian Blue (pH 2.5) Solution	250mL
	MER HMM125	Hematoxylin, Mayer's (Lillie's)	2 x 125mL
	MER PAQ250	Periodic Acid (1%) Solution	250mL
	MER SRF250	Schiff's Solution	250mL
AMYLOID (CONGO RE	ED)		
Kit	MER AMY1 (1,600 SL	DES) OR MER AMY2 (100 SLIDES)	
Component	MER BRT500	Bluing Reagent	500mL
	MER CRA500	Congo Red Solution	500mL
	MER HMM500	Hematoxylin, Mayer's (Lillie's Modification)	500mL
BIELSCHOWSKY'S ST	AIN KIT (MODIFIED)		
Kit	MER BSK1 (400 SLID	ES)	
Component	MER CASO08	Citric Acid Solution (Bielschowsky's)	8mL Dropper
	MER FRL008	Formalin (20%) Solution	8mL Dropper
	MER NASO08	Nitric Acid Solution (Bielschowsky's)	8mL Dropper
	MER SNZ500	Silver Nitrate (20%) Solution	500mL
	MER STB125	Sodium Thiosulfate (5%) Solution	125mL
CALCIUM (MODIFIED	VON KOSSA)		
Kit	MER CVK1 (400 SLID	ES) OR MER CVK2 (100 SLIDES)	
Component	MER NFS125	Nuclear Fast Red Solution	125mL
	MER SNV125	Silver Nitrate (5%) Solution	125mL
	MER STB125	Sodium Thiosulfate (5%) Solution	125mL
COLLOIDAL IRON			
Kit	MER CIK1 (400 SLIDE	S) OR MER CIK2 (100 SLIDES)	
	MER AAK500	Acetic Acid (12%) Solution	500mL
Component			
Component	MER CIS125	Colloidal Iron Stock Solution	125mL
Component	·	Colloidal Iron Stock Solution  Hydrochloric Acid (1N) Solution	125mL 125mL
Component	MER CIS125		






Cit	MER CSK1 (400 SLIDI	ES) OR MER CSK2 (100 SLIDES)	
Component	MER SAB500	Acetate Buffer (pH 8.0) Solution	2 x 500mL
	MER HMM125	Hematoxylin, Mayer's (Lillie's)	125mL
	MER RSS030	Rhodanine (Stock) Solution	30mL
LASTIC (MODIFIED	VERHOFF'S)		
Cit	MER ETS1 (400 SLIDE	ES) OR MER ETS2 (100 SLIDES)	
Component	MER FCC125	Ferric Chloride (10%, Aqueous) Solution	125mL
	MER FCB125	Ferric Chloride (2%) Differentiating Solution	125mL
	MER HSV250	Hematoxylin (5%) Solution	250mL
	MER LIS125	Lugol's lodine Solution	125mL
	MER STB125	Sodium Thiosulfate (5%) Solution	125mL
	MER VGS125	Van Gieson's Solution	125mL
OSINOPHIL - MAST	CELL		
lit	MER CEM1 (400 SLID	ES) OR MER CEM2 (100 SLIDES)	
Component	MER AST125	Astra Blue Solution	125mL
	MER HMM125	Hematoxylin, Mayer's (Lillie's)	125mL
	MER VRN125	Vital New Red Solution	125mL
ITE'S (LEPROSY & N	NOCARDIA PROCEDURES)		
(it	MER FLS1 (400 SLIDE	ES) OR MER FLS2 (100 SLIDES)	
Component	MER AAM500	Acid Alcohol (1%) Solution	500mL
	MER CFZ125	Carbol Fuchsin Solution	125mL
	MER MBS125	Methylene Blue Solution	125mL
	MER XPO125	Xylene-Peanut Oil Solution	125mL
ONTANA-MASSON	(FOR ARGENTAFFIN CELLS &	MELANIN)	
(it	MER FMS1 (5 ASSAYS	5/35ML)	
Component	MER GCB125	Gold Chloride (0.2%) Solution	125mL
	MER NFS125	Nuclear Fast Red Solution	125mL
	MER SNX009	Silver Nitrate (10%) Solution	5 x 9mL
	MER STB125	Sodium Thiosulfate (5%) Solution	125mL
GIEMSA (MAY-GRUN	WALD)		
Cit	MER GMG1 (1,600 SL	IDES) OR MER GMG2 (100 SLIDES)	
Component	MER GGS500	Giemsa Stock Solution	500mL
	MER MAY500	May-Grunwald Stock Solution	500mL
	MER PBM500	Phosphate Buffer (pH 6.8) Solution	500mL
GMS			
Kit	MER KAA1 (400 SLID	ES) OR MER KAA2 (100 SLIDES)	
Component	MER BOR015	Borax Solution	15mL
	MER CHR125	Chromic Acid Solution	125mL
	MER GCB125	Gold Chloride (0.2%) Solution	125mL
	MER LGA125	Light Green Solution	125mL
	MER MET125	Methenamine Solution	125mL
	MER SNS125	Silver Nitrate (0.2%) Solution	125mL
	MER SBC125	Sodium Bisulfite Solution	125mL
	MER STB125	Sodium Thiosulfate (5%) Solution	125mL

GRAM			
Kit	MER GSK1 (400 SLID	ES) OR MER GSK2 (100 SLIDES)	
Component	MER CFX125	Carbol Fuchsin Counterstain	125mL
·	MER GVS125	Gentian Violet Solution	125mL
	MER GDS125	Gram's Decolorizer Solution	125mL
	MER LIS125	Lugol's Iodine Solution	125mL
	MER TZQ125	Tartrazine Solution	125mL
GRAM (MODIFIED BR	OWN & BRENN)		
Kit	MER BBS1 (400 SLID	ES) OR MER BBS2 (100 SLIDES)	
Component	MER GVS125	Gentian Violet Solution	125mL
	MER GDS125	Gram's Decolorizer Solution	125mL
	MER LIS125	Lugol's lodine Solution	125mL
	MER PAB125	Picric Acid - Acetone (0.1%) Solution	125mL
	MER SOH125	Safranin O Solution	125mL
H. PYLORI			
Kit	MER AYH1 (400 SLID	ES) OR MER AYH2 (100 SLIDES)	
Component	MER ANY125	Alcian Yellow Solution	125mL
	MER PAQ500	Periodic Acid Solution	500mL
	MER SHC008	Sodium Hydroxide (3%) Solution	8mL
	MER SME125	Sodium Metabisulfite Solution	125mL
	MER TQF006	Toluidine Blue Solution	2 x 6mL
H. PYLORI RAPID (1-S	TEP)		
Kit	MER HPS999 (10,000	SLIDES) OR MER HPS500 (5,000 SLIDES)	
Component	MER HPS999 (1,000N	1L), MER HPS500 (500ML)	
HEMATOXYLIN & EOS	SIN		
Kit	MER HAE1 (1,600 SLI	DES) OR MER HAE2 (100 SLIDES)	
Component	MER BRT500	Bluing Reagent	500mL
	MER EYB500	Eosin Y (Modified Alcohol) Solution	500mL
	MER HMM500	Hematoxylin, Mayer's (Lillie's Modification)	500mL
IRON			
Kit			
IXIL	MER IRN1 (400 SLIDE	S) OR MER IRN2 (100 SLIDES)	
	MER IRN1 (400 SLIDE MER HQB500	S) OR MER IRN2 (100 SLIDES)  Hydrochloric Acid (2%) Solution	500mL
			500mL 125mL
	MER HQB500	Hydrochloric Acid (2%) Solution	
Component	MER HQB500 MER NFS125 MER PFB500	Hydrochloric Acid (2%) Solution Nuclear Fast Red Solution	125mL
Component  JONES (FOR BASEME	MER HQB500 MER NFS125 MER PFB500	Hydrochloric Acid (2%) Solution  Nuclear Fast Red Solution  Potassium Ferrocyanide Solution	125mL
Component JONES (FOR BASEME Kit	MER HQB500  MER NFS125  MER PFB500  NT MEMBRANE)	Hydrochloric Acid (2%) Solution  Nuclear Fast Red Solution  Potassium Ferrocyanide Solution	125mL
Component JONES (FOR BASEME Kit	MER HQB500  MER NFS125  MER PFB500  NT MEMBRANE)  MER JSK1 (11 ASSAYS	Hydrochloric Acid (2%) Solution  Nuclear Fast Red Solution  Potassium Ferrocyanide Solution  6/55ML)	125mL 500mL
Component JONES (FOR BASEME Kit	MER HQB500  MER NFS125  MER PFB500  NT MEMBRANE)  MER JSK1 (11 ASSAYS  MER BOR125	Hydrochloric Acid (2%) Solution  Nuclear Fast Red Solution  Potassium Ferrocyanide Solution  5/55ML)  Borax Solution	125mL 500mL 125mL
Component JONES (FOR BASEME Kit	MER HQB500 MER NFS125 MER PFB500 NT MEMBRANE) MER JSK1 (11 ASSAYS MER BOR125 MER GCB125	Hydrochloric Acid (2%) Solution  Nuclear Fast Red Solution  Potassium Ferrocyanide Solution  6/55ML)  Borax Solution  Gold Chloride (0.2%) Solution	125mL 500mL 125mL 125mL
Component JONES (FOR BASEME Kit	MER HQB500 MER NFS125 MER PFB500 INT MEMBRANE) MER JSK1 (11 ASSAYS MER BOR125 MER GCB125 MER MET250	Hydrochloric Acid (2%) Solution  Nuclear Fast Red Solution  Potassium Ferrocyanide Solution  (5/55ML)  Borax Solution  Gold Chloride (0.2%) Solution  Methenamine Solution (2 Bottles)	125mL 500mL 125mL 125mL 250mL
Component JONES (FOR BASEME Kit	MER HQB500 MER NFS125 MER PFB500 NT MEMBRANE) MER JSK1 (11 ASSAYS MER BOR125 MER GCB125 MER MET250 MER NFS125	Hydrochloric Acid (2%) Solution  Nuclear Fast Red Solution  Potassium Ferrocyanide Solution  6/55ML)  Borax Solution  Gold Chloride (0.2%) Solution  Methenamine Solution (2 Bottles)  Nuclear Fast Red Solution (Enhanced)	125mL 500mL 125mL 125mL 250mL 125mL
Component JONES (FOR BASEME Kit	MER HQB500 MER NFS125 MER PFB500 NT MEMBRANE) MER JSK1 (11 ASSAYS MER BOR125 MER GCB125 MER MET250 MER NFS125 MER PAQ250	Hydrochloric Acid (2%) Solution  Nuclear Fast Red Solution  Potassium Ferrocyanide Solution  6/55ML)  Borax Solution  Gold Chloride (0.2%) Solution  Methenamine Solution (2 Bottles)  Nuclear Fast Red Solution (Enhanced)  Periodic Acid Solution	125mL 500mL 125mL 125mL 250mL 125mL 250mL
Component  JONES (FOR BASEME Kit  Component	MER HQB500 MER NFS125 MER PFB500 NT MEMBRANE) MER JSK1 (11 ASSAYS MER BOR125 MER GCB125 MER MET250 MER NFS125 MER NFS125 MER PAQ250 MER SNV030	Hydrochloric Acid (2%) Solution  Nuclear Fast Red Solution  Potassium Ferrocyanide Solution  6/55ML)  Borax Solution  Gold Chloride (0.2%) Solution  Methenamine Solution (2 Bottles)  Nuclear Fast Red Solution (Enhanced)  Periodic Acid Solution  Silver Nitrate (5%) Solution	125mL 500mL 125mL 125mL 250mL 125mL 250mL 30mL
Component  JONES (FOR BASEME Kit  Component	MER HQB500 MER NFS125 MER PFB500  INT MEMBRANE) MER JSK1 (11 ASSAYS MER BOR125 MER GCB125 MER MET250 MER NFS125 MER PAQ250 MER SNV030 MER STB125	Hydrochloric Acid (2%) Solution  Nuclear Fast Red Solution  Potassium Ferrocyanide Solution  6/55ML)  Borax Solution  Gold Chloride (0.2%) Solution  Methenamine Solution (2 Bottles)  Nuclear Fast Red Solution (Enhanced)  Periodic Acid Solution  Silver Nitrate (5%) Solution	125mL 500mL 125mL 125mL 250mL 125mL 250mL 30mL
Component  JONES (FOR BASEME Kit  Component	MER HQB500 MER NFS125 MER PFB500  INT MEMBRANE) MER JSK1 (11 ASSAYS MER BOR125 MER GCB125 MER MET250 MER NFS125 MER PAQ250 MER SNV030 MER STB125	Hydrochloric Acid (2%) Solution  Nuclear Fast Red Solution  Potassium Ferrocyanide Solution  8/55ML)  Borax Solution  Gold Chloride (0.2%) Solution  Methenamine Solution (2 Bottles)  Nuclear Fast Red Solution (Enhanced)  Periodic Acid Solution  Silver Nitrate (5%) Solution  Sodium Thiosulfate (5%) Solution	125mL 500mL 125mL 125mL 250mL 125mL 250mL 30mL
Component  JONES (FOR BASEME Kit  Component  LUXOL FAST BLUE Kit	MER HQB500 MER NFS125 MER PFB500  INT MEMBRANE) MER JSK1 (11 ASSAYS MER BOR125 MER GCB125 MER MET250 MER NFS125 MER PAQ250 MER SNV030 MER STB125  MER LBC1 (400 SLIDI	Hydrochloric Acid (2%) Solution  Nuclear Fast Red Solution  Potassium Ferrocyanide Solution  6/55ML)  Borax Solution  Gold Chloride (0.2%) Solution  Methenamine Solution (2 Bottles)  Nuclear Fast Red Solution (Enhanced)  Periodic Acid Solution  Silver Nitrate (5%) Solution  Sodium Thiosulfate (5%) Solution  Sodium Thosulfate (5%) Solution	125mL 500mL 125mL 125mL 250mL 125mL 250mL 30mL
Component  JONES (FOR BASEME Kit  Component  LUXOL FAST BLUE Kit	MER HQB500 MER NFS125 MER PFB500  INT MEMBRANE) MER JSK1 (11 ASSAYS MER BOR125 MER GCB125 MER MET250 MER NFS125 MER PAQ250 MER SNV030 MER STB125  MER LBC1 (400 SLIDI MER EAS500	Hydrochloric Acid (2%) Solution  Nuclear Fast Red Solution  Potassium Ferrocyanide Solution  8/55ML)  Borax Solution  Gold Chloride (0.2%) Solution  Methenamine Solution (2 Bottles)  Nuclear Fast Red Solution (Enhanced)  Periodic Acid Solution  Silver Nitrate (5%) Solution  Sodium Thiosulfate (5%) Solution	125mL 500mL 125mL 125mL 250mL 125mL 250mL 30mL 125mL





	OLUTION (LOEFFLER'S)		
Kit	MER MBS125 (125ML	.), MER MBS500 (500ML) OR MER MBS999 (1,000ML)	
MOVAT PENTACHEO	OME (MODIFIED MOVAT)		
Kit		PES) OR MER MPS2 (100 SLIDES)	
Component	MER AAE250	Acetic Acid (1%) Solution	250mL
	MER AAG125	Acetic Acid (3%) Solution	125mL
	MER ANC125	Alcian Blue (pH 2.5) Solution	125mL
	MER BSU125	Biebrich Scarlet - Acid Fuchsin Solution	125mL
	MER FCB125	Ferric Chloride (2%) Differentiating Solution	125mL
	MER FCC125	Ferric Chloride (10%) Aqueous Solution	125mL
	MER HSV250	Hematoxylin (5%) Solution	250mL
	MER LIS125	Lugol's lodine Solution	125mL
	MER MYQ125	Metanil Yellow Solution	125mL
	MER PGC250	Phosphotungstic Acid (5%) Solution	250mL
	MER STB125	Sodium Thiosulfate (5%) Solution	125mL
MUCICARMINE (MO	DIFIED SOUTHGATES'S)	· ,	
(it		ES) OR MER SMS2 (100 SLIDES)	
Component	MER BRT125	Bluing Reagent	125mL
·	MER HMM125	Hematoxylin, Mayer's (Lillie's)	125mL
	MER SMG125	Mucicarmine Solution	125mL
	MER TZQ125	Tartrazine Solution	125mL
NUCLEAR FAST RED	SOLUTION (ENHANCED SOLU	UTION)	
Kit		), MER NFS500 (500ML) OR MER NFS999 (1,000ML)	
DIL RED O (FOR FAT	·)		
		ES) OR MER ORK2 (100 SLIDES)	
(it		ES) OR MER ORK2 (100 SLIDES)  Hematoxylin, Mayer's (Lillie's)	125mL
(it	MER ORK1 (400 SLID		125mL 125mL
(it	MER ORK1 (400 SLID MER HMM125	Hematoxylin, Mayer's (Lillie's)	
Cit Component	MER ORK1 (400 SLID MER HMM125 MER ORG125	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution	125mL
Component  DRCEIN (FOR HEPA	MER ORK1 (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution	125mL
Cit Component DRCEIN (FOR HEPA Cit	MER ORK1 (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution Propylene Glycol	125mL
Cit Component DRCEIN (FOR HEPA Cit	MER ORK1 (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)  MER HBK1 (400 SLID	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution Propylene Glycol ES) OR MER HBK2 (100 SLIDES)	125mL 500mL
Cit Component DRCEIN (FOR HEPA Cit	MER ORK1 (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)  MER HBK1 (400 SLID  MER DSA125	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution Propylene Glycol  ES) OR MER HBK2 (100 SLIDES) Differentiating Solution	125mL 500mL 125mL
(it Component DRCEIN (FOR HEPA (it	MER ORK1 (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)  MER HBK1 (400 SLID  MER DSA125  MER OAA125	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution Propylene Glycol  ES) OR MER HBK2 (100 SLIDES) Differentiating Solution Orcein Solution	125mL 500mL 125mL 125mL
(it Component DRCEIN (FOR HEPA (it	MER ORK1 (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)  MER HBK1 (400 SLID  MER DSA125  MER OAA125  MER OQB125	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution Propylene Glycol  ES) OR MER HBK2 (100 SLIDES) Differentiating Solution Orcein Solution Oxalic Acid (2%) Solution	125mL 500mL 125mL 125mL 125mL
(it Component DRCEIN (FOR HEPA (it Component	MER ORK1 (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)  MER HBK1 (400 SLID  MER DSA125  MER OAA125  MER OAB125  MER OQB125  MER PPG030  MER SAQ030	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution Propylene Glycol  ES) OR MER HBK2 (100 SLIDES) Differentiating Solution Orcein Solution Oxalic Acid (2%) Solution Potassium Permanganate (5%) Solution	125mL 500mL 125mL 125mL 125mL 30mL
Component  DRCEIN (FOR HEPA  Kit  Component	MER ORKI (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)  MER HBK1 (400 SLID  MER DSA125  MER OAA125  MER OAB125  MER OQB125  MER PPG030  MER SAQ030	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution Propylene Glycol  ES) OR MER HBK2 (100 SLIDES) Differentiating Solution Orcein Solution Oxalic Acid (2%) Solution Potassium Permanganate (5%) Solution	125mL 500mL 125mL 125mL 125mL 30mL
Cit  Component  Component  Cit  Component  CAPANICOLAOU (PA	MER ORKI (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)  MER HBK1 (400 SLID  MER DSA125  MER OAA125  MER OAB125  MER OQB125  MER PPG030  MER SAQ030	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution Propylene Glycol  ES) OR MER HBK2 (100 SLIDES) Differentiating Solution Orcein Solution Oxalic Acid (2%) Solution Potassium Permanganate (5%) Solution Sulfuric Acid (3%) Solution	125mL 500mL 125mL 125mL 125mL 30mL
Cit  Component  Component  Cit  Component  CAPANICOLAOU (PA	MER ORKI (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)  MER HBKI (400 SLID  MER DSA125  MER OAA125  MER OAB125  MER OGB125  MER PPG030  MER SAQ030  AP)  MER PAP1 (1,600 SLID	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution Propylene Glycol  ES) OR MER HBK2 (100 SLIDES) Differentiating Solution Orcein Solution Oxalic Acid (2%) Solution Potassium Permanganate (5%) Solution Sulfuric Acid (3%) Solution  DES) OR MER PAP2 (100 SLIDES)	125mL 500mL 125mL 125mL 125mL 30mL 30mL
Cit  Component  Component  Cit  Component  CAPANICOLAOU (PACITY)	MER ORKI (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)  MER HBKI (400 SLID  MER DSA125  MER OAA125  MER OQB125  MER PPG030  MER SAQ030  AP)  MER PAP1 (1,600 SLII  MER EAC500	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution Propylene Glycol  ES) OR MER HBK2 (100 SLIDES) Differentiating Solution Orcein Solution Oxalic Acid (2%) Solution Potassium Permanganate (5%) Solution Sulfuric Acid (3%) Solution  DES) OR MER PAP2 (100 SLIDES) EA-50 Stain Solution	125mL 500mL 125mL 125mL 125mL 30mL 30mL
Component  Component  Component  Component  Component  Component  Component  Component	MER ORK1 (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)  MER HBK1 (400 SLID  MER OAA125  MER OAA125  MER OQB125  MER PPG030  MER SAQ030  AP)  MER PAP1 (1,600 SLII)  MER EAC500  MER HMM500  MER OGA500	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution Propylene Glycol  ES) OR MER HBK2 (100 SLIDES) Differentiating Solution Orcein Solution Oxalic Acid (2%) Solution Potassium Permanganate (5%) Solution Sulfuric Acid (3%) Solution  DES) OR MER PAP2 (100 SLIDES) EA-50 Stain Solution Hematoxylin, Mayer's (Lillie's Modification)	125mL 500mL 125mL 125mL 125mL 30mL 30mL
Cit  Component  Component  Cit  Component  Cit  Component  Cit  Component  Cit  Component	MER ORKI (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)  MER HBKI (400 SLID  MER OAA125  MER OAA125  MER OQB125  MER PPG030  MER SAQ030  AP)  MER PAP1 (1,600 SLID  MER HMM500  MER HMM500  MER OGA500	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution Propylene Glycol  ES) OR MER HBK2 (100 SLIDES) Differentiating Solution Orcein Solution Oxalic Acid (2%) Solution Potassium Permanganate (5%) Solution Sulfuric Acid (3%) Solution  DES) OR MER PAP2 (100 SLIDES) EA-50 Stain Solution Hematoxylin, Mayer's (Lillie's Modification)	125mL 500mL 125mL 125mL 125mL 30mL 30mL
Cit  Component  CRCEIN (FOR HEPA  Cit  Component  Cit  Component  Cit  Component  Cit  Component  Cit  Component	MER ORKI (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)  MER HBKI (400 SLID  MER DSA125  MER OAA125  MER OQB125  MER PPG030  MER SAQ030  AP)  MER PAP1 (1,600 SLID  MER EAC500  MER HMM500  MER HMM500  MER OGA500  MIFF (PAS)  MER PAS1 (400 SLID	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution Propylene Glycol  ES) OR MER HBK2 (100 SLIDES) Differentiating Solution Orcein Solution Oxalic Acid (2%) Solution Potassium Permanganate (5%) Solution Sulfuric Acid (3%) Solution  DES) OR MER PAP2 (100 SLIDES) EA-50 Stain Solution Hematoxylin, Mayer's (Lillie's Modification) OG-6 Solution  ES) OR MER PAS2 (100 SLIDES)	125mL 500mL 125mL 125mL 125mL 30mL 30mL 500mL 500mL
Cit Component  CRCEIN (FOR HEPA  Cit Component  Cit Component  Cit Component  Cit Component  Cit Component  Cit Component	MER ORKI (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)  MER HBKI (400 SLID  MER DSA125  MER OAA125  MER OQB125  MER PPG030  MER SAQ030  AP)  MER PAP1 (1,600 SLID  MER EAC500  MER HMM500  MER HMM500  MER OGA500  MER PAS1 (400 SLID)  MER PAS1 (400 SLID)  MER BRT125	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution Propylene Glycol  ES) OR MER HBK2 (100 SLIDES) Differentiating Solution Orcein Solution Oxalic Acid (2%) Solution Potassium Permanganate (5%) Solution Sulfuric Acid (3%) Solution  DES) OR MER PAP2 (100 SLIDES) EA-50 Stain Solution Hematoxylin, Mayer's (Lillie's Modification) OG-6 Solution  ES) OR MER PAS2 (100 SLIDES) Bluing Reagent	125mL 500mL 125mL 125mL 30mL 30mL 500mL 500mL
OIL RED O (FOR FAT  Kit  Component  ORCEIN (FOR HEPAT  Kit  Component  PAPANICOLAOU (PAT  Kit  Component  PERIODIC ACID SCH  Kit  Component	MER ORKI (400 SLID  MER HMM125  MER ORG125  MER PRG500  TITIS B & ELASTIC FIBERS)  MER HBKI (400 SLID  MER DSA125  MER OAA125  MER OQB125  MER PPG030  MER SAQ030  AP)  MER PAP1 (1,600 SLID  MER EAC500  MER HMM500  MER HMM500  MER OGA500  MIFF (PAS)  MER PAS1 (400 SLID	Hematoxylin, Mayer's (Lillie's) Oil Red O Solution Propylene Glycol  ES) OR MER HBK2 (100 SLIDES) Differentiating Solution Orcein Solution Oxalic Acid (2%) Solution Potassium Permanganate (5%) Solution Sulfuric Acid (3%) Solution  DES) OR MER PAP2 (100 SLIDES) EA-50 Stain Solution Hematoxylin, Mayer's (Lillie's Modification) OG-6 Solution  ES) OR MER PAS2 (100 SLIDES)	125mL 500mL 125mL 125mL 125mL 30mL 30mL 500mL 500mL

PERIODIC ACID SCH	IFF (PAS FOR FUNGUS)		
Kit		DES) OR MER PASF2 (100 SLIDES)	
Component	MER LGA125	Light Green Solution	125mL
	MER PAQ250	Periodic Acid Solution	250mL
	MER SRF250	Schiff's Solution	250mL
PICRO-SIRIUS RED (F	FOR CARDIAC MUSCLE)		
Kit	MER SRC1 (800 SLID	ES) OR MER SRC2 (100 SLIDES)	
Component	MER AAD250	Acetic Acid (0.5%) Solution	250mL
	MER PQA250	Phosphomolybdic Acid (0.2%) Solution	250mL
	MER SRS250	Picro-Sirius Red Solution	250mL
PICRO-SIRIUS RED (F	FOR COLLAGEN)		
Kit	MER PSR1 (800 SLID	ES) OR MER PSR2 (100 SLIDES)	
Component	MER AAD250	Acetic Acid (0.5%) Solution	2 x 250mL
	MER SRS250	Picro-Sirius Red Solution	250mL
PNEUMOCYSTIS			
Kit	MER PCS1 (400 SLID	ES) OR MER PCS2 (100 SLIDES)	
Component	MER CEA030	Cresyl Echt Violet (0.1%) Solution	125mL
	MER NYS030	Naphthol Yellow S Solution	125mL
		Staining Jar	60mL
PTAH (PHOSPHOTUN	NGSTIC ACID HEMATOXYLIN)		
Kit	MER PTA1 (400 SLID	ES) OR MER PTA2 (100 SLIDES)	
Component	MER FAS125	Ferric Ammonium Sulfate	125mL
	MER HPA125	PTAH Solution	125mL
	MER ZCS500	Zinc Chloride (10%) Solution	500mL
RETICULUM			
Kit	MER GRT1 (5 ASSAYS		
Component	MER FAT125	Ferric Ammonium Sulfate (3%) Solution	125mL
	MER FRL125	Formalin (20%) Solution	125mL
	MER GCS125	Gold Chloride (0.1%) Solution	125mL
	MER LGA125	Light Green Solution	125mL
	MER PMS125	Potassium Metabisulfite (3%) Solution	125mL
	MER PPE125	Potassium Permanganate (1%) Solution	125mL
	MER SNX010	Silver Nitrate (10%) Solution	5 Vials x 10mL
	MER SHC125	Sodium Hydroxide (3%) Solution	125mL
DETICULUM (MODIE	MER STB125	Sodium Thiosulfate (5%) Solution	125mL
RETICULUM (MODIFI		C/OFMI )	
Kit	MER GRS1 (6 ASSAYS	<u> </u>	1251
Component	MER FAS125	Ferric Ammonium Sulfate Solution	125mL
	MER FRL125	Formalin (20%) Solution	125mL
	MER GCB125	Gold Chloride (0.2%) Solution	125mL
	MER NFS125	Nuclear Fast Red Solution	125mL
	MER PHC015	Potassium Pormanganata Solution	15mL
	MER PPD250	Potassium Permanganate Solution	250mL
	MER PMS125	Potassium Metabisulfite (3%) Solution	125mL
	MER SNX065	Silver Nitrate (10%) Solution	65mL
	MER STB125	Sodium Thiosulfate (5%) Solution	125mL
	MER SAQ015	Sulfuric Acid (1N) Solution	15mL





STEINER /EOR SPIRO	CHETES)		
STEINER (FOR SPIRO Kit	MER SSK1 (6 ASSAYS/	/40ML)	
Component	MER GUM125	Gum Mastic Solution	125mL
•	MER HFS001.5	Hydroquinone	1.5gm
	MER SOS125	Oxidizer Solution	125mL
	MER SNS125	Silver Nitrate (0.2%) Solution	125mL
	MER SNY009	Silver Nitrate (1%) Solution	9mL
	MER ZFN125	Zinc Formalin Solution	125mL
TRICHROME (MODIFI	ED MASSON'S)		
Kit	MER TRM1 (400 SLIDE	ES) OR MER TRM2 (100 SLIDES)	
Component	MER AAE125	Acetic Acid (1%) Solution	125mL
	MER ABP125	Aniline Blue Solution	125mL
	MER BSU125	Biebrich Scarlet/Acid Fuchsin Solution	125mL
	MER BNF125	Bouin's Fluid	125mL
	MER PPA125	Phosphomolybdic/Phosphotungstic Acid Solution	125mL
	MER HWIA125	Weigert's Iron, Hematoxylin (A)	125mL
	MER HWIB125	Weigert's Iron, Hematoxylin (B)	125mL
TRICHROME (MODIFI	ED GOMORI'S)		
Kit	MER TRG1 (400 SLIDE	S) OR MER TRG2 (100 SLIDES)	
Component	MER AAD250	Acetic Acid (0.5%) Solution	250mL
	MER AAL250	Acid Alcohol (0.5%) Solution	250mL
	MER BNF125	Bouin's Fluid	125mL
	MER TGB125	Trichrome Stain (Blue) Solution	125mL
	MER HWIA125	Weigert's Iron, Hematoxylin (A)	125mL
	MER HWIB125	Weigert's Iron, Hematoxylin (B)	125mL
TWORT'S COUNTERS	TAIN		
Kit	MER TCK1 (600 SLIDE	S)	
Component	MER AAF250	Acetic Alcohol (2%) Solution	250mL
	MER TCB030	Twort's Fast Green (Part B)	30mL
	MER TCA125	Twort's Neutral Red (Part A)	125mL
WARTHIN-STARRY			
Kit	MER WSS1 (8 ASSAYS	/20ML)	
Component	MER GSA065	Gelatin (4%), Acidulated	2 x 65mL
	MER HSA030	Hydroquinone (0.1%) Solution, Acidulated	2 x 30mL
	MER SNB125	Silver Nitrate (0.5%) Solution, Acidulated	125mL
	MER SNA030	Silver Nitrate (2%) Solution, Acidulated	30mL
	MER SSE125	SpiroPrep	125mL
WRIGHT-GIEMSA			
Kit		DES) OR MER WGK2 (100 SLIDES)	
Component	MER PBM500	Phosphate Buffer (pH 6.8) Solution	2 x 500mL
	MER WGS500	Wright-Giemsa Solution	500mL





#### **CARBOYS**

- Produced from high density polyethylene (HDPE)
- · Standardized industrial packaging design
- Store, mix and dispense bulk quantities of liquids, solvents and cleaning solutions
- · Features comfortable handles and secure seal caps
- Larger sizes (5L, 10L, 20L) conform to UN regulations on transport of dangerous liquids

GLO 601780	2L Bottle	Each
GLO 601781	5L Bottle	Each
GLO 601782	10L Bottle	Each
GLO 601783	20L Bottle	Each

#### WITH SPIGOT

- · Features very thick walls but retain transparency
- · Can be used in sub-zero temperatures
- DESIGNED FOR LIFETIME SERVICE
- Included spigot (3/4") mates with 12mm deep thread in aspirator for total leak-proof performance

GLO 601660	5L Bottle	Each
GLO 601662	10L Bottle	Each
GLO 601664	25L Bottle	Each
GLO 601666	50L Bottle	Each





### PLASTIC GRADUATED CYLINDERS

- · Made of superior quality translucent polypropylene
- Safe and accurate measurement of all types of liquids
- Wide rim and tapered pour spout for filling and pouring
- Chemically resistant to most acids, bases and solvents
- Non-wetting interior surface eliminates concave menisci
- Heavy-duty plastic withstands repeated autoclaving
- Pentagonal shaped base for stability and roll resistance
- · Greatly reduce breakage and the potential for injury

ITEM #	CAPACITY	HEIGHT	QTY
GLO 601075	10mL	140mm	50/Case
GLO 601077	25mL	195mm	30/Case
GLO 601078	50mL	199mm	30/Case
GLO 601079	100mL	249mm	30/Case
GLO 601080	250mL	315mm	12/Case
GLO 601081	500mL	361mm	12/Case
GLO 601082	1000mL	439mm	6/Case
GLO 601094	2000mL	531mm	6/Case

#### **LOW-FORM BEAKERS**

- · Produced from high clarity polypropylene (PP)
- Space saving design with molded graduations
- Stackable for convenient storage and will not stick together when stacked
- Chemically resistant to most acids, bases and solvents
- Tapered sides for easy handling
- · Precise no-drip pour spout

ITEM #	CAPACITY	HEIGHT	QTY
GLO BPM0025P	25mL	51mm	10/Box
GLO BPM0050P	50mL	61mm	10/Box
GLO BPM0100P	100mL	73mm	10/Box
GLO BPM0250P	250mL	96mm	10/Box
GLO BPM0400P	400mL	111mm	10/Box
GLO BPM0500P	500mL	121mm	10/Box
GLO BPM0600P	600mL	126mm	5/Box
GLO BPM1000P	1000mL	148mm	5/Box
GLO BPM2000P	2000mL	188mm	Each





#### **DISPOSABLE FUNNELS**

- Economical, disposable funnels are available in either transparent polystyrene or translucent polypropylene
- · Lightweight with reinforced rims and inside fluting
- Polypropylene Funnels: Tolerate most strong acids, alkalis and organic solvents
- Polystyrene Funnels: Suitable for aqueous solutions, weak acids and alkalis
- Heat-resistant and withstand temperatures up to 121°C
- · Polypropylene funnels are autoclavable
- Tapered at a 60° angle
- PACKAGED 100/BOX (25/SLEEVE, 4 SLEEVES/BOX)





ITEM #	FUNNEL SIZE	MATERIAL	TOP ID	HEIGHT	STEM LENGTH	STEM OPENING	FILTER PAPER SIZE
GLO 8301	55mm	Polystyrene (PS)	55mm	69mm	25mm	4.3mm	11cm
GLO 8302	65mm	Polystyrene (PS)	65mm	78mm	30mm	4.3mm	12.5cm
GLO 8303	55mm	Polypropylene (PP)	55mm	69mm	25mm	4.3mm	11cm
GLO 8304	65mm	Polypropylene (PP)	65mm	78mm	30mm	4.3mm	12.5cm

#### **ANALYTICAL FUNNELS**

- Heavy-duty, translucent polypropylene (PP) funnels are perfect for use in analytical chemistry applications
- Long, narrow stems are ideal for use with liquids
- Internal and external ribbing helps to prevent air lock
- Angled exactly at 60° for rapid filtration
- · Excellent chemical resistance
- Can withstand temperatures up to 121°C
- Autoclavable at 121°C for 20 minutes
- Available in 9 different sizes
- Designed to accommodate standard filter paper sizes
- PACKAGED 5 PER PACK, 2 PER PACK OR BY THE EACH DEPENDING ON WHICH SIZE FUNNEL





ITEM #	FUNNEL SIZE	CAPACITY	TOP DIAMETER	STEM LENGTH	STEM DIAMETER	QTY
GLO 600145	27mm	3.5mL	25mm	37mm	4mm	5/Pack
GLO 600146	37mm	10mL	35mm	37mm	5mm	5/Pack
GLO 600147	46mm	20mL	50mm	43mm	5mm	5/Pack
GLO 600148	66mm	50mL	65mm	62mm	10mm	5/Pack
GLO 600150	81mm	100mL	80mm	70mm	11mm	2/Pack
GLO 600152	100mm	200mL	100mm	82mm	11mm	2/Pack
GLO 600153	120mm	350mL	120mm	86mm	11mm	2/Pack
GLO 600155	150mm	700mL	150mm	115mm	14mm	Each
GLO 600156	183mm	1,250mL	180mm	140mm	14mm	Each

**CONTACT SALES REPRESENTATIVE TO ORDER OR FOR MORE INFORMATION!** 

















#### **PRE-PRINTED WASH BOTTLES**

- Produced for low density polyethylene (LDPE)
- GHS compliant wash bottle to be only printed in 3 languages (English, French and Spanish)
- Solvent labels include the NEW globally harmonized signal words and pictograms, NFPA diamond, risk phrases and safety information
- · Specially designed jet-tip for precise dispensing control
- DripLok® vapor venting valve prevents solvent drips
- ALL BOTTLES COME FULLY ASSEMBLED
- Comply with the latest GHS and OSHA requirements
- · Color coded for quick solvent recognition
- · Wide mouth opening for easy filling

#### SAFETY-VENTED, 500ML

ITEM #	SOLVENT	CAP	QTY
GLO WGW537VTML	Acetone	Red	Each
GLO WGW538VTML	Methanol	Green	Each
GLO WGW539VTML	Isopropanol	Yellow	Each
GLO WGW541VTML	Ethanol	Orange	Each
GLO WGW601PML	Distilled Water	White	Each
GLO WGW500AST	1 of Each Bottle	Assorted	5/Set

NOTE: DISTILLED WATER BOTTLES ARE NOT SAFETY VENTED

#### **PLAIN WASH BOTTLES**

- Produced from low density polyethylene (LDPE)
- Fitted with color-coded one-piece polypropylene (PP) closures that maintain their shape in use
- Square shoulders and wide mouth for spill-free filling
- Closures have precision jet-tip aperture for fine control
- PACKAGED 5 PER BOX

ITEM # (250ML)	ITEM (500ML)	ITEM # (1000ML)	CAP
GLO WGW230P	GLO WGW236P	GLO WGW443P	White
GLO WGW233P	GLO WGW239P	GLO WGW447P	Green
GLO WGW232P	GLO WGW238P	GLO WGW444P	Yellow
GLO WGW231P	GLO WGW237P	GLO WGW445P	Red
GLO WGW234P	GLO WGW235P	GLO WGW446P	Blue

#### **DO-IT-YOURSELF WASH BOTTLES**

- Do-It-Yourself DripLok® safety-vented wash bottles provided with blank write-on panels
- Indelibly printed in English, French and German
- Supplied with adhesive hazard labels and white closures

GLO WGW001VTML	250mL	5/Box
GLO WGW002VTML	500mL	5/Box
GLO WGW020VTML	1000mL	5/Box



#### **COREDISH® BIOPSY CONTAINERS**

- Multi-compartment containers in the shape of a dish and half prefilled with 10% Neutral Buffered Formalin (NBF), for holding and transporting biopsies
- Supplied with a leakproof closure with "o-ring" ensuring total protection of contents
- Measures only 15 x 95mm in diameter
- Each compartment is clearly identified to allow proper placement and visualization of the biopsy being inserted
- Label permits essential information to be written such as patient identification, doctor, date and time

#### COREDISH® COMPONENTS

- 1) Well labeled, easy-to-read identification of biopsy origin
- 2) Each biopsy is isolated from each other by divider
- 3) Molded ridges around lid make it easy to open and close
- 4) Fast thread design allows lid removed with half turn
- 5) "O-Ring" provides secure and positive leakproof seal
- 6) Warning label has space for patient identification

SIM M970 CoreDish™ Biopsy Containers 10/Case



#### **GRADUATED CLASS A CYLINDERS**

- · Manufactured from borosilicate glass
- · Hexagonal base for extra stability
- · Spout for easy, mess free pouring
- Includes protective guard and printed graduations
- Capacity in accordance with ASTM™ E1272 standards

ITEM #	CAPACITY	SUBDIVISION	QTY
EIS CH0341ABL	10mL	0.2mL	Each
EIS CH0341BBL	25mL	0.5mL	Each
EIS CH0341CBL	50mL	1.0mL	Each
EIS CH0341DBL	100mL	1.0mL	Each
EIS CH0341EBL	250mL	2.0mL	Each
EIS CH0341FBL	500mL	5.0mL	Each
EIS CH0341GBL	1000mL	10.0mL	Each
EIS CH0341HBL	2000mL	20.0mL	Each

500ML









- Low form with spout made from borosilicate glass
- Graduated as per DIN 12231 ISO 3819

**SIZES AVAILABLE:** 5mL | 10mL | 25mL | 50mL | 100mL | 150mL | 250mL | 400mL | 500mL | 600mL | 800mL | 1,000mL | 2000mL | 3000mL | 5,000mL | 10,000mL

EIS CH0126

Low Form Glass Beaker

Each



#### **CONICAL ERLENMEYER FLASKS**

· Made of borosilicate glass

#### WITH TEFLON LINER SCREW CAP

SIZES AVAILABLE: 100mL | 250mL | 500mL | 1,000mL

#### NARROW-NECK

**SIZES AVAILABLE:** 25mL | 50mL | 100mL | 150mL 200mL | 250mL | 300mL | 500mL | 1,000mL | 2,000mL 3,000mL | 5,000mL

#### **WIDE-NECK**

**SIZES AVAILABLE:** 25mL | 50mL | 100mL | 250mL 500mL | 1,000mL | 2,000mL

EIS CH0430	Flask w/Teflon Liner Cap	Each
EIS CH0424	Narrow-Neck Flask	Each
EIS CH0426	Wide-Neck Flask	Fach



#### **REAGENT BOTTLES**

- · Made from high-quality borosilicate glass
- Autoclavable screw cap (GL45) and pouring ring for drip free pouring and clean, safe working
- Complies with DIN ISO 4796-1 standards







#### **CLEAR GLASS**

EIS CH0164A	100mL	Each
EIS CH0164B	250mL	Each
EIS CH0164C	500mL	Each
EIS CH0164D	1000mL	Each
EIS CH0164E	2000mL	Each







#### AMBER GLASS

EIS CH0165A	100mL	Each
EIS CH0165B	250mL	Each
EIS CH0165C	500mL	Each
EIS CH0165D	1000mL	Each
EIS CH0165E	2000mL	Each



#### **WASTE CANS**

- · Step-on foot pedal opens lid
- · Ideal for exam rooms, restrooms, and front office areas
- Meets OSHA standard on blood-borne pathogens
- · Hazardous waste label included
- · AVAILABLE IN RED. WHITE, AND BEIGE

#### **RECTANGULAR PLASTIC SERIES**

- · Plastic can and lid
- LARGE: 52QT (13GL) Capacity, 12"D x 15"W x 26"H
- SMALL: 32QT (8GL) Capacity, 12"D x 15"W x 19"H

#### **SQUARE STEEL SERIES:**

- Steel outer can and lid with removable plastic inner liner
- LARGE: 32QT (8GL) Capacity, 111/2"D x 121/4"W x 181/4"H
- **SMALL:** 20QT (5GL) Capacity, 11"D x 11<sup>3</sup>/<sub>4</sub>"W x 17<sup>1</sup>/<sub>4</sub>"H

#### **ROUND STEEL SERIES:**

- Steel outer can and lid with removable plastic inner liner
- STAINLESS STEEL FINISH AVAILABLE
- SIZE: 32QT (8GL) Capacity, 111/2"D x 24"H

BRE 252XX	Rectangular Plastic	Each
BRE 352XX	Square Steel	Each
BRE 452XX	Round Steel	Each

#### **CHAIRS & STOOLS**

- Heavy duty steel, lightweight aluminum and upholstered components provide stability and strength
- Optional high resistance casters for free chair movement when weighted or occupied, but resist unnecessary chair movement when unweighted or exiting the chair
- Each stool or chair has a capacity of 250lbs

#### **AVAILABLE IN 25+ VINYL COLORS**

Examples from Classic Vinyl Colors (Group A):

02: Royal Blue06: Hunter Green14: Burgundy03: Taupe07: Gray16: Mauve04: Moss Green08: Black17: Blueridge



ITEM #	DESCRIPTION	DESCRIPTION	SEAT HEIGHT	WEIGHT
WIN 4300	Gas-Lift Task Stool	5-Star Base, Comfortable Padded Seat Upholstered in High-Grade Vinyl, "One-Touch" Height Adjustment	16½" - 21½"	14lbs
WIN 4400	Deluxe Gas-Lift Task Stool	5-Star Polished Aluminum Base, Comfortable Padded Seat Upholstered in High-Grade Vinyl, "One-Touch" Height Adjustment	17³/4" - 23"	15.5lbs
WIN 4350	Gas-Lift Task Chair	Seat-Back that Raises or Lowers, 5-Star Base, Comfortable Padded Seat Upholstered in High-Grade Vinyl, Gas-Lift Height Adjustment	173/4" - 223/4"	20lbs
WIN 4450	Deluxe Gas-Lift Task Chair	Adjustable Back Height, Polished Aluminum 5-Star Base, Fingertip Height Adjustment	181/4" - 231/4"	21.5lbs
WIN 4380	Tall Gas-Lift Task Chair	Adjustable Back, Includes Foot Ring, 16" Diameter Seat Quickly Raises and Lowers with Fingertip Height Adjustment	23½" - 33¾"	22lbs

#### **FLAMMABLE STORAGE CABINETS**

- Sturdy, all-welded double wall steel construction
- 2" wide louvers support shelves, adjustable on 4" centers
- · Safely contains hazardous chemicals to reduce fires
- Electrical ground on bottom right hand side aids in the elimination of an electrostatic discharge
- Improved security with recessed lockable paddle latch and 3-point locking system
- Bunged air vents with flame arrestors on each side of the cabinet allows for ventilation
- · SAFETY YELLOW POWDER COAT FINISH

DHM 1004M50	<sup>4</sup> GL, Yellow, <sup>1</sup> Shelf, <sup>173</sup> / <sub>8</sub> "W x <sup>181</sup> / <sub>8</sub> "D x <sup>221</sup> / <sub>8</sub> "H	Each
DHM 1012M50	<sup>12</sup> GL, Yellow, <sup>1</sup> Shelf, <sup>23</sup> "W x <sup>18</sup> "D x <sup>35</sup> "H	Each
DHM 1016M50	<sup>16</sup> GL, Yellow, <sup>1</sup> Shelf, <sup>23</sup> "W x <sup>18</sup> "D x <sup>44</sup> "H	Each
DHM 1030M50	<sup>30</sup> GL, Yellow, <sup>1</sup> Shelf, <sup>43</sup> "W x <sup>18</sup> "D x <sup>44</sup> "H	Each
DHM 1045M50	$^{45}$ GL, Yellow, $^2$ Shelves, $^{43}$ "W x $^{18}$ "D x $^{65}$ "H	Each
DHM 1060M50	<sup>60</sup> GL, Yellow, <sup>2</sup> Shelves, <sup>34</sup> "W x <sup>34</sup> "D x <sup>65</sup> "H	Each

## ACID & CORROSIVE STORAGE CABINETS

- Leak tight sump and shelves are equipped with removable corrosive-resistant polyethylene trays to contain accidental spills
- Finished in corrosive-resistant blue urethane paint
- LARGE WARNING LABEL: "CAUTION CORROSIVES"

AAS PE3045	<sup>24</sup> GL, <sup>18</sup> "H x <sup>38</sup> "W x <sup>13</sup> "D,	Each
	High Density Polyethylene	
AAS C330	<sup>30</sup> GL, <sup>46</sup> "H x <sup>43</sup> "W x <sup>18</sup> "D,	Each
	One Shelf, Self-Close,	
	Self-Latch Safe-T-Door	
AAS C345	<sup>45</sup> GL, <sup>67</sup> "H x <sup>43</sup> "W x <sup>18</sup> "D,	Each
	Two Shelves, Self-Close,	
	Self-Latch Safe-T-Door	
AAS C360	<sup>60</sup> GL, <sup>67</sup> "H x <sup>31</sup> "W x <sup>31</sup> "D,	Each
	Two Shelves, Self-Close,	
	Self-Latch Safe-T-Door	





#### **ABSORBANCE PADS**

- · Impermeable to strong solvents (i.e. xylene, alcohols)
- · Clean labware, absorb spills and contain leaks

CUR BH32000	3" x 3", Heavy	2000/Case
CUR BH44200	4" x 4", Heavy	2000/Case
CUR BH81020	8" x 9", Heavy	400/Case
CUR BH10600	16" x 100', Heavy	2 Rolls/Case
CUR BH93016	16" x 100', Super	2 Rolls/Case
CUR BH16100	16" x 100', Standard	2 Rolls/Case

#### ABSORBANCE PAD HOLDERS (4" X 4")

CUR BH10000	White Cardboard	10/Pack
CUR AC9512	Clear Acrylic	Each

#### **WOVEN COTTON GAUZE SPONGES**

- Made from 100% woven cotton
- Ideal for applications including general cleaning, wound dressings, prepping, packing and debriding wounds
- · Non-sterile gauze suits multiple market needs

DUK 2085	2" x 2", 8-Ply	200/Bag
DUK 2128	2" x 2", 12-Ply	200/Bag
DUK 3124	3" x 3", 12-Ply	200/Bag
DUK 4084	4" x 4", 8-Ply	200/Bag
DUK 4122	4" x 4", 12-Ply	200/Bag

#### BENCHKOTE® SURFACE PROTECTOR

- Coated with an impermeable layer of polyethylene to protect from liquid spills, pathogens and bacteria
- · Strong and tear resistant when wet or dry
- Smooth surface may be written on with ink or pencil
- · Suitable for saturation with disinfectant
- Can be safely incinerated after use

**FIS 12007182A** 18" x 164', White Roll





#### **NON-STERILE GAUZE SPONGES**

- Made from 100% premium bleached cotton gauze
- · Bleached for maximum absorbency
- Applicable for general cleaning, wound dressing, prepping, packing and debriding wounds


PRA P157115	2" x 2", 8-Ply	200/Sleeve
PRA P157112	2" x 2", 12-Ply	200/Sleeve
PRA P157116	3" x 3", 8-Ply	200/Sleeve
PRA P157117	4" x 4", 8-Ply	200/Sleeve
PRA P157114	4" x 4", 12-Ply	200/Sleeve
PRA P157118	4" x 4", 16-Ply	200/Sleeve

#### WHITE COUNTERTOP PAPER

- Polybacking provides extra strength, leak-proof protection, softness and reduced clean-up time
- · Perforation adds convenience and neatness
- 99.9% protection against bacteria pass-through
- Covers exam table sizes from pediatrics to bariatrics

TID 980980	20" x 150', Rolls	4/Case
TID 980982	12" x 20", Sheets	600/Case
TID 980984	18" x 20", Sheets	400/Case











#### SUPER SANI-CLOTH® WIPES

- · Ideal for daily use in fast-paced environments
- Test effective against 30 microorganisms
- · Thick, strong, pliable material with high fluid capacity reaches crevices and difficult corners of equipment
- 2 MINUTE CONTACT TIME

PDI Q55172	Large, 6" x 6 <sup>3</sup> / <sub>4</sub> "	160/Can
PDI Q86984	X-Large, 8" x 14"	65/Can
PDI H04082	Large, 5" x 8"	50/Box
PDI U87295	X-Large, 11½" x 11¾"	50/Box

#### SANI-CLOTH® HB WIPES

- · Alcohol free disinfection formulation
- · Test effective against 100 microorganisms
- 10 MINUTE CONTACT TIME

PDI Q08472	Large, 6" x 63/4"	160/Can
PDI Q85484	X-Large, 8" x 14"	65/Can

#### SANI-CLOTH® PLUS WIPES

- Ideal for alternate care settings including physician offices on hard nonporous surfaces and equipment
- · Low alcohol disinfection formulation
- · Test effective against 16 microorganisms
- 3 MINUTE CONTACT TIME

PDI Q89072	Large, 6" x 6¾"	160/Can
PDI Q85084	X-Large, 8" x 14"	65/Can

#### SANI-CLOTH® BLEACH WIPES

- · Ideal for disinfecting high risk areas endemic with Clostridium difficile spores and Norovirus
- · Sodium hypochlorite disinfection formulation
- · Thick, strong wipes with high fluid capacity and pliable to get into crevices and corners of equipment
- · Test effective against 50 microorganisms
- 4 MINUTE CONTACT TIME

PDI P54072 X-Large, 6" x 101/2" 75/Can



#### **KIMWIPES®**

- Wiper easily picks up liquid, dust and tiny particles
- Antistatic dispensing reduces lint and electrostatic discharge for sensitive environments like labs
- Features one-at-a-time dispensing to reduce waste
- Ideal for cleaning laboratory instruments, polishing glassware and other delicate surfaces
- Push-Up Box holds 5" x 5" Kimwipes®

KIM 34155	White, 41/2" x 81/2"	280/Box
AME TWDO0101E	E" v E" Disponsor	Fach



#### KIMWIPES® BOX HOLDERS

- **Dimensions:** 45/8" x 5" x 31/8" (Small)
- May be mounted to wall or lab cabinets for easy access with screws or double sided tape
- One-handed removal of tissues
- Rubber feet prevent holders from sliding on counters

MIT VP170G	Small, Neon Green	Each
MIT VP170R	Small, Neon Red	Each
MIT VP180	2 Sided, Clear	Each



#### PROADVANTAGE® HAND SANITIZERS

- Helps condition, soften and disinfect with 62% alcohol
- Kills over 99% of germs in as little as 15 seconds
- · Rubs quickly for clean, sanitized hands and skin

PRA P779114	4oz Bottle, Flip-Cap	24/Case
PRA P779118	10oz Pump Bottle	12/Case



#### **PURELL® INSTANT HAND SANITIZERS**

- Refreshing gel formulation Leaves hands feeling soft
- Kills more than 99.99% of illness causing germs
- · Clinically-proven to maintain skin health

GOJ 960524	2oz Flip-Cap Bottle	24/Case
GOJ 960624	2oz Pump Bottle	24/Case
GOJ 9651	4oz Pump Bottle	12/Case
GOJ 9652	8oz Pump Bottle	12/Case
GOJ 962504	2L Pump Bottle	4/Case
GOJ 302312	20oz Pump Bottle	12/Case



#### WYPALL® L40 WIPERS

- · Ideal for routine industrial cleaning and maintenance
- · Removes soils and liquids from face and hands
- · Easily absorbs large liquid spills and leaks
- 56 Sheets/Box, 18 Boxes/Case

**KIM 05701** White, 12.5" x 12" 1,008/Case

#### WYPALL® X50 WIPERS

- · Reinforced for strength when wet or dry
- · Smaller sheet size means less waste
- · 26 Sheets/Pack, 32 Packs/Case

**KIM 35025** White, 10" x 12.5" 832/Case





## WYPALL® X60 WIPERS

- · Reinforced for extra scrubbing and cleaning power
- Tough enough for big jobs, soft enough for face or hands
- · No adhesives or binders, leaves no residue

**KIM 34955** White, 12.5" x 13.4" 1,100/Roll

#### WYPALL® X80 SHOP TOWELS

- · Absorbs more water and oil than textile rags
- · Resistant to solvents and other harsh cleansing agents

KIM 41055	Red, 12.5" x 13.4"	475/Roll
KIM 41025	White, 12.5" x 13.4"	475/Roll



#### WYPALL® MICROFIBER CLOTHS

- · Designed with Microban® protection
- Launderable and reusable for up to 300 washings

KIM 83610	Gold, 15 <sup>3</sup> / <sub>4</sub> " x 15 <sup>3</sup> / <sub>4</sub> "	24/Case
KIM 83620	Blue, 15¾" x 15¾"	24/Case
KIM 83630	Green, 15 <sup>3</sup> / <sub>4</sub> " x 15 <sup>3</sup> / <sub>4</sub> "	24/Case

#### WYPALL® DRY-UP™ L40 TOWELS

- Multi-purpose, absorbent bath-sized towels reduce cost and hassles of laundering and pilferage
- Ideal for health clubs, manufacturing and schools
- Packaged in pop-up box

**KIM 05860** White, 19.5" x 42" 200/Box

the



- Excellent supplement to existing safety equipment
- Packaged in an instant access container that can be hung on the wall or stored in any accessible location
- Made of woven fiberglass with a textured surface; all edges are serged with fireproof Kevlar\* thread
- Close weave greatly reduces the amount of atmospheric oxygen available to feed the fire
- May be wrapped around the head and body to protect against flames and embers when fleeing a burning room, or to protect from splashing molten metal

**BEL H248690000** 36" x 40" Each



MDL FA024300001Z

#### **EMERGENCY EYEWASH STATIONS**

- · Quick and easy access to eye wash for emergencies
- · Includes sterile, irrigating solution for flushing the eyes
- Plastic squeeze bottle and nozzle help create a fluid stream to direct solution into the affected areas

MDL FAO24000001Z	Eye Wash Station w/(1) 16oz Solution Bottle	Each
MDL FAO24102001Z	Eye Wash Station w/(2) 16oz Solution Bottles	Each
MDL FAO24300001Z	Eye Wash Station w/(2) 32oz Solution Bottles	Each
MDL FAO24101001	16oz Solution Refill	12/Case
MDL FAO24201001	32oz Solution Refill	12/Case











#### **OPTI-KLENS® EMERGENCY EYEWASH**

- Chrome-plated, solid brass eyewash unit easily installs on most existing sink or basin faucets
- Unlimited water supply unlike portable models
- Frees both hands to hold eyes open to water flow
- Delivers more than twice the required water flow
- · Reduces high water pressure to prevent injury
- · Two of the most common adapters included
- Eliminator valve presets faucet water temperature, activate water flow by use of push through plunger
- · Meets ANSI and OSHA standards

DES 329200	Opti-Klens™ I Eyewash Station	Each
DES 469022	Opti-Klens™ Wall Mounted Face & Eyewash Station	Each
DES 269422	Opti-Klens™ II Eyewash Fountain	Each
DES 269622	Opti-Klens™ II Eyewash Fountain with Eliminator Valve	Each









VHI BH2004

WHI BH2008

WHI BH2005

WHI BH2003 WHI BH2009

#### BIOHAZARD KEEPER™ - FLOOR

- Sturdy corrugated paperboard construction for strength with 2 mil red plastic liner for leak-resistance
- Marked with orange biohazard waste warnings printed in English and Spanish (NOT A SHARPS CONTAINER)
- Unique wide-mouth opening with reclosable, removable cover for disposing larger items and reducing odors
- Separate flap can be folded over the top to permanently close the container for final disposal

**WHI BH2009** 12"W x 12"D x 27"H 6/Case

#### **BROKEN GLASS KEEPER™ - FLOOR**

- Sturdy corrugated container lined with a clear, thick plastic bag and marked with broken glass symbols
- Unique wide-mouth opening with reclosable, removable cover for disposing larger glass or sharp items
- Inner flap seals contents for final disposal when full

**WHI BH2008** 12"W x 12"D x 27"H 6/Case

#### BIOHAZARD BENCHTOP KEEPER™

- Free standing paperboard container stores contaminated waste including petri dishes, test tubes, applicator sticks, pipets and other non-sharp items at your workstation
- · Plastic lining holds moderate amounts of liquid
- · Self-locking tops keep contents safely enclosed
- Autoclave contaminated materials then dispose of according to your local, state or federal regulations

WHI BH2002	4"L x 4"W x 7"H	100/Case	
WHI BH2003	6"L x 6"W x 10"H	50/Case	

#### BROKEN GLASS BENCHTOP KEEPER™

- Sturdy paperboard container has a plastic lining to contain moderate amounts of liquid
- Safely collect non-hazardous broken glass at workstation
- Securely close using the interlocking tabs on the lid and dispose of according to your facility's instructions

**WHI BH2004** 6"L x 6"W x 10"H 50/Case





















**SHARPS® CONTAINERS** 

- Made of autoclavable polypropylene to resist puncture
- Stackable design provides a convenient stable platform during use, economical storage and lower shipping costs
- Built-in phlebotomy needle port for easy disposal
- · Wide horizontal drop lid models accommodate larger syringes and allow full utilization of fill line capacity
- Transparent lids facilitate visual identification of fill level
- Translucency facilitate easy side viewing of content level
- Compact containers with locking lids for easy transport

#### STACKABLE

Biohazard Labeling, Biohazard Symbol

ITEM #	CAPACITY	DIMENSIONS	COLOR
MEX 8702	1QT	3½ x 7 x 3½"	Red
MEX 8702T	1QT	3½ x 7 x 3½"	Translucent
MEX 8703	1GL	10 x 5 x 7"	Red
MEX 8703T	1GL	10 x 5 x 7"	Translucent
MEX 8704	2GL	10 x 9½ x 7"	Red
MEX 8704T	2GL	10 x 9½ x 7"	Translucent
MEX 8705	8GL	13 <sup>3</sup> / <sub>4</sub> x 14 x 13 <sup>3</sup> / <sub>4</sub> "	Red
MEX 8705T	8GL	13 <sup>3</sup> / <sub>4</sub> x 14 x 13 <sup>3</sup> / <sub>4</sub> "	Translucent
MEX 8716	16GL	13½ x 25½ x 13¾"	Red

#### STACKABLE

Tortuous Path Lid, Biohazard Labeling, Biohazard Symbol

ITEM #	CAPACITY	DIMENSIONS	COLOR
MEX 8706	1GL	10 x 6 <sup>3</sup> / <sub>4</sub> x 7"	Red
MEX 8706T	1GL	10 x 6 <sup>3</sup> / <sub>4</sub> x 7"	Translucent
MEX 8707	2GL	10 x 11½ x 5"	Red
MEX 8707T	2GL	10 x 11½ x 5"	Translucent
MEX 8710	3GL	12½ x 13½ x 6"	Red

#### **CHEMO WASTE CONTAINERS**

CYTA Symbol, Chemotherapy Waste

ITEM #	CAPACITY	DIMENSIONS	COLOR	
MEX 9751	8GL	10 x 6 <sup>3</sup> / <sub>4</sub> x 7"	White	
MEX 932	23QT	12 x 9 x 18"	White	

#### NONSTACKABLE

Biohazard Labeling, Biohazard Symbol

ITEM #	CAPACITY	DIMENSIONS	COLOR
MEX 182	6.2QT	10 x 6 x 11"	Red
MEX 184	4.8QT	10 x 6 x 9"	Red
MEX 185	1.7QT	3½ x 3½ x 10"	Red
MEX 186S	12QT	12 x 9 x 11"	Translucent
MEX 186W	12QT	12 x 9 x 11"	Red
MEX 187	23QT	12 x 9 x 18"	Red
MEX 187S	23QT	12 x 9 x 18"	Translucent
MEX 195	2QT	5 x 4 x 10"	Translucent
MEX 199	1.3QT	5 x 4 x 7½"	Translucent

141













#### Item # ALC 2001

#### **CITRAJET**®

- · Low foaming, phosphate-free, concentrated, liquid cleaner and metal brightener for high performance removal of metal oxides, scales, salts and inorganic residues in power-spray systems, washers, and manual and ultrasonic cleaning
- · Corrosion inhibited and free rinsing provides reliable results and eliminates interfering residues
- · Ideal for acid washing and rinsing

#### Item # ALC 1715

#### DETERGENT 8°

- · Nonionic, low-foaming, phosphate-free detergent for labware washers, PC board washers, power-spray systems, ultrasonic and
- · Removes oils, resins rosins, and fluxes from PC circuit boards, electronic parts, phosphate sensitive labware and nuclear contaminated equipment
- · Excellent replacement for solvent cleaning

#### Item # ALC 1801

#### **CITRANOX**®

- Phosphate free, concentrated cleaner and metal brightener for high performance removal of metal oxides, scale, salts and inorganic residues in manual and ultrasonic cleaning
- · Corrosion inhibited and free rinsing provides reliable results and eliminates interfering residues
- Ideal for trace metal analysis

#### Item # ALC 1201

#### LIQUI-NOX<sup>®</sup>

- · Concentrated, anionic phosphate free, detergent for manual and ultrasonic cleaning
- · Free rinsing provides reliable results and eliminates interfering residues
- · Extremely mild and completely soluble in hard and soft water

#### Item # ALC 1104

#### **ALCONOX**<sup>®</sup>

- · Concentrated, anionic detergent for manual and ultrasonic cleaning
- Free rinsing provides reliable results and eliminates interfering residues
- Ideal for cleaning contaminants from glassware, metals, plastics, ceramic, porcelain, rubber and fiberglass
- Excellent for corrosive acids and hazardous solvents

#### Item # FIS 16000141

#### **LUMINOX**<sup>®</sup>

- · Low foaming, neutral pH detergent for use in manual, ultrasonic, and machine cleaning
- Safe-handling detergent used for non-corrosive cleaning without chelation or alkalinity waste treatment problems
- Recommended for aluminum and soft metal, medical devices, pharmaceutical equipment, alkaline sensitive labware, and medical instruments in aluminum trays

#### Item # ALC 21011

#### **SOLUJET®**

- Nonionic, low-foaming, phosphate-free liquid cleaner for use in labware washers, manual and ultrasonic cleaning
- · Corrosion inhibited formula recommended for glassware, metals, plastic, ceramic, porcelain, rubber and fiberglass
- Free rinsing provides reliable results and eliminates interfering residues

#### Item # ALC 1500

#### **ALCOTABS®**

- · Effervescent detergent tablet for maximum cleaning in syphon pipette washers to clean pipettes and tubes
- · Slow release bubble action
- Free rinsing provides reliable results and eliminates interfering residues

#### CONTACT SALES REPRESENTATIVE FOR SIZES, QUANTITIES OR MORE INFO!













Item # ALC 1325

#### TERG-A-ZYME®

- Concentrated, anionic detergent with protease enzyme for manual and ultrasonic cleaning
- Excellent for removal of proteinaceous soils, tissue, blood and body fluids from glassware, metals, plastic, ceramic, porcelain, rubber and fiberglass
- Ideal as a cleaning agent in Reverse Osmosis and Ultra-Filtration Systems

Item # ALC 2001

#### **TERGAJET®**

- Low foaming, phosphate-free powdered cleaner for use in labware washers, manual and ultrasonic cleaning
- Corrosion inhibited formula recommended for glassware, metals, plastic, ceramic, porcelain, rubber and fiberglass
- Free rinsing provides reliable results and eliminates interfering residues

Item # ALC 1425

#### **ALCOJET®**

- Nonionic, low foaming powder cleaner for use in labware washers, parts washers, power-spray systems, manual and ultrasonic cleaning
- Corrosion inhibited formula recommended for glassware, metals, plastic, ceramic, porcelain, rubber and fiberglass
- Free rinsing provides reliable results and eliminates interfering residues

#### **SPARKLEEN® 1 DETERGENT**

- Universal detergent for manual or hand washing laboratory glassware and equipment
- Removes proteins, fats, oils, media, wax, dried blood and other organic-inorganic deposits from glassware

**FIS 043204** 3.25lb Carton Each

#### SPARKLEEN® 2 DETERGENT

- Effectively cleans the most tenacious soils such as grease, petrolatum, dried blood and organic deposits
- Used for mechanical spray washing of glassware, stainless steel and most plastics

**FIS 043206** 3.25lb Carton Each



**CONTACT SALES REPRESENTATIVE FOR SIZES, QUANTITIES OR MORE INFO!** 



SAFETY









#### **CIDEX®**

#### **OPA SOLUTION**

- Glutaraldehyde-free (0.55% ortho-phthalaldehyde) high-level disinfecting solution effective at room temperature
- · Effective on various endoscopes and instruments
- · Requires no activation or mixing
- · Reusable for up to 14 days, 2-year shelf life

**JOH 20390** 1GL Bottle Each

#### **ACTIVATED DIALDEHYDE SOLUTION**

- · 2.4% alkaline glutaraldehyde solution provides superior microbicidal and anticorrosion properties
- · Works quickly to disinfect metal and plastic devices
- Reduces costs linked to instrument damage and repair
- Activated solution offers up to 14-day reuse life

**JOH 2266** 4.7L Bottle Each

#### **HYDROGEN PEROXIDE 3%**

- Topical solution USP, first aid antiseptic and oral debriding agent
- · Treats and reduces risk of infection in minor cuts, scrapes, burns and abrasions
- · Removes mucus, phlegm or other secretions associated with occasional sore mouth
- · May be used as a gargle or a rinse
- · Keep tightly closed at room temperature





- · Convenient, ready-to-use, intermediate-level surface disinfectant and non-aerosol decontaminant
- · Effective against TB, HCV, HBV, viruses, bacteria (including MRSA and VRE) and fungi
- Safe for all facility spaces including operating rooms, isolation rooms, patient care areas and laboratories

MET 133324 24oz Spray Bottle 12/Case

#### **CAVICIDE®**

- · All-in-one disinfectant, cleaner and decontaminant
- · Ideal for environmental and medical device surfaces
- · Effective against TB, MRSA, HIV-1, HCV and HBV
- · Safe for use on hard, non-porous surfaces

MET 131000	1GL Bottle	Each
MET 131008	8oz Spray Bottle	Each
MET 131024	24oz Spray Bottle	Each
MET 131025	2.5GL Bottle	Each

#### **ERADO-SOL®**

- · Fast, effective biological and chemical stain remover
- · Removes most medical and histological stains from hands, glassware, lab coats, equipment and surfaces
- · Biodegradable, pH neutral and phosphate-free

CAM 20716	16oz Bottle	Each
CAM 20732	32oz Bottle	Each
CAM 2071G	1GL Bottle	Each









#### **ACTIVATE® BLEACH SPRAY**

- · Dual-liquid sprayer automatically combines bleach and water to create a fresh 10:1 diluted solution to disinfect
- · Special lemon-scented bleach formula cleans, spreads, and adheres to surfaces better than ordinary bleach
- · Bleach formula is highly effective against a broad spectrum of pathogens, including C. difficile spores
- Use on nonporous surfaces in hospitals and laboratories







TAK 77062	Includes: (2) Sprayer Heads, (6) 11oz Bleach Cartridges, (2) 11oz Water Cartridges	Kit
TAK 77012	11oz Bleach Cartridge Refills	12/Pack
TAK 77006	(6) Sprayers Heads, (6) 11oz Water Cartridges	6/Pack

#### LYSOL®

· Multipurpose, hospital-grade cleaner cleans, disinfects and deodorizes on hard, nonporous surfaces

#### LYSOL® QUARTERNARY DISINFECTANT CLEANER

- · Concentrated formula with a neutral pH in dilution
- · Effective against HIV-1, staph, strep, pseudomonas and salmonella bacteria
- For floors, walls, porcelain, plastic and metal surfaces

#### LYSOL® DISINFECTANT SPRAY

- · Kills 99.99% of germs of viruses, bacteria and fungi
- · Effective against over 50 microorganisms
- · Control valve produces a wetter spray that enhances surface coverage
- · Ideal for showers, sinks, door handles, waiting rooms, sick rooms, laboratories and around toilet areas

SUL 95029	19oz Spray Can	12/Case
SUL 74983	1GL Bottle	4/Case

#### CLOROX® GERMICIDAL BLEACH

- · Pre-mixed formulated disinfectant kills a broad spectrum of microorganisms including norovirus and C. difficile
- · Excellent in removing molds, mildews and stains
- · Anti-corrosive ingredients are compatible with most hard, nonporous surfaces
- · Meets OSHA bloodborne pathogen standards
- EPA REGISTERED & CDC RECOMMENDED

CLO 30966	Ultra Clorox® Bleach, 121oz	3/Case
CLO 68967	22oz Trigger Spray Bottle	Each
CLO 68970	32oz Trigger Spray Bottle	Each








- WITH PREPAID ANALYSIS
- Can be worn as a badge near the breathing zone to measure personal exposure, or placed in a room to measure area concentration
- Air diffuses through a micro-porous membrane and collects vapor on special prepared adsorbents
- Meet and exceed OSHA requirements
- No need to hire expensive consultants to for testing Indoor Air Quality or Area Monitoring
- SHOULD BE SENT WITHIN TWO WEEKS (14 DAYS) AFTER USE FOR ANALYSIS

ACS F50	Formaldehyde	Each
ACS G10	Glutaraldehyde	Each
ACS OV115	Xylene	Each
ACS OV35	Ethyl Alcohol	Each
ACS OV57	Isopropyl Alcohol	Each
ACS OV60	Methyl Alcohol	5/Case
ACS OV85	Phenol	Each

#### **HOW IT WORKS:**





# TN60 WITH UV-C DISINFECTION CRYOSTAT

#### item # TNR TN60

Dimensions (WxDxH): 26¾" x 23-5/8" x 43¾"

Weight (w/o accessories): 280 lbs Electrical Parameters: 115V; 50/60 Hz



#### **FIRST AID KITS**

- Contain a variety of first aid supplies suited to the size and specific hazards of your work environment
- · Kits come in sturdy plastic, wall mountable container
- Meet and exceed ANSI Z308.1 and OSHA 1910.151 standards (except NSC FAK01)
- · CALL SALES REP FOR DETAILS & OTHER SIZES!

NSC FAK01	Personal	Each
NSC FAKLPEK25	Economy, 25 People	Each
NSC FAKLPEK50	Economy, 50 People	Each
NSC FAKXXV	Deluxe, 25 People	Each
NSC FAKXXL	Deluxe, 50 People	Each
NSC FAKTRMK1	Trauma Kit w/Case	Each

#### **FIRST AID CABINETS**

- · Metal cabinets are wall mountable
- · Meets and exceeds ANSI and OSHA standards

#### **CALL SALES REP FOR DETAILS & OTHER SIZES!**

NSC FAC4	4-Shelf, Deluxe, Large	Each
NSC FAKLXXV	3-Shelf, 75 People	Each

#### **FIRST AID REFILLS**

 Completely replace exhausted or expired first aid supplies originally included in their respective kits

#### **CALL SALES REP FOR DETAILS & OTHER SIZES!**

NSC FAC4R First Aid Refill (FAC4) Each

#### NSC FAKTRMK1

Trauma Kit Includes: Skin Closures, Gauze Bandages, Gauze Rolls, Gauze Pads, Tape, Finger Splint, Eye Wash, Bite Stick, Scissors, Forceps, Burn Spray, Lancet, Thermometer, Rescue Blanket, Triangular Bandage, Absorbent Compress, Ammonia Inhalant, Eye Pads, Insect Sting Wipes, Tourniquet, Biohazard Bag, Burn Cream, Elastic Bandages, Eye Cup, Cold Packs, Alcohol Wipes, Non-Stick Pads, Antiseptic Towelettes, Antibiotic Cream, Iodine Wipes, Hydrocortisone Cream, Water Jel Burn Dressing, Nitrile Gloves, First Aid Facts Guide

#### NSC FAKXXV

- 16 Certi-Strip Bandages (1" x 3")
- 4 Certi-Gauze Pads (3" x 3")
- 10 Certi-Sporyn Antibiotic Cream (1g)
- 10 Antiseptic BZK Towelettes (5" x 7")
- 6 Hand Sanitizer (1g)
- 1 Triangular Bandage (40" x 40" x 56")
- 2 Pair of Nitrile Gloves
- 1 Eye Wash (10z)
- 1 Water Gel Burn Dressing (4" x 4")

- 10 Certi-Burn Cream (1g)
- 1 Certi-Tape (1/2" x 2.5 yds)
- 4 Eye Pads w/Tape
- 1 Certi-Gauze Roll (2" x 6 yds)
- 2 Trauma Pad/Compress (5" x 9")
- 1 Cold Pack, Certi-Cool (Junior)
- 1 CPRotector Plastic Bag
- 1 Scissors w/Red Handle
- 1 First Aid Facts Guide Card







#### **NEUTRA-FORM®**

- Spill control agent quickly diminishes harmful vapors, and neutralizes hazardous aldehydes
- Absorbs large spill of 10% formalin or 2% glutaraldehyde

**SCG 4304** 1.4kg Bottle 2/Case

#### **NEUTRA-WIPES®**

- · Rapidly absorbs and neutralizes small formalin spills
- Wipes (12" x 121/2") come in convenient pop-up dispenser

SCG 4278 60 Wipes Each

#### **DID YOU KNOW?**

Neutra-Pads® neutralize 10% formalin while reducing vapors during grossing, keeping work areas clean and dry. The blue surface provides an excellent contrast to tissues.

#### **NEUTRALEX®**

- Neutralizes 4 liters (1GL) of 10% formalin or 8 liters (2GL) of 2% glutaraldehyde in 15 minutes
- · Safe and effective drain disposal with no clogging solids
- Final neutralization is less than 1 ppm

#### STARTER PACK

 Includes: (4) Neutralex\* Pouches, (1) Reaction Container, (10) 8" x 10" Neutra-Pads\*, (1) Neutra-Form\* Container, (60) Neutra-Wipes\*, (1) Aldehyde Test Kit

#### **REACTION CONTAINER**

- · Waste collection and formalin treatment container
- Capacity: 9.5L (2.5GL)

SCG 4046	Starter Pack Eac	
SCG 4047	0.75kg Pouch	16/Case
SCG 4278	Reaction Container	Each

#### **HYDEAWAY**<sup>M</sup>

- Neutralizes formaldehydes, formalin and glutaraldehydes into a non-toxic, non-hazardous end product
- · Liquid allows for drain disposal
- · Crystals solidify waste for trash, excellent for spills
- Environmentally friendly, makes for a safer workspace
- Highest Neutralization Ratios: Liquid 1:15, Crystals 1:30

DEC 2201	1GL Bottle (Liquid)	Each
DEC 2224	32oz Bottle (Crystals)	Each











#### FORMALEX® "GREEN"

- Liquid control agent for safe drain disposal of formalin, glutaraldehyde and other aldehyde waste solutions
- Improved formula eliminates any sludgy residue and does not require any pH adjustments
- · Neutralizes into a nonhazardous, non-toxic solution
- · Contains no phosphates and is non-vapor forming
- · Spray bottle cleans small countertop or bench spills
- · Environmentally friendly to personnel and laboratories

ABS FXG32E	Spray Bottle, 32oz	Each
ABS FXG01	Cubi-Tainer, 1GL, Spigot	Each
ABS FXG250	Cubi-Tainer, 2.5GL, Spigot	Each
ABS FXG05	Cubi-Tainer, 5GL, Spigot	Each

#### **FORM-ZERO®**

#### FORMALDEHYDE REDUCING AGENT

- Unique, powdered reagent blend safely and quickly neutralizes formaldehyde and glutaraldehyde waste
- From formalin to drain in 20 minutes
- Simply add one bottle of Form-Zero to each gallon of 10% formalin waste or 4% glutaraldehyde waste
- · No pH adjustment required

ABS FZ55012	500cc Bottle	12/Case
ABS FZ55032	500cc Bottle	32/Case

#### **SAFECUBE DISPENSER™**

- 5-gallon cubitainer hangs over counter edge
- · Control drips, drops and accidental spills
- Liquid detection alarm immediately signals working personnel of any accidental spills
- · Dispenser tray balances container during filling
- Enhances regulatory compliance

ABS SCD100A 5GL Capacity Each

#### SASCO SOLIDIFIERS/NEUTRALIZERS

- Acid Handler: flowable powder designed to raise the pH of an acid solution and solidify the liquid providing for a safe nonhazardous clean up
- Base Control: flowable powder designed to lower the pH of a caustic solution and solidify the liquid providing for a safe nonhazardous clean up
- Solvent Handler: free flowing granular material reduces the flammability and solidify flammable hydrocarbon solutions providing a safe method of clean up
- **Polyform-F:** very unique granular material designed to destroy formaldehyde, formaldehyde-based fixatives, glutaraldehyde and other aldehyde solutions
- Chlorine Control: specially formulated powder designed for the destruction and control of various types of liquid chlorine spills and leaks

ABS AH32	Acid Handler™, 32oz	Each
ABS BC32	Base Control™, 32oz	Each
ABS PF32	Polyform-F™, 32oz	Each
ABS SH32	Solvent Handler™, 32oz	Each





#### SASCO SPILL CONTROL PRODUCTS

- A = ACID HANDLER
- B = BASE CONTROL
- P = POLYFORM-F
- S = SOLVENT HANDLER
- C = CHLORINE CONTROL

#### SASCO ACID HANDLER™ | USE ON THE FOLLOWING ACIDS:

Acetic Acid • Acetic Anhydride • Acetyl Chloride • Sulfonic Acid Aluminum Chloride • Chlorosulfonic Acid • Citric Acid Formic Acid • Glacial Acetic Acid • Hydrochloric Acid Hydrofluosilic Acid • Muriatic Acid • Nitric Acid • Sulfuric Acid Perchloric Acid • Phosphoric Acid • Phosphoric Anhydride Phosphorous Pentoxide • Phosphorous Trichloride Dodecylbensylsulfonic Acid • 54% Hydrofluoric Acid Solution Chromic Acid Solutions • Chromium Waste Hazard

#### DO NOT USE ON THE FOLLOWING:

Chlorine • Oxidizers • Sodium Amide • Picric Acid Hydrogen Peroxide • Iodic Acid • Sulfurous Fluoride Antimony Pentrafluoride ("Super-Acids") • Concentrated Hydrofluoric Acid

#### SASCO BASE CONTROL™ | USE ON THE FOLLOWING BASES:

Triethylalamine • Monoethanolamine • Morpholine
Most Alkali Detergents • Anhydrous Ammonia
Sodium Metasilicate Solution
Sodium Hydroxide = Caustic Soda
Potassium Hydroxide = Caustic Potash
Ammonium Hydroxide = Aqua Ammonia

#### DO NOT USE ON THE FOLLOWING:

Chlorine • Oxidizers • Sodium Hypochloride Hydrogen Peroxide • Sulfurous Fluoride • Sodium Amide

#### **GENERAL SPILL RESPONSE KITS**

 Includes: Safety Goggles, Gloves, Waste Bags, Broom and Pan for Safe Waste Collection

ABS SRKGP3204	Choose Any <u>FOUR</u> Spill Control Products from Product List	Kit
ABS SRKGP3202	Choose Any <b>TWO</b> Spill Control Produ	ucts Kit

from Product List

**Note:** Limited to <u>ONE</u> Solvent Handler<sup>™</sup> per kit. Each additional Solvent Handler<sup>®</sup> per kit will be an additional charge.

#### SASCO POLYFORM-F™ | USE ON THE FOLLOWING:

Formaldehyde • Formaldehyde-Based Embalming Solution 10% Neutral Buffered Formalin • Glutaraldehyde Solution; e.g. (Cidex\*, Cidex-OPA\*, Metracide\*, Wavacide\*, OmniCide\*)

**Bouins Fixative:** Add Polyform-F™ to destroy formaldehyde. Neutralize pH with bicarbonate or caustic soda. This will convert picric acid to sodium picrate, which can be disposed of safely.

**B-5 Fixative:** Contains heavy metals, which constitutes a hazardous waste.

#### SASCO SOLVENT HANDLER™ | USE ON THE FOLLOWING:

Alcohols • Chloroform • Methyl Ethyl Ketone • 111 Trichlorethane N,N-Dimethylformamide • Most Hydrocarbons • Xylene Acetone • Ethyl Acetate • Gasoline • Diesel • Most Halogens

DO NOT USE in large volume in enclosed area.

**BE EXTREMELY CAREFUL** using **SOLVENT HANDLER**™ on solvents with low auto-ignition temperature like "Nitromethane."

The lists given herein are general and do not necessarily include all the materials Acid Handler<sup>™</sup>, Polyform-F<sup>™</sup>, Base Control<sup>™</sup>, and Solvent Handler<sup>™</sup> can or cannot be used on. If you would like an exotic species tested, or have questions as to application of the products, contact the S&S Company of Georgia, Inc. or American Bio-Safety, Inc. for further information. They will try applying their products to your substance and determine a suitable procedure for spill situations.





#### **MERCEDES SCIENTIFIC®**

#### **Q-SORB® FORMALDEHYDE ABSORBENT PADS**

- Designed to absorb and neutralize formaldehyde vapors in a laboratory environment
- · Light blue colored sheets provide excellent contrast for tissue grossing
- · Reduce harmful vapors while transferring, dissecting, storing or transporting tissue in formalin containers
- · Tear-resistant, durable and disposable
- · Useful for areas prone to leakage

MER AP0811	8" x 11"	25/Pack
MER AP1115	11" x 15"	25/Pack
MER AP1520	15" x 20"	25/Pack





#### **ABSORBENT UNDERPADS**

- 3-ply, highly absorbent tissue-filled or fluff-and-tissue-filled center with a waterproof polyethylene backing
- Soft, nonwoven material protects non-irritating skin from wetness for enhanced comfort

TID 16650	17" × 24"	300/Case
TID 16651	24" x 36"	150/Case
TID 16652	23" x 24"	200/Case
TID 16653	12" x 17"	500/Case



#### **FORMALIN PADS**

- · Super absorbent pads neutralize and eliminate formalin, formaldehyde, glutaraldehyde and other aldehydes
- Convert the aldehydes to a nonhazardous polymer to neutralize harmful vapors once the solution is absorbed

ABS FAN55	Mini Pad, 11" x 11"	6/Case
ABS FAN614	Small Pad, 6" x 14"	48/Case
ABS FAN825	Large Pad, 8" x 25"	30/Case



#### **QUANTOFIX® FORMALDEHYDE**

Test sticks and reagents for the semi-quantitative determination of formaldehyde in solutions. Gradation of scale: 0-10-20-40-60-100-200 mg/l HCHO beige to blue-violet. Other aldehydes such as e. g. acetaldehyde or glutardialdehyde only react at higher concentrations, resulting in a discoloration which deviates from the colour scale. Acetone does not react. Strong oxidizing and reducing agents cause low results.

CSS 91328 100 Test Strips + Reagent Kit



## ACID-NEUTRAL-ALKALINE PH TEST PAPER

- · Includes color comparison chart
- Simply immerse strip into solution under test and compare to chart
- Determines if a solution is acidic (red), neutral (yellow-green), basic (green) or alkaline (blue)

**FIS S63129** 100 Strips/Pack 12/Pack



#### **MQUANT® PH-INDICATOR STRIPS**

- Different pH ranges from the universal range 0 to 14
- Contain special indicator dyes that are covalently bound to the cellulose of the indicator paper
- Suited for all media in environmental analysis and industrial in-process controls
- · Non-bleeding and highly accurate even in weak buffers
- pH value is measured semiquantitatively by visual comparison of the reaction zone on the test strip with the fields of the color scale

EMD 1095400001	0-2.5	100/Pack
EMD 1095310001	0-6.0	100/Pack
EMD 1095350001	0-14.0	100/Pack
EMD 1095410001	2.5-4.5	100/Pack
EMD 1095420001	4.0-7.0	100/Pack
EMD 1095840001	2.0-9.0	100/Pack
EMD 1095330007	5.0-10.0	600/Pack
EMD 1095430001	6.5-10.0	100/Pack
EMD 1095320001	7.5-14.0	100/Pack
EMD 1095450001	11.0-13.0	100/Pack

WEB MercedesScientific.com



153

#### **OUTFRONT REVERSIBLE MAT**

- Flexible scrapers and quick-dry drainage ports remove dirt and moisture in more moderate weather while stiff ridges scrape heavier debris and packed snow
- MOST VERSATILE SCRAPER MAT AVAILABLE
- · Outdoor mat accommodates all weather conditions
- · Provides an all-season line of defense for entry ways
- UV RESISTANT

Overall Thickness: 7/16" (11mm)
 SIZE AVAILABLE: 3' x 6'

· COLOR: Black

**TEN 227** Outdoor Reversible Mat Each

#### **TILE-TOP AM SPONGECOAT MATS**

- Doubles down on combating fatigue even while providing a top to bottom anti-microbial solution
- Ideal for laboratories, pharmacies, biotechnology and pharmaceutical manufacturing facilities
- Absorbs repeated compression and bounces back to its original form better than any PVC sponge alternative
- THE ONLY FULLY CUSTOMIZABLE TILE-TOP MATTING PRODUCT AVAILABLE
- · Machine cut and beveled edges prevent tripping
- Tough but attractive abrasion and chemical resistant PVC surface
- Anti-microbial top surface provides top level defense and is backed with an anti-microbial Nitricell sponge
- Overall Thickness: 1/2" (Standard), 7/8" (UltraSoft)
- SIZES AVAILABLE: 2' x 3' | 2' x 60' | 3' x 5' | 3' x 10' 3' x 60' | 4' x 60' | Custom
- · COLORS: Black, Blue, Charcoal, Gray

TEN 419UltraSoft Tile-Top MatsEachTEN 420Standard Tile-Top MatsEach

#### **TILE-TOP SELECT MATS**

- Mid-range offering within the popular family of marbleized antifatigue mats
- Anti-microbial marble top surface is a chemical-resistant and durable PVC that is lightly textured, supple to the touch and easy-to-clean
- Special "Flex-Link" sponge base or closed-cell PVC foam provides solid, stable performance to combat fatigue
- · Resistant to many common chemicals and cleaners
- Overall Thickness: 1/2" (Standard), 7/8" (UltraSoft)
- SIZES AVAILABLE: 2' x 3' | 2' x 60' | 3' x 5' | 3' x 60' 4' x 60' | Custom
- COLORS: Black, Blue, Charcoal, Gray

TEN 494 Tile-Top Select Mats Each


#### 227 MAT SERIES






419/420 MAT SERIES









494 MAT SERIES













WEB MercedesScientific.com

#### **COMFORTPRO MATS**

- · Special "Flex-Link" PVC sponge base is a budget friendly fatigue fighting solution and is easy-to-clean
- · Wide corrugated pattern provides good traction
- Best solution when a corrugated surface is a must
- Overall Thickness: 1/2" (Standard), 7/8" (UltraSoft)
- **SIZES AVAILABLE:** 2' x 3' | 3' x 4' | 3' x 5' | 2' x 75' 3' x 75' | 4' x 75' | Custom
- · COLOR: Black

**TEN 433** ComfortPro Mats Each

#### **AUTOCLAVABLE MAT**


- · Ideal for sterile environments such as hospital operating rooms, surgery centers, laboratories, pharmaceutical and life sciences
- Features 100% latex-free, non-porous rubber compound
- · Promotes health and wellness for professionals working in sterile applications
- Rubber surface domes (flexible and resilient underfoot) provide relief during long shifts and extended standing
- · 100% anti-microbial and USP 797 compliant
- Overall Thickness: 1/2" (13mm) • SIZE AVAILABLE: 2' x 3'
- · COLOR: Gray


**TEN 580** Autoclavable Mat Each

#### **ABRASIVE COATED KUSHION WALK SLOTTED MATS**

- · Industrial tough, and abrasion and highly cut resistant
- · Durable textured surface guards against slipping
- Unique GRITSHIELD surface is specifically designed for high COF (coefficient of friction) in the slickest locations
- · Effective standing solution in areas requiring high traction to combat slips and falls
- · Combats the dangerous nuisance of slips and falls in today's workplace
- · Excellent choice for covering grated walkways
- Overall Thickness: 3/8" (9mm)
- SIZES AVAILABLE: 2' x 60' | 3' x 60' | Custom
- **COLOR:** Black

**TEN 480** Abrasive Coated Kushion Mats Roll











#### **ULTRA-TRED ARMORCOTE MATS**

- More durable urethane coating for greater abrasion resistance than standard blown PVC sponge products
- Tough outer skin increases service life, especially in applications where abrasion can cause the mat to fail
- · Economic comfort for medium traffic areas
- Overall Thickness: 3/8" (9mm)
- SIZES AVAILABLE: 2' x 3' | 2' x 60' | 3' x 5' | 3' x 12' 3' x 60' | 4' x 60' | Custom
- · COLOR: Black w/Yellow Borders

TEN 440 Ultra-Tred ArmorCote Mats Each

#### **REJUVENATOR MATS**

- Large modular tiles interlock on all 4 sides to customize on-site virtually any flooring configuration
- · Deep positive lock keeps seams secure
- · 100% urethane with thick outer wear layer
- Ability to customize in ONE FOOT (30cm) increments
- Choose from beveled  $\ensuremath{\mathbf{BLACK}}$  or  $\ensuremath{\mathbf{YELLOW}}$  ramp edging
- Overall Thickness: 5/8" (16mm)
- SIZES AVAILABLE: 1' x 3' | 3' x 3'

· COLOR: Black

TEN 502 Rejuvenator Square Mats Each

#### **CLEAN ROOM MATS**

- Trap impurities from shoe soles for dust and dirt control
- · Sheets of "tacky" film guards against contamination
- Tear-resistant sheets are numbered to indicate when a new mat is needed by simply peeling off the old sheet
- Use on low-pile carpet, tile or concrete surfaces
- · Adhesive backing holds mat securely in place
- SIZES AVAILABLE: 18" x 36" | 18" x 45" | 24" x 45" | 36" x 45" | 2' x 3' | 3' x 5'
- COLORS: Blue, White

TEN 095 Clean Room Mats 4/Case

#### **COMFORT2GO MAT**

- · Long-lasting, portable, versatile, anti-fatigue mat
- · Provides standing and kneeling support
- · Built-in handle for easy carrying and storing
- Double thick Nitricell\* sponge base infused with nitrile rubber for exceptional comfort
- Overall Thickness: 1/8" (22mm)
- SIZE AVAILABLE: 17" x 30"
- · COLOR: Black

TEN 518 Comfort2Go Mat Each

## WEIGHTS & MEASURES



#### **ENCLOSED CHAMBER THERMOMETERS**

ITEM #	RANGE	DIVISIONS	APPLICATION	HEIGHT	BOTTLE MEDIA	
THM ACCR0201S	-5° to 15°C	0.5°C	Refrigerator	125mm	Non-Mercury Spirit Filled Bio-Safe Glycol/Water, 30mL	
THM ACCR20001S	-2° to 10°C	0.1°C	Refrigerator	240mm	Non-Mercury Spirit Filled Bio-Safe Glycol/Water, 60mL	
THM ACCAM0401S	10° to 30°C	0.5°C	Ambient/Room	140mm	Non-Mercury Spirit Filled, Bio-Safe Glycol/Water, 30mL	

#### **TISSUE FLOTATION BATH THERMOMETERS**

ITEM #	RANGE	DIVISIONS	IMMERSION	LENGTH	BOTTLE MEDIA
THM ACC5351S	-15° to 105°C	0.1°C	Total	145mm	Non-Mercury Red Spirit
THM ACC5351SSC	-15° to 105°C	0.2°F	Total	145mm	Non-Mercury Red Spirit, Safety Coated

#### **OVEN THERMOMETERS**

ITEM #	RANGE	DIVISIONS	IMMERSION	LENGTH	BOTTLE MEDIA
THM ACCOV070S	20° to 130°C	1.0°C	-	135mm	Non-Mercury Spirit Filled, Sand, 30mL
THM ACCOV700S	50° to 115°C	0.5°C	-	240mm	Non-Mercury Spirit Filled, Sand, 60mL
THM ACCOV080S	35° to 200°C	1.0°C	-	145mm	Non-Mercury Spirit Filled, Sand, 100mL
THM ACC637S	0° to 200°C	1.0°C	76mm	280mm	Non-Mercury Spirit Filled
THM ACC637SSC	0° to 200°C	1.0°C	76mm	280mm	Non-Mercury Spirit Filled, Safety Coated
THM ACC639S	0° to 200°C	1.0°C	125mm	405mm	Non-Mercury Spirit Filled
THM ACC639SSC	0° to 200°C	1.0°C	125mm	405mm	Non-Mercury Spirit Filled, Safety Coated

#### **LARGE DIGIT - TRIPLE DISPLAY ALARM**

#### **Current-Min-Max Temperatures**

The Triple Display Digital Thermometer simultaneously shows the current, minimum, and maximum temperatures and updates continuously. An audible alarm will sound when the temperature rises above or goes below the set temperatures. The minimum/maximum feature monitors critical temperature for any period of time, such as overnight or weekend. Reset buttons on front panel.

Each thermometer is certified and traceable to NIST. Each thermometer is serialized. The thermometer probe is housed in an unbreakable plastic bottle or vial filled with a Bio-Safe Glycol / Water or Sand buffering media which maintains the thermometer reading when the instrument door is opened.

The thermometer and the bottle housing the probe have dual magnets allowing attachment to doors or walls of the instrument.

Other features include a flip-open stand, wall mounting key hole, and supplied with a AAA alkaline battery. The instrument is °C and °F switchable.



ITEM #	APPLICATION	RANGE	BOTTLE MEDIA	CERTIFIED TEMPERATURE
THM ACC895DIG	General Laboratory	-50° to 70°C		
THM ACC895FRE	Freezers	-50° to 70°C	Bio-Safe Glycol/Water	-20.0°C
THM ACC895REF	Refrigerators	-50° to 70°C	Bio-Safe Glycol/Water	+4.0°C
THM ACC895INC	Incubators	-50° to 70°C	Bio-Safe Glycol/Water	+37.0°C
ТНМ АСС895АМВ	Ambient / Room	-50° to 70°C	Bio-Safe Glycol/Water	+22.0°C
THM ACC895FREV	Freezers	-50° to 70°C	Bio-Safe Glycol/Water	-20.0°C
THM ACC895REFV	Refrigerators	-50° to 70°C	Bio-Safe Glycol/Water	+4.0°C
THM ACC895BLK	Block Heaters	-50° to 70°C	Mineral Oil 1ml	+37.0°C
THM ACC895WB	Water Baths	-50° to 70°C	Sand 30ml	+37.0°C


#### PARAFFIN MELTING POINT THERMOMETERS

ITEM #	RANGE	DIVISIONS	IMMERSION	LENGTH	BOTTLE MEDIA	
THM ACC014C	38° to 82°C	0.1°C	79mm	380mm	Mercury Filled	
THM ACC014F	100° to 180°F	0.2°F	79mm	380mm	Mercury Filled	













#### **JUMBO BLACK TIMER**

- Simultaneous programming of two different activities from 24 hours to 1 second with 0.01% accuracy
- · Second display mode shows time of day and a stopwatch
- LCD display can be read up to 30 feet away

CCM 5022CC Jumbo Black Timer Each

#### **MULTI-COLORED TIMER**

- 3-bright colorful supplied frames permit personalization
- Unique memory returns channel to previously programmed setting at touch of a key
- Counts up during alarming indicating elapsed time since alarming, 12/24 hour clock
- · 4 channels count up in stopwatch mode

CCM 5012CC Multi-Colored Timer Each

#### THREE-CHANNEL ALARM TIMER

- Triple display shows 3 different times at once, tests can be started independently or simultaneously
- Adjustable alarm volume can be set to sound for 60 seconds or to sound continuously until turned off
- Loud decibel for noisy areas and low for quiet areas
- Features: Count Up/Down, Time In/Out, Memories, Alarm, Clock, Adjustable Alarm Duration/Volume

CCM 5000CC 3-Channel Alarm Timer Each

WEB MercedesScientific.com

#### **MONITORING HYGROMETER**

- · Simultaneously displays time, temperature and humidity
- Internal sensors are ideal for use in hoods, storerooms, incubators, drying chambers and environmental cabinets
- Memory recalls temperature and humidity readings

CCM 4040 Monitoring Hygrometer Each

#### **FOUR-CHANNEL ALARM TIMER**

- · Massive high display allows viewing from across the lab
- Remarkable memory returns display to previously programmed countdown time at the touch of a button
- Extra-loud alarm for 1 min or can be silenced manually
- Features: Count Up/Down, Time In/Out, Memories, Alarm, 12/24 Hour Clock

CCM 5004 4-Channel Alarm Timer Each

#### **DIGITAL TIMER PLUS™**

- Six-digit LCD display shows time remaining in hours, minutes, seconds with quartz-crystal accuracy
- Alarm sounds for 60 seconds or may be silenced manually when zero is reached
- Handiest timer designed to clip to any surface such as lab coat, belt, notebook or clipboard
- Features: Count Up/Down, Time In/Out, Keychain, Alarm

CCM 5045CC Digital Timer Plus Each

EACH UNIT INCLUDES TRACEABLE® CERTIFICATE SUPPORTED BY NIST













#### **WEIGHING BOATS**

- Flexible polystyrene boats feature smooth surfaces to provide accurate pour-outs with minimal sample loss
- Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bases
- · Sloping sides and flat bottoms resist tipping
- · Rounded corners simplify transfer

#### **HEXAGONAL (WHITE)**

HEA HD14251A	Small	500/Pack	
HEA HD14251B	Medium	500/Pack	
HEA HD14251C	Large	500/Pack	
HEA HD14251D	X-Large	500/Pack	

#### STANDARD

- AVAILABLE IN REGULAR & ANTISTATIC
- SIZES AVAILABLE: Small, Medium, Large

HEA HD1420	White	500/Pack
HEA HD1421	Blue	500/Pack
HEA HD1423	Black	500/Pack

#### **ALUMINUM**

- Sturdy, lightweight boats feature smooth, flat bottom for sample handling and flanged, crimped wall for rigidity
- Finger tab handle can be etched for identification

HEA HD14521A	Small	500/Pack
HEA HD14521B	Medium	500/Pack
HEA HD14521C	Large	500/Pack

#### **ANTISTATIC WEIGHING DISHES**


- Made from a special grade of polystyrene
- Resist acids, aqueous solutions, alcohols and bases
- · Safe, contaminant-free, biologically inert and economical
- · Ideal for weighing, dispensing or storing
- Smooth surface for accurate pours with tiny sample loss
- · Non-reactive with most substances
- · Stackable for easy storage
- · Excellent for handling solids or liquids during weighing
- Facilitate dispensing of powdered and liquid materials
- Suitable for weighing static affected samples
- · Withstands temperatures up to 80°C


#### **POUR BOATS**

GLO 3623	20mL, Small	250/Pack
GLO 3624	140mL, Medium	250/Pack
GLO 3625	270mL, Large	250/Pack
SQUARE		
GLO 3620	20mL, Small	500/Pack
GLO 3621	100mL, Medium	500/Pack
GLO 3622	330mL, Large	500/Pack
HEXAGONAL		
GLO 3615	20mL, Small	500/Pack
GLO 3616	50mL, Medium	500/Pack
GLO 3617	200mL, Large	500/Pack
GLO 3618	350mL, X-Large	500/Pack













Aprimproproprogramment of the control of the contro

FIS 09016

#### **RULERS**

#### PAPER

- · White paper ruler in centimeters and inches
- · Features specimen and date fields

#### **PLASTIC**

- · Red plastic ruler in centimeters and inches
- · Cork and stopper size guide included on back
- Temperature conversion scale included on back featuring -20° to +100°C and -10° to +212°F ranges

#### VINY

- Blue vinyl ruler in centimeters and inches
- Dual-scale, 0.062" and 1mm increments

MER 1000	Paper, 8"	250/Pack
MER RULER	Plastic, 6"	100/Box
FIS 09016	Vinvl. 6"	100/Pack







## **Microscope Slide Label Printer**


A complete, ready-to-print microscope slide label printer package that comes with everything you need to print histology grade slide labels — printer, software, labels, and ribbon.

The slide labels and thermal transfer ribbon are optimized for harsh chemicals and will not absorb eosin or hematoxylin coloring.

SlideID printers function with or without LIS. They can be set up for batch printing or with on-demand printing for "laboratory best practices."

#### **FEATURES:**

- Ready-to-print
- Scalable
- · Resistant to stains, solvents and xylene
- · Will not absorb eosin or hematoxylin coloring
- · Function with or without LIS
- Can be setup for batch printing or with on-demand printing




QRL SLIDEIDPRINTERPK	BASIC PACKAGE: SlideID CL-E331 Slide Label Printer, SlideID Slide Printing Software, SlideID Microscope Slide Labels, SlideID Thermal Transfer Ribbon	Each
QRL SLIDEIDPROPK	PRO PACKAGE: Includes All Contents of the Basic Package Above Plus Denso AT30Q Barcode Scanner	Each
QRL SLIDEIDLABEL	SlideID Microscope Slide Labels, 0.875" x 0.875", 1" core (1,750 Per Roll)	Roll
QRL SLIDEIDRIBBON	SlideID Special Resin Ribbon, 1.57" W x 984' (prints 12,000 +\- Slide labels)	Each
QRL SLIDEIDRIBBONZ	SlideID Thermal Transfer Ribbon For ZEBRA® Style Printers, Xylene and solvent resistant, histology grade	Each

#### **SPECIMEN & PATHOLOGY LABELS**

- Pressure-sensitive, non-smudge labels accept written notations in pen, pencil or marker
- Adhere instantly to glass, plastic or waxed surfaces
- · Identify specimens more easily
- · Holds test requisitions to specimen containers

SHA SL3	1" x <sup>21</sup> / <sub>4</sub> "	500/Roll
SHA UPCR9404	2" x 3"	500/Roll
SHA UPCR6052	1" x <sup>2</sup> 1/4"	500/Roll
SHA UPCR6054	1" x <sup>21</sup> / <sub>4</sub> "	500/Roll




#### THERMAL PRINTER LABELS

FOR COMPUTER GENERATED CONSECUTIVE NUMBER PRINTING

SHA THERM755R	7/8" x <sup>11</sup> /4"	3,400/RoII
SHA THERM628R	7/8" x 11/2"	3,400/Roll
SHA THERM123R	1" x <sup>15</sup> /16"	3,000/Roll
SHA THERM125R	1" x 1½"	3,000/Roll
SHA THERM566R	3/4" x 11/2"	3,900/Roll
SHA THERM100R	3/8" x 11/4"	5,000/Roll
SHA THERM617R	5/8" x 1"	4,500/Roll
SHA THERM504R	5/8" x <sup>7</sup> /8"	4,500/Roll
SHA THERM305R	3/4" X 3/4"	3,900/Roll




SHA THERM1213R	<sup>15</sup> / <sub>16</sub> " X <sup>15</sup> / <sub>16</sub> "	12,800/RoII	
SHA THERM1318R	7/8" x 7/8"	12,800/Roll	



#### MICROSCOPE SLIDE FORM LABELS

• RC = Round Corners

SHA SLR15	Standard, <sup>15</sup> / <sub>16</sub> " x <sup>15</sup> / <sub>16</sub> "	1,000/RoII
SHA SLR15RC	Standard, RC, <sup>15</sup> / <sub>16</sub> " x <sup>15</sup> / <sub>16</sub> "	1,000/RoII
SHA SLRP15	Pathology, <sup>15</sup> / <sub>16</sub> " x <sup>15</sup> / <sub>16</sub> "	1,000/RoII
SHA SLRP15RC	Pathology, RC, <sup>15</sup> / <sub>16</sub> " x <sup>15</sup> / <sub>16</sub> "	1,000/RoII
SHA SER16	Standard, <sup>15</sup> / <sub>32</sub> " x <sup>15</sup> / <sub>16</sub> "	1,000/RoII
SHA SER16RC	Standard, RC, <sup>15</sup> / <sub>32</sub> " x <sup>15</sup> / <sub>16</sub> "	1,000/RoII
SHA SERP16	Pathology, <sup>15</sup> / <sub>32</sub> " x <sup>15</sup> / <sub>16</sub> "	1,000/Roll
SHA SERP16RC	Pathology, RC, <sup>15</sup> / <sub>32</sub> " x <sup>15</sup> / <sub>16</sub> "	1,000/Roll





#### LASER PRINTER LABELS

FOR MICROSCOPE SLIDES


- Supplied on 8½" x 11" sheets
- RC = Round Corners
- CALL FOR COLOR OPTIONS!

ITEM #	SIZE, SHEETS	QTY
SHA SLSLASER	<sup>15</sup> / <sub>16</sub> " x <sup>15</sup> / <sub>16</sub> ", 88	5,000/Pack
SHA SLSLASERRC	<sup>15</sup> / <sub>16</sub> " x <sup>15</sup> / <sub>16</sub> ", 80	5,000/Pack
SHA SESLASER	15/32" x 15/16", 168	5,000/Pack









SHA SBH103 PRESSURE-SENSITIVE SHA SBH104 NON-ADHESIVE



SHA SBH101 PRESSURE-SENSITIVE SHA SBH102 NON-ADHESIVE



SHA SBH5

#### **BIOHAZARD WARNING COMMUNICATION LABELS**

ITEM #	DESCRIPTION	SIZE	COLOR (LABEL, INK)	QTY
SHA SBH101	Biohazard Symbol, Permanent Adhesive	8" x 10"	Fluorescent Red, Black	10/Pack
SHA SBH102	Biohazard Symbol, Non-Adhesive, Cardstock	8" x 10"	Fluorescent Red, Black	10/Pack
SHA SBH103	Biohazard, Infectious Agent, Permanent Adhesive	8" x 10"	Fluorescent Red, Black	10/Pack
SHA SBH104	Biohazard, Infectious Agent, Non-Adhesive, Cardstock	8" x 10"	Fluorescent Red, Black	10/Pack
SHA SBH2	Biohazard, Hazard Identity, Non-Adhesive, Cardstock	8" x 10"	Fluorescent Red, Black	10/Pack
SHA SBH3	Biohazard, Hazard Identity, Permanent Adhesive	8" x 10"	Fluorescent Red, Black	10/Pack
SHA SBH1	Biohazard Symbol	3" x 21/4"	Fluorescent Red, Black	500/Roll
SHA SBH4	Biohazard, Hazard Identity	1" x 3"	Fluorescent Red, Black	500/Roll
SHA SBH5	Biohazard, Hazard Identity	3½" x 2½"	Fluorescent Red, Black	100/Pack
SHA SBH60	Biohazard, Etiologic Agents, Krome Material	2" x 4"	White, Red	100/Pack
SHA SBH61	Biohazard, Clinical Specimens Biological Products	2" x 4"	Fluorescent Red, Black	100/Pack
SHA SBH62	Biohazard, Contaminated, Blood/Body Fluid Presence	1½" x 3"	Fluorescent Red, Black	500/Roll
SHA SBH8D	Biohazard Symbol	1" Diameter	Fluorescent Red, Black	500/Roll
SHA SBH8L	Biohazard Symbol	<sup>5</sup> / <sub>16</sub> " × 1½"	Fluorescent Red, Black	500/Roll



SHA UPCR9700



SHA UPCR9701



SHA SB372H



SHA SB375H



SHA UPCR9250



SHA SB370H



SHA SB371H



SHA SB373H

CAUTION ACID ADDED

SHA UPCR6028

**CAUTION: ACID**Prevent Skin Contact

SHA UPCR6029

CAUTION CONTAINS CAUSTIC OR TOXIC AGENT

SHA UPCR6031



SHA UPCR7006



SHA UPCR6032



SHA UPCR6034

CAUTION CARCINOGEN

SHA UPCR6030



SHA UPCR7007



SHA UPCR1300



SHA UPCR1307

CLEAN: NO GLOVES

SHA UPCR3006



SHA UPCR7008

#### SPECIALTY CAUTION LABELS

SHA SB370H	Caution Caustic, 1" x 21/4"	SHA UPCR6030	Caution Carcinogen, 1" x 21/4"	500/Roll
SHA SB371H	Caution Contaminated, 1" x 21/4"	SHA UPCR6031	Caustic or Toxic Agent, 1" x 21/4"	500/Roll
SHA SB372H	Caution Contains Acid, 1" x 21/4"	SHA UPCR6032	Caution Flammable, 1" x 21/4"	500/Roll
SHA SB373H	Caution Poison, 1" x 21/4"	SHA UPCR6034	Poison, 1" x 21/4"	500/Roll
SHA SB375H	Hazardous Materials, 1" x 21/4"	SHA UPCR7006	Antineoplastic Material, 1" x 3"	500/Roll
SHA UPCR9250	Danger Acid, 1¾" x 2¾"	SHA UPCR7007	Corrosive, 1" x 3"	500/Roll
SHA UPCR9700	Caution This Bottle, 2½" x 4"	SHA UPCR7008	Flammable, 1" x 3"	500/Roll
SHA UPCR9701	Danger Flammable, 2½" x 4"	SHA UPCR1300	Poison, <sup>5</sup> / <sub>16</sub> " x 11/4"	500/Roll
SHA UPCR6028	Caution Acid Added, 1" x 21/4"	SHA UPCR1307	Carcinogen, 5/16" x 11/4"	500/Roll
SHA UPCR6029	Prevent Skin Contact, 1" x 21/4"	SHA UPCR3006	Clean: No Gloves, 3/4" x 11/2"	500/Roll

## 10% FORMALIN

SHA SB374H

#### 10% FORMALIN

CAUTION CONTAINS FORMALDEHYDE

VISE CANCER REPEATED OR PROLONGED EXPO-RE INCREASES THE RISK. CAT. NO. SB-381H

SHA SB381H

#### CAUTION

CONTAINS FORMALDEHYDE

SHA SB385H



SHA UPCR6042



SHA UPCR6043

#### 10% FORMALIN POISON 😣 CAT, NO. SB-382H

SHA SB382H

### 10% FORMALIN @ CAUTION: CONTAINS FORMALDEHYDE

SHA UPCR6046

## ALCOHOL

WARNING: FLAMMABLE LIQUID

CAT, NO. SB-389H

SHA SB389H

#### CAUTION RELEASES FORMALDEHYDE

TOXIC BY INHALATION, IRRITATING TO THE TOJIC BY INMALATOR, BHITIATING TO THE EYES, RESPIRATORY SYSTEM AND SKIN. MAY CAUSE SENSITIZATION BY INHALATION OR SKIN CONTACT, MAY CAUSE CANCER. REPEATED OR PROLONGED EXPOSURE INCREASES THE RISK.

CAT, NO. 58-387H

SHA SB387H



SHA UPCR6044

## DO NOT DRINK

SHA UPCR5009

PAP C	YTOLOGY FIXATIVE
PATIENT	
DOCTOR .	
SOURCE .	
DATE	HOSP. #

SHA UPCR6047

WARNING: FLAMMABLE LIQUID

CAUTION: CONTAINS ETHYL BENZENE, TOKIC BY INHALATION AND IF SWALLOWED, IRRITATING TO EYES. RESPIRATORY SYSTEM, AND SKIN, WASH THOROUGHLY AFTER HANDLING.


CAT, NO. SB-388H

SHA SB388H

#### 10% FORMALIN CAUTION CONTAINS FORMALDEHYDE

TODG BY INHALATION AND IF SWALLOWED, INHITEAT-ING TO THE EYES, RESPIRATORY SYSTEM, AND SKIN MAY CAUSE SENSITIZATION BY INHALATION OR SKIP CAT, NO. SB-380H


SHA SB380H



SHA UPCR6045

#### 10% FORMALIN POISON S

SHA UPCR1319



SHA UPCR1320



SHA UPCR1318

CAUTION: Contains FORMALDEHYDE. Toxic by inhalation and if swallowed. Irritating to the eyes, respiratory system and skin. May cause sensitization by inhalation or skin contact. Risk of serious damage to eyes. Potential cancer hazard. Repeated or prolonged exposure increases the risk.

SHA UPCR5008

## **GLUTARALDEHYDE**

SHA UPCR7017

95% ALCOHOL FOR PAP SMEARS

SHA UPCR5007




SHA SB500HT

#### CAUTION LABELS FOR PATHOLOGY SPECIMEN STORAGE

ITEM #	DESCRIPTION	SIZE	ITEM #	DESCRIPTION	SIZE
SHA UPCR1319	10% Formalin/Poison	5/16" x 11/4"	SHA SB385H	Caution / Contains Formaldehyde	1" x 2"
SHA UPCR1320	Pathologist Review Needed	5/16" x 11/4"	SHA SB388H	Xylene / Flammable Liquid	1" x 2"
SHA UPCR1318	<sup>10</sup> % Formalin	5/16" x 11/4"	SHA SB389H	Alcohol / Flammable Liquid	1" x 2"
SHA UPCR7017	Glutaraldehyde	1" x 3"	SHA SB381H	10% Formalin / Formaldehyde	1" x 2"
SHA UPCR5009	Do Not Drink	<sup>3</sup> / <sub>4</sub> " x <sup>2</sup> "	SHA SB380H	10% Formalin / Formaldehyde	13/4" × 21/2"
SHA UPCR5007	95% Alcohol For Pap Smears	<sup>3</sup> / <sub>4</sub> " x <sup>2</sup> "	SHA UPCR6046	10% Formalin / Formaldehyde	1" x 21/4"
SHA UPCR6044	Fixative	1" x 21/4"	SHA UPCR6042	Caution / Contains Formaldehyde	1" x 21/4"
SHA SB382H	10% Formalin / Poison	1" x 2½"	SHA UPCR6045	Sputum Cytology / Poison	1" x 21/4"
SHA UPCR6047	PAP Cytology Fixative	1" x 21/4"	SHA SB500HT	<sup>10</sup> % Neutral Buffered Formalin	³⁄₄" x 60yd
SHA UPCR6043	Cytology Fixative Poison	1" x 21/4"	SHA SB387H	Caution / Releases Formaldehyde	1 <sup>3</sup> / <sub>4</sub> " × 2 <sup>1</sup> / <sub>2</sub> "
SHA SB374H	10% Formalin	1" x 21/4"	SHA UPCR5008	Caution / Contains Formaldehyde	<sup>3</sup> / <sub>4</sub> " x <sup>2</sup> "






SHA UPCR30



SHA UPCR37



CHEMICAL HAZARD IDENTIFICATION LABELS




SHA UPCR36








SHA UPCR33



**SHA UPCR32** 



SHA UPCR7002







SHA SBH31



SHA UPCR34



SHA UPCR40



SHA UPCR7003



#### **ACETONE**

SHA UPCR42



SHA UPCR43



SHA UPCR44



SHA UPCR45



SHA UPCR46



**SHA UPCR47** 



SHA UPCR48



SHA UPCR49

ITEM #	DESCRIPTION	SIZE	ITEM #	DESCRIPTION	SIZE
SHA UPCR30	Health / Flammability / Reactivity	2" x 2"	SHA UPCR45	Isopropyl Alcohol	1" x 3"
SHA UPCR31	Reagent / Prep. Date / Exp. Date	1½" x 3"	SHA UPCR46	Methyl Alcohol	1" x 3"
SHA UPCR32	Caution / Contains Formaldehyde	1" × 3"	SHA UPCR47	Nitric Acid	1" x 3"
SHA UPCR33	Caution / Contains Formaldehyde	1" x 21/4"	SHA UPCR48	Sodium Hydroxide	1" x 3"
SHA UPCR34	2.5% Glutaraldehyde / EM Fixative	1" x 21/4"	SHA UPCR49	Sulfuric Acid	1" x 3"
SHA UPCR36	Health / Fire / Specific Hazard	2" x 3"	SHA UPCR50	Xylene	1" x 3"
SHA UPCR37	Chemical Name / MSDS # / Date	21/4" x 21/4"	SHA UPCR7002	Chemical Hazard Symbol	1" x 3"
SHA UPCR40	Chemical Name / Common Name	1" x 3"	SHA UPCR7003	Chemical / Owner / Date	1" x 3"
SHA UPCR41	Acetic Acid	1" x 3"	SHA UPCR9403	Hazardous Materials	2" x 3"
SHA UPCR42	Acetone	1" x 3"	SHA SBH30*	Emergency Guide	2½" x 4"
SHA UPCR43	Ethyl Alcohol	1" x 3"	SHA SBH31*	Chemical Hazard Symbol	2" x 2"
SHA UPCR44	Hydrochloric Acid	1" x 3"	SHA SBH32*	Chemical Hazard Symbol	1" × 1"

LABELS PACKAGED 500 PER ROLL UNLESS NOTED. (\*) = 1,000 LABELS PER ROLL.



## TANNER SCIENTIFIC® BLUTOUCH® NITRILE EXAM GLOVES

- · Powder-Free, Latex-Free, Non-Sterile
- Made of Synthetic Rubber
- Disposable (Single Use Only)
- · Comfortable Fit and Feel
- · High Barrier Protection
- · Superior Solvent Resistance
- Textured Fingertips for Better Grip
- Exceeds ASTM Standards

TNR 700S	Small (6.5-7)	100/Box
TNR 700M	Medium (7.5-8)	100/Box
TNR 700L	Large (8.5-9)	100/Box



#### **NITRILE**



Product was **FENTANYL TESTED** 

#### **NITRIDERM® ULTRA BLACK NITRILE EXAM GLOVES**

- · Powder-Free, Non-Latex, Non-Sterile
- · Named a true medical grade exam glove
- #1 choice for law enforcement and corrections
- · Medium modulus and high tensile strength
- Popular black color masks a variety of stains
- · Textured for excellent grip

INN 187050	X-Small	100/Box
INN 187100	Small	100/Box
INN 187200	Medium	100/Box
INN 187300	Large	100/Box
INN 187350	X-Large	1,000/Case



#### NITRIDERM® ULTRA BLUE

- Powder-Free, Non-Latex, Non-Sterile
- Exceptionally soft
- · Very comfortable for extended periods
- · Low modulus formulation
- COLOR: LIGHT BLUE

INN 157050	X-Small	100/Box
INN 157100	Small	1,000/Case
INN 157200	Medium	100/Box
INN 157300	Large	100/Box
INN 157350	X-Large	100/Box



#### PRODERM™ XR

- Powder-Free, Non-Latex, Non-Sterile
- · Thinner and lightweight
- Soft, comfortable, easy to don
- Extraordinary tactile sensitivity yet equally extraordinary strength
- · COLOR: BLUE

INN 158050	X-Small	250/Box
INN 158100	Small	250/Box
INN 158200	Medium	250/Box
INN 158300	Large	250/Box
INN 158350	X-Large	250/Box

#### **NITRILE**





- · Powder-Free, Non-Latex, Non-Sterile
- · Exceptionally soft, thin-film is very comfortable for extended periods
- Slightly tacky, textured fingers provide excellent gripping properties
- · COLOR: WHITE

#### **PULSE® NITRILE**

- · Powder-Free, Non-Latex, Non-Sterile
- · Thinner glove for incredible elasticity
- Finely textured fingertips for unequaled tactile sensitivity and excellent gripping
- Automated flat packing allows gloves to dispense easily
- **COLOR: DARK LAVENDER**

#### NITRIDERM® ULTRA ORANGE

- · Powder-Free, Non-Sterile, Textured
- Ideal for EMS, police and fire, hazmat functions, spills, etc
- Extremely elastic and comfortable while providing the tensile strength
- · COLOR: ORANGE

#### KIMBERLY-CLARK® PURPLE

- · Powder-Free, Non-Latex, Non-Sterile
- · Maximum coverage and protection
- · Excellent dexterity and sensitivity
- Cleared for use in chemotherapy
- · Textured fingertips
- · Beaded cuff
- COLOR: PURPLE

INN 167050	X-Small	100/Box
INN 167100	Small	100/Box
INN 167200	Medium	1,000/Case
INN 167300	Large	100/Box

INN 177052	X-Small	200/Box
INN 177102	Small	2,000/Case
INN 177202	Medium	200/Box
INN 177302	Large	200/Box
INN 177352	X-Large	200/Box

INN 187100	Small	100/Box
INN 187200	Medium	100/Box
INN 187300	Large	100/Box

KIM 55080	X-Small	100/Box
KIM 55081	Small	100/Box
KIM 55082	Medium	100/Box
KIM 55083	Large	100/Box
KIM 55084	X-Large	100/Box

#### **LATEX**



#### **PULSE® LATEX**


- · Powder-Free, Non-Sterile, Textured
- Engineered to be slightly thinner for superior tactile sensitivity
- · Online chlorination creates superior gripping ability for a very soft feel
- · Excellent donning properties and a more comfortable fit
- · COLOR: WHITE

#### INN 151200 Medium 100/Box INN 151300 Large 100/Box INN 151350 X-Large 1.000/Case

Small

INN 151100

#### **PRODERM™**



- · Engineered to be slightly thinner for incredible tactile sensitivity
- Polymer bonded interior and lightly textured surface
- Unique polymer process creates a glove easy to don without using chlorination
- · Very comfortable for extended wear
- · COLOR: WHITE

INN 155050	X-Small	100/Box
INN 155100	Small	100/Box
INN 155200	Medium	100/Box
INN 155300	Large	100/Box
INN 155350	X-Large	100/Box



100/Box



#### OMNIMED® GLOVE BOX HOLDERS

- Ideal for medical, laboratory or any area where cross contamination occurs or protective gloves are needed
- Mount vertically or horizontally on a variety of surfaces
- · Easy to clean, load and mount
- All models accommodate single glove boxes
- Allow easy access to disposable apparel in exam rooms, hospitals, laboratories and more
- Available in a variety of sizes and styles to suit any volume of usage

#### METAL

- Stainless steel construction built for heavy duty use
- Spot welded for a sleek finish to withstand years of use

OMN 305303	Triple,	10"H x 16"W x 33/4"D	2/Pack

#### WIRE

· Heavy gauge, powder coated, steel wire construction

OMN 305375	Triple, 16"H x 11½"W x 4½"D	2/Pack

## HEA HD23458





HEA HD23457

#### **WIRE GLOVE BOX HOLDERS**

- Free up benchtop space by mounting holders on wall
- Sturdy, epoxy-coated steel wire racks
- Fit most major brands of glove boxes
- · Vertical or horizontal mounting
- · Mounting screws included

HEA HD23457	Single, 140 x 209 x 108mm	Each
HEA HD23458	Triple, 250 x 454 x 108mm	Each

#### **SHOE COVERS**

- · Footwear protection with fluid repellant, slip-on covers
- Sewn seams provide increased tear strength
- Elastic opening for easy on-and-off convenience
- · Non-conductive and non-sterile

AMD 8008	Skid	100/Box
AMD 8007	Non-Skid	100/Box



#### **UNISEX LAB COATS**

- 40" unisex lab coat made of polyester/cotton
- (2) Roomy Patch Pockets, (1) Chest Pocket with Pencil Stitch, Vented Back
- **SIZES:** XS | S | M | L | XL | 2XL | 3XL | 4XL | 5XL

DLU 1346 Unisex White Lab Coat Each

#### **UNISEX V-NECK TUNIC TOPS\***

- · Made of polyester/cotton material
- 1 Chest/Breast Pocket with Pencil Slot, Dolman Sleeve, Pull Over Head, Center Back Length: 27½"
- **SIZES:** XS | S | M | L | XL | 2XL | 3XL | 4XL | 5XL

**DLU 4777** Unisex V-Neck Tunic Tops Each

#### **UNISEX DRAWSTRING PANTS\***

- · Made of polyester/cotton material
- Adjustable Webbed Drawstring Waist, (1) Back Pocket, (1) Cargo Pocket with Cell Phone Pocket, 31" Inseam
- **SIZES:** XS | S | M | L | XL | 2XL | 3XL | 4XL | 5XL

**DLU 4100** Unisex Drawstring Pants Each

\*COLORS: Black, Navy, Pewter, Royal, White, Ceil, Wine, Khaki, Red, Teal, Grey, Chocolate, Galaxy Blue, Hunter, Caribbean Blue, Turquoise, Raspberry, Shocking Pink, Grape, Olive, Eggplant, Orchid, Surgical Green, Sage





#### **PROTECTION PLUS® LAB COATS**

#### UNISEX LONG LENGTH LAB COATS

- · Lightback design promotes air flow for user comfort
- · Maximum fluid resistance front and back
- Pockets: double-stitched pockets (1 Breast, 2 Lower) with BarTack reinforced at top corners for durability
- Cuffs: stain resistant, low pill, spun woven polyester, comfortable, 41" long
- Tapered Sleeves prevent catching when reaching, cuffs overgloved provides protection for forearms
- · Heavy Duty E-Z Snap Closures
- Side Vents for access to pant pockets
- · Hang Loop is located on the back of neckline
- Tested 200+ Commercial Wash and Dry Cycles
- SIZES: XS | S | M | L | XL | 2XL | 3XL | 4XL | 5XL
- · COLORS: Black, Ceil, Emerald Green, Royal, White

DLU DL170BJ ProtectionPlus® Lab Coat Each



#### **KNEE LENGTH LAB COATS**

#### **PROTECTION LAB COAT (3460)**

- Constructed of spunbonded polypropylene (SPP) with exceptional cloth-like feel and breathability
- **SIZES:** S | M | L | XL | 2XL | 3XL
- · COLORS: Sky Blue, White

#### **EASY-BREATHE SMS LAB COAT (3560)**

- **SIZES:** S | M | L | XL | 2XL | 3XL | 4XL | 5XL
- COLORS: Ceil Blue, Medical Blue, Purple, Raspberry, Teal, White

#### EXTRA-SAFE SMS LAB COAT (3660)\*\*

- · Anti-static treatment prevents clinging
- **SIZES:** XS | S | M | L | XL | 2XL | 3XL | 4XL | 5XL

#### **LIQUIDGUARD LAB COAT (3960)**

- · Microporous breathable film material
- Impervious yet breathable, highly anti-static treated
- **SIZES:** S | M | L | XL | 2XL | 3XL | 4XL
- · COLOR: White

VAL 3460	Protection Lab Coat	10/Pack
VAL 3560	Easy-Breathe SMS Lab Coat	10/Pack
VAL 3660	Extra-Safe SMS Lab Coat	10/Pack
VAL 3960	LiquidGuard Lab Coat	10/Pack

#### HIP LENGTH LAB JACKETS

#### **EASY-BREATHE LAB JACKET (3530)**

- Durable medium weight unisex lab jacket made for protection and comfort with 3-layer SMS fabrics
- **SIZES:** S | M | L | XL | 2XL | 3XL | 4XL | 5XL
- COLORS: Ceil Blue, Medical Blue, Purple, Raspberry, Teal, White (No Pockets)

#### EXTRA-SAFE LAB JACKET (3630)\*\*

- Durable unisex lab jacket made for protection and comfort with soft 3-layer SMS fabrics
- · Anti-static treatment prevents clinging
- **SIZES:** XS | S | M | L | XL | 2XL | 3XL | 4XL | 5XL

VAL 3530	Easy-Breathe Lab Jacket	10/Pack
VAL 3630	Extra-Safe Lab Jacket	10/Pack

\*\*COLORS: Aqua, Black, Blueberry, Ceil Blue, Coffee, Coral Pink, Cranberry, Deep Purple, Deep Sea Blue, Hot Pink, Lavender, Light Pink, Medical Blue, Multicolor, Navy Blue, Olive, Orange, Pink, Purple, Plum, Raspberry, Red, Royal Blue, Sky Blue, Steel Blue, Tan, Teal, True Blue, Violet Purple, White





#### **DUPONT™ TYVEK™ BIB APRON**

- Nonwoven Tyvek™ fabric helps resist penetration by particles and other contaminants
- Inherent barrier stops submicron particles even after abrasion
- · Provides enhanced frontal protection
- · Bound seams around neck and ties at waist
- ONE SIZE FITS MOST
- For use with lead and asbestos abatement/remediation, general maintenance/operations, spray painting, and general clean-up
- Should not be used around heat, flames, sparks or in potentially flammable or explosive environments

FIS 0136183

Bib Apron, Universal Size, White, Length: 36" (69cm) 100/Case

## WHITE APRONS DIMENSIONS: 27" X 43"

Disposable, non-woven fabric impervious to fluids

VAL 1550 Polycoated 25/Pack



#### ANTI-FOG SURGICAL MASK

- Nosepieces help conform to the contours of the face
- · Pleat-style with ties and film strip

KIM 49235 Surgical Mask, Green

300/Case

#### SPLASH SHIELD° FACE SHIELDS

- Anti-fog and optically clear
- · Soft strap for comfort
- · Latex free, disposable
- **Dimensions:** 13"W x 5"L
- Impenetrable foam barrier for wide coverage
- Maximum face protection hospitals, labs, dental and EMT use



SPL 4511 Half Face 24/Pack



#### **SAFETY GOGGLES**

- Scratch-resistant, polycarbonate lenses provide 99% protection against harmful UV rays
- Polyvinyl chloride (PVC) frame and nose piece
- Elastic headband made of polyester and rubber thread
- Clear, anti-fog lens for indoor applications
- Lens Dimensions: 144mm x 68mm x 2.0mm
- Ventless: Protects against odor vapors
- Chemical: Vent caps restrict influx of liquids

PYR G200T	Ventless Goggles	Each
PYR G204T	Chemical Goggles	Each

#### VENTURE











PYR SB1810S

10S PYR SB1890S

PYR SB1830S

PYR SB1860S

PYR SB1840S

#### ONIX











PYR SB4910S

PYR SGR4910ST

PYR SP4910S

PYR SB4920S

PYR SR4980S

#### **ONIX® EYEWEAR**

- · Cushioned temple tips provide non-slip comfort
- Adjustable nose piece and ratcheting temples create different angles for perfect view and fit
- Lightweight frame and suspended lens design provide a perfect fit for a variety of users
- Suspended lens provides stylish look, comfort fit, panoramic view and maximum protection
- Scratch resistant polycarbonate lens for UV protection

PYR SB4910ST

Clear, Anti-Fog

Each

#### **VENTURE II® EYEWEAR**

- · Temples adjust to four different lengths
- · Nose buds secure the glass and prevent slipping
- Curved lens provides full side protection without distracting seams or traditional type side shields
- Scratch resistant polycarbonate lens for UV protection
- Colored lenses offer different properties gray stops excessive glare, orange gives low light image resolution
- Passes MIL-PRF 32432 high velocity impact standards

PYR SB1810ST

Clear, Anti-Fog

Each

#### MINI ZTEK







PYR S2520SN





INTRUDER











PYR S4110SUC

PYR S4120S

PYR SN4110S

PYR SR4110S

PYR SB4110S

#### MINI ZTEK® EYEWEAR

- Design is better suited for slimmer facial sizes
- Integrated nose piece
- Soft, non-slip rubber temple tips for non-binding fit
- Wrap-around single lens provides full panoramic view
- Scratch resistant polycarbonate lens for UV protection

PYR S2510SNT Clear, Anti-Fog

#### INTRUDER® EYEWEAR

- Lightweight, frameless protection
- Superior comfort and fit
- Scratch resistant polycarbonate lens for UV protection
- · Mini Intruder versions are suited for smaller facial size
- · Perfect for work areas where color coding is necessary

PYR S4110ST

Clear, Anti-Fog

Each

#### **CONTACT SALES REPRESENTATIVE FOR EYEWEAR OPTIONS!**

Each

#### A

<del>-</del> -	
Absorbance Pads	135
Absorbent Pads, Formaldehyde	152
Acetic Acid	115
Acetone	105
Acid Fast Bacteria Stain Kit	119
Acid, Acetic	115
Activated Charcoal	24
Adapter, Chuck	57
Adsorber, Fume	22
Agent, Mold Release	55
Agents, Clearing	110
Air Filtration System	25
Alcian Blue Stain Kit	119
Alcohol Recycling System	25
Alcohols	105
Aligner, Microtome	32
Aluminum Ammonium Sulfate	115
Aluminum Chloride	115
Amyloid Stain Kit	119
Applicator Sticks	49
Aprons	175
Ascorbic Acid	115
Automatic Multistainer	16
Azure A Powder	115
В	
Bags, Biohazard	34
Bags, Biopsy	39
Bags, Specimen	34
Bags, Tea	41
Base Molds	55

Basic Fuchsin	115
Baskets, Histology Cassette	77
Bath, Tissue	26
Beakers	126, 130
Biohazard Bags	34
Biohazard Labels	165
Biopsy Bags	39
Biopsy Capusles	41
Biopsy Containers	129
Biopsy Foam Pads	40
Biopsy Paper	40, 41
Biopsy Wraps	40
Blade Removal	44
Blades, Microtome	53, 54
Blades, Razor	46
Blades, Scalpel	43, 44
Blanket, Fire	139
Bleach	145
Block Cooling System	61
Bloodbank Saline	118
Board, Cutting	41, 42
Board, Dissection	42
Boats, Weighing	160
Boric Acid	115
Bottles, Reagent	131
Bottles, Wash	128
Boulin Fixative	115
Boxes, Lock	36
Boxes, Slide	94, 95
Brilliant Cresyl Blue	115
Brushes, Cytology	102
Brushes, Frozen Section	57

Brushes, Paraffin59	Coats, Lab	173, 174
Buffer Solution115	Collection Kits	34
	Colloidal Iron Stain Kit	119
C	Congo Red	115
Cabinets, Slide93	Containers Prefilled	35
Cabinets, Storage133	Containers, Biopsy	129
Calcium Stain Kit119	Containers, Sharps	141
Cans, Waste131	Control Slides	87, 88
Capsules, Biopsy41	Convection Ovens	30
Carbol Fuchsin Stain115	Cooling System, Block	61
Carbon Filters23	Coplin Jar	65
Carboys126	Copper Stain Kit	120
Carono Fixative115	Copper Sulfate	115
Cassette Printer19	Countertop Paper	135
Cassette Storage Drawers77	Coverglass	85, 86
Cassette Storage System78	Covers, Shoe	172
Cassettes, Histology69-76	Coverslipping Film	85
Caution Labels166, 167	Cryostats	2-6
Cedarwood Oil115	Crystal Violet	115
Centrifuge Tubes103	Cutting Board	41, 42
Cervical Scrapers102	Cylinders	126, 130
Chairs	Cytology Brushes	102
Charcoal24	Cytology Funnels	101
Chemical Hazard Labels167, 168	Cytology Preservatives	100
Chromic Acid115	Cytology Slide Processor	17, 99
Chuck Adapter 57	Cytology Stain	102
Chuck, Specimen31		
Chucks	D	
Circles, Qualitative63	Davidson Fixative	115
Cleaners144, 145	Decalcifiers	111, 115
Clearing Agents110	Detergents	142, 143
Clocks159	Dipper, Staining	98

Disc, Leveling	32	Flasks	130
Discs, Specimen	31	Folders, Slide	91, 92
Dishes, Staining	64, 65, 98	Fontana Stain Kit	120
Disinfectants	144, 145	Forceps	38
Dispenser, Paraffin	29, 30	Formaldehyde	116
Dissection Board	42	Formaldehyde Absorbent Pads	152
Drawers, Cassette Storage	77	Formalin Fixative	109
Drawstring Pants	173	Formalin Pads	152
Dye, Tissue Marking	47-50	Formalin Recycling System	25
		Formic Acid	116
E		Freeze Spray	61
Embedding Centers	10-13	Freezing Media	60, 61
Embedding Well Bars	56, 57	Frozen Section Brushes	57
Eosin Y	113, 116	Fume Adsorber	22
Eosinophil Stain Kit	120	Fume Hoods	22
Ethyl Acetate	116	Funnels	127
Extractor, Heat	32	Funnels, Cytology	101
Eyewash Station	139		
Eyewear	175, 176	G	
		Gallego Solution	116
F		Gauze Sponges	135
Face Shields	175	Gelatin	116
Ferric Chloride	116	Gentian Violet	116
Film, Coverslipping	85	Giemsa Stain	116
Filter Paper	63	Giemsa Stain Kit	120
Filters, Carbon	23	Glove Box Holders	172
Filtration System, Air	25	Gloves	170, 171
Fire Blanket	139	Glutaraldehyde	116
First Aid	147	Glycerin	116
Fixative Spray	102	GMS Stain Kit	120
Fixative, Formalin	109	Goggles	175
Fixatives	114	Gold Chloride	116

Gomori Trichome116	J
Gram Stain116	Jackets, Lab174
Gram Stain Kit121	Jar, Coplin65
Gravity Convection Ovens30	Jars, Pap103
Grossing Tissue Paper39	Jenner Stain116
н	Jones Stain Kit121
H. Pylori Stain Kit121	K
Hand Sanitizers137	Karnovsky Fixative116
Handles, Scalpel43	Kimwipes®137
Hazardous Disposal140	Kit, Slide Staining67
Heat Extractor32	Kits, Pap Smear102
Hematoxylin112, 116	Kits, Slide Filter100
Hemotoxylin & Eosin Stain Kit121	Kits, Specimen Collection34
Histology Cassett Baskets77	Knives46
Histology Cassettes69-76	
Histo-Wrap39	L Control of the Cont
Holders, Glove Box172	Lab Coats173, 174
Holders, Slide92	Labels, Biohazard165
Hoods, Fume22	Labels, Caution166, 167
Hot Foil Tape31	Labels, Chemical Hazard167, 168
	Labels, Printer18, 163, 164
Hydrochloric Acid116	Latex Gloves171
Hydrogen Peroxide144	Leishman Stain116
	Leveling Disc32
	Linear Stainers21
Immersion Oil116	Lithium Carbonate116
Ink49	Lock Boxes36
Iodine116	Lubricating Oil31
Iron Stain Kit121	Luxol Fast Blue Stain Kit121

## M

Market College In	117
Magnesium Carbonate	
Mailers, Paraffin Block	
Mailers, Slide	
Malachite Green	117
Markers	51
Markers, Slide Index	95
Markers, Tissue Index	77
Masks	175
Mats	154-156
Media, Mounting	106-109
Mesh Strainers	65
Metanil	117
Methanol	105
Methylene	117
Methylene Blue Stain Kit	122
Michel	117
Microscope Slides	80-85
Microscopes	27, 28
Microtome Aligner	32
Microtome Blades	53, 54
Microtome Leveling Disc	32
Microtomes	7-9
Mold Release Agent	55
Molds, Base	55
Monitors Vapor	146
Mounting Media	106-109
Mounting Media, Toulene	107
Mounting Media, Water	108
Mounting Media, Xylene	106
Movat Stain Kit	122
Multistainer, Automatic	16

### N

Needles, Teasing	43
Neutral Buffered Formalin Prefilled	Containers 35
Neutral Red	117
Neutralizers	149, 150
Nitric Acid	117
Nitrile Gloves	170, 171
Nuclear Fast Red	117
Nuclear Fast Red Stain Kit	122

#### C

Oil Red O	117
Oil Red O Stain Kit	122
Oil, Lubricating	31
Optical Lens Paper	41
Orange G	117
Orcein	122
Ovens	30

#### -

Pads, Absorbance	135
Pads, Biopsy Foam	40
Pads, Formalin	152
Pants, Drawstring	173
Pap Jars	103
Pap Smear Kits	102
Papanicolaou Stain	117
Papanicolaou Stain Kit	122
Paper, Biopsy	40, 41
Paper, Countertop	135
Paper, Filter	63
Paper, Grossing	39

Paper, Optical Lens4	1 Processor, Slide	17, 99
Paraffin58	Processor, Tissue	14
Paraffin Block Mailers55	5	
Paraffin Brushes59	9 Q	
Paraffin Dispenser29, 30	Qualitative Circles	63
Paraffin Repellant114	4 Quik-Dip® Stain	98
Paraffin Scrapers59	9	
Paraffin Wax Trimmer20	R	
Parafilm59	Rack, Slide Staining	67
Peanut Oil Xylenes118	Racks, Slide Tray	91
Pencils5	1 Razor Blades	46
Pens5	1 Reagent Bottles	131
Periodic Acid117	Recycling System, Alcohol	25
Periodic Acid Schiff Stain Kit122, 123	Recycling System, Formalin	25
pH Tests153	Repellant, Paraffin	114
Phenol117	7 Reticulum Stain Kit	123
Phloxine B117	7 Rhodanine	117
Phosphomolybdic Acid117	7	
Picric Acid117	7 <b>S</b>	
Picro-Sirius Red Stain Kit123	Saccomanno Fixative	118
Pins43	Safranine O	118
Pipets126	S Sample Storage Trays	103
Pneumocystis Stain Kit123	3 Sanitizers	137
Potassium Bromide117	7 Scalpel Blades	43, 44
Potassium Ferricyanide117	7 Scalpel Handles	43
Potassium Hydroxide117	7 Scalpels	45
Potassium Iodide117	7 Schiff Stain	118
Prefilled Containers	5 Scissors	38
Preservatives, Cytology100	Scott's Tap Water Substitute	112, 118
Printer Labels18, 163, 164	Scrapers, Cervical	102
Printer, Cassette19	Scrapers, Paraffin	59
Printer, Slide18, 19, 163	3 Scrubs	173

Sets, Staining	98	Softener	114
Sharps Containers	141	Solidifiers	150
Shields, Face	175	Spacers, Sponge Block	93
Shoe Covers	172	Specimen Bags	34
Silver Nitrate	118	Specimen Chuck	31
Slide Boxes	94, 95	Specimen Collection Kits	34
Slide Cabinets	93	Specimen Discs	31
Slide Filter Kits	100	Spill Control Kits	151
Slide Folders	91, 92	Sponge Block Spacers	93
Slide Holders	92	Sponges, Gauze	135
Slide Index Markers	95	Spray, Fixative	102
Slide Mailers	89, 90	Spray, Freeze	61
Slide Printer	18, 19, 163	Stain, Cytology	102
Slide Processor	17, 99	Stain, Quik-Dip®	98
Slide Stainer	15, 16, 21	Stainer, Slide	15, 16, 21
Slide Staining Kit	67	Stainers, Linear	21
Slide Staining Rack	67	Staining Dipper	98
Slide Staining System	67	Staining Dishes	64, 65, 98
Slide Storage	91-96	Staining Sets	98
Slide Tray	66	Staining Trays	65, 66
Slide Tray Racks	91	Steiner Stain Kit	124
Slide Warmer	20	Sticks, Applicator	49
Slides	80-85	Stools	132
Slides, Control	87, 88	Storage Cabinets	133
Sodium Bicarbonate	118	Storage, Slide	91-96
Sodium Carbonate	118	Strainers, Mesh	65
Sodium Chloride	118	Substitute, Water	112
Sodium Citrate	118	Substitute, Xylene	110
Sodium Hydroxide	118	Sudan	118
Sodium Phosphate	118	Sufosalicylic Adid	118
Sodium Sulfate	118	Sulfuric Acid	118
Sodium Thiosulfate	118		

## **LOOKING TO SET UP A NEW LAB?**

## Mercedes Scientific is here to help!

We currently provide a full line of **NEW and Refurbished Laboratory Equipment!**On-site installation, training and extended warranties are available for most models.

## TYPES OF EQUIPMENT:

- Tissue Embedding Units
- Microtomes
- Cryostats
- Autostainers
- Manual Staining Systems
- Coverslippers
- Paraffin Dispensers
- Tissue Flotation Baths
- Slide Warmers/Dryers
- Microscopes
- Slide/Cassette Labelers

Time-WVVI

• Fume Hoods

AND MUCH MORE!

#### **INCLUDES:**

- Professional Packaging and Full-Value Insured Shipping
- Full Warranty on ALL Equipment\*
- Extended Service Contracts
- Experienced Techs for Timely and Quality Repairs if needed
- Lease Options



TN55 CRYOSTAT (NEW)

FOR MORE INFORMATION CALL 800.331.2716

\*PLEASE NOTE: Not all items come with parts and labor warranties. Ask for details on specific refurbished items. Financing Options Available.

TISSUE-TEK° VIP PROCESSOR (REFURBISHED)



12210 Rangeland Pkwy Lakewood Ranch, FL 34211

PHONE 800.331.2716
FAX 800.359.8807

**WEB MercedesScientific.com** 



FREE HENRY, HAZEL, OR POP-CORN

With Orders of \$99 ore more



MONTHLY RAFFLES ENTER-TO-WIN

Yes. Every month!



**DOG-FRIENDLY WORKPLACE** 

Why? Because they are awesome.

CATALOG # RH12

**Mercedes Scientific** offers full product lines and alternate product options. Call today to ask about the best product selection that suit your needs.